01

ALYA

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
CONTENT

1	SAFETY REQUIREMENTS	36
1.1	SYMBOLS USED.....	36
1.1.1	Symbols used in the manual	36
1.1.2	Symbols on the labels	36
1.1.3	Symbols on the packaging.....	37
1.2	PLANNED USE.....	37
1.3	PLANNED USER.....	37
1.3.1	Professional qualification:	37
1.3.2	Minimum skills.....	37
1.3.3	Experience	37
1.3.4	Possible user handicaps	37
1.4	GENERAL STANDARDS AND MAIN WARNINGS.....	38
1.5	STORAGE AND USE: ENVIRONMENTAL PROVISIONS.....	39
1.6	REQUIREMENTS FOR ELECTROMAGNETIC COMPATIBILITY	39
2	GENERAL FEATURES.....	41
2.1	DESCRIPTION OF THE PRODUCT	41
2.2	DESCRIPTION OF THE PARTS	42
2.3	DEVICE IDENTIFICATION.....	42
3	DEVICE INSTALLATION.....	44
3.1	DIMENSIONS	44
3.2	COMPLETE DENTAL LIGHT	45
3.2.1	Electrical requirements.....	45
3.2.2	Complete combined version light assembly	46
3.2.3	Remote cable connection.....	46
3.2.4	Complete dental light connection to the ceiling, wall or floor	47
3.2.5	Installation applications.....	47
3.3	HEAD.....	56
3.3.1	Mechanical requirements.....	56
3.3.2	Electrical requirements.....	57
4	INSTRUCTIONS FOR USE.....	58
4.1	SWITCH ON AND OFF	58
4.1.1	Adjustment:	58
4.1.2	Light / Head WITH PROXIMITY.....	58
4.1.3	Light / Complete light with Theia Tech / “ALYA” head WITH REMOTE CONTROL.....	59
4.1.4	SYNCRONIZATION command.....	59
4.2	JOYSTICK ALYA JOYSTICK LEVER ASSEMBLY	61
5	ORDINARY MAINTENANCE.....	61
6	CLEANING	61
6.1	PARABOLA CLEANING	61
6.2	CLEANING THE HEAD.....	62
6.3	ARM CLEANING	62
7	STERILISING THE HANDLES	62
7.1	Removing the handles	62
7.2	Decontamination and disinfectant.....	62
7.3	Sterilisation.....	63
8	PERIODIC CHECKS.....	63
9	ACOUSTIC ALARM SIGNALS.....	63
9.1	Acoustic Signals.....	63
9.2	Guide to problems	64
10	Technical features.....	65
10.1	WIRING DIAGRAMS.....	66

1. SAFETY REQUIREMENTS

Dear Client,

FARO hopes you enjoy your work with the new high quality light. For safe work and to take full advantage of the performance of the product, carefully read this manual before using the device. In particular, follow all the warnings and the notes given.

1.1 SYMBOLS USED

1.1.1 Symbols used in the manual

	WARNING
The paragraphs marked with this symbol contain instructions that must be carefully followed to avoid damaging the device, harming the operator or the patient.	
	WARNING
These instructions warn you that you must pay attention to avoid situations that could damage the device.	
	FORBIDDEN
This icon highlights what you should not do to avoid damaging the device.	
	NOTES
This icon supplies information that allows you to use the device more efficiently.	

1.1.2 Symbols on the labels

The data plate is fixed:

- for the complete light: on the rear arm
- for the head: under the heat sink cover

and outlines the following data:

Serial Number (SN): year (YY) / range of origin (LD for dental light - TE for head only) plus the progressive number (NNNNNN) e.g.: SN14LD000001 for the complete light SN 14TE000001 for the head.

The following standardised symbols are also present:

	Read the use instructions. Supply electricity
	Manufacturing symbol pursuant to Directive 93/42/EEC
	The use instructions include safety warnings
	WEEE equipment in compliance with the Directive 1012/19/EC. Dispose of the product in compliance with this directive.
	Double insulation. Class 2 device against electrical risk
	Serial Number
	Can be sterilised with humid heat at 134°C
	Use the device at a temperature between 10°C and 40°C
	Use the device at pressure between 800mbar and 1060mbar
	Use the device at relative humidity between 30 RH and 75RH
	Symbol to adjust light intensity
	Symbol to switch on the light
	[Symbol to switch on/off the light on the rear arm]

1.1.3 Symbols on the packaging

HIGH

FRAGILE

DO NOT WET

DO NOT ROLL

DO NOT USE HOOKS


500 kg max


MAXIMUM STACKABLE WEIGHT

**STORAGE TEMPERATURE
TRANSPORT CONDITIONS**

RELATIVE HUMIDITY

ATMOSPHERIC PRESSURE

RECYCLABLE CARDBOARD

1.2 PLANNED USE

The device is intended for use solely in dental practices by dentists and hygienists to light the operative and intervention site in the treatment of oral cavity and dental treatments. In normal use, the device is positioned over the patient's body at a distance of 700mm, the distance for which the lighting features were designed.

The treated patients can be of all ages with typical dental pathologies.

1.3 PLANNED USER

The planned user is a dentist or dental hygienist.

1.3.1 Professional qualification:

- Degree in medicine with dentistry specialisation
- Degree in dentistry
- Degree in dental hygiene

1.3.2 Minimum skills

- Those planned for the professional qualification
- Understanding of language: Those acquired for the professional qualification

1.3.3 Experience

- Those outlined to conduct the profession

1.3.4 Possible user handicaps

- For use, complete use is necessary of an upper limb.
- Visual faculty compatible with the profession

1.4 GENERAL STANDARDS AND MAIN WARNINGS

- The device can be applied to the dental unit, but also be installed on specific applications. The device can be powered both by the dental unit and by a power supply unit connected directly to the mains. See the specific installation paragraph.
- The device does not have Essential Performances for which inadequacy of the device performance does not prejudice patient safety.
- The device does not provide life support.
- The device must be clean before use (see Device Cleaning paragraph).
- The packaging of the light is suitable to adequately protect it from penetration of external agents.

	Warnings against electrical danger and fire Do not use the light in the event parts are damaged. Installation of the device must only be carried out by qualified staff. The dental light must be installed on a specific control and power supply device, such as dental units, or with an electrical system that meets standard IEC 60364-1 and “national installation regulations for electrical systems in premises for medical use”. The appliance must be installed with an omnipolar separation device from the mains and compliant with Standard IEC 61058-1. Installation and maintenance of device conformity with the standard IEC 60601-1 is the responsibility of the installation technician or the manufacturer of the combined units. Check the power supply voltage, indicated on the data plate, corresponds to the mains voltage. Do not carry out any maintenance on the light when the power supply is inserted: then disconnect the power supply from the mains before intervening.
	Warning against danger of wear of the mechanical parts and falling suspended weights Do not use detergents containing the following to clean plastic parts: AMMONIUM HYDROXIDE - SODIUM HYDROXIDE - METHYLENE CLORIDE – METHYL ALCOHOL. Non-compliance with this specification could cause: RISK OF WEAR ON PLASTIC PARTS RESULTING IN BREAKAGE. Do not spray any kind of chemical agent directly on the light. In particular, use is forbidden of abrasive substances, acids, containing chlorine.
	Warnings for danger of suspended loads falling Strictly comply with compliance for maximum loads planned. Do not knock against or overload the limit switches on the arms and heads.
	Warnings for biological danger and glare Do not fasten or focus the light strip directly in the patient’s eyes, especially patients at greater risk of eye injury (e.g. children with eye diseases). In this case, always use appropriate guards and precautions. The light is classified as exempt from photo-biological risk in compliance with EN 62471. However, it cannot be excluded that particularly photo-sensitive patients or those who have taken photo-sensitising medicine can develop a rash or allergic reactions to light. In this case, suspend the treatment and use very low lighting levels. The articulated arm and the joints of the light allow correct positioning of the light strip.
	Warning for danger of damaging the electrical parts Do not overcharge the arms and the joints with end of stroke knocks. Rotation of the head and arms as well as the limit switches can damage the conductor insulations.
	Warning for danger of explosion The device is not suitable for installation in environments with the presence of inflammable gas or risks of oxygen.
	Warning for danger of patient-patient cross contamination The dentist must use disposable protection on the handles of the light and sterilise them after each patient. To disinfect the surfaces use water-alcohol mixed disinfectant (see maintenance/cleaning paragraph).
	Warning for danger of wrong maintenance Do not carry out maintenance operations or replacements of parts other than those outlined in the manual. Any intervention not indicated in the manual could compromise the safe appearance of the device. Only carry out the maintenance operations in the manual; in any other case, contact technical support.

The product is covered by WEEE Directive 2012/19/EU

To scrap and dispose of the materials, comply with the standard in force in your country, if necessary contacting recognised and authorised specialist companies.

At the end of the life cycle, divide the materials based on their type (ferrous, rubber, plastic).

Do not leave small parts of the equipment unguarded or within reach of exposed people (children) because they are a potential sources of danger.

The company FARO does not allow any changes to the product not expressly authorised in writing, otherwise conformity to the safety standards and the warranty expires.

Other warnings are outlined in the titles of this manual.

1.5 STORAGE AND USE: ENVIRONMENTAL PROVISIONS


The appliance in the original packaging can be transported or kept in a warehouse for a period of 15 weeks if the following environmental conditions are met:

- Environmental temperature from -20°C to + 70°C
- Relative humidity from 10% to 90%
- Atmospheric pressure from 500 to 1060 mbar


The appliance must be used in the following environmental conditions:

- Temperature from 10° to 40°C
- Max altitude: 2000 m
- Relative humidity from 30% to 75%


ATMOSPHERIC PRESSURE

RELATIVE HUMIDITY

USE TEMPERATURE

1.6 REQUIREMENTS FOR ELECTROMAGNETIC COMPATIBILITY

This medical device requires particular precautionary measures to ensure electromagnetic compatibility, and must be installed and used in compliance with the information provided in the accompanying documentation.

Manufacturer's guide and declaration – Electromagnetic emissions		
The ALYA lamp is made to function in the electromagnetic environment specified below. The client or user should make sure that it is effectively used in this environment.		
Emission tests	Conformity	Guide – Electromagnetic environment
RF Emission CISPR15	Yes	The ALYA lamp uses RF energy only for internal functioning. Therefore, its RF emissions are very low and probably do not cause any interference with other nearby electronic equipment.
RF Emission CISPR15	Yes	The ALYA lamp is suitable for use in all buildings, including domestic buildings and those connected directly to the low voltage public electricity mains that power domestic buildings.
Harmonic emission	Class C	
RF Emission CISPR11 / EN 55011	Yes	The ALYA lamp is not suitable for interconnections with other devices (ceiling model).

ELECTROMAGNETIC IMMUNITY

Manufacturer's guide and declaration – Electromagnetic immunity		
The ALYA lamp is made to function in the electromagnetic environment specified below. The client or user should make sure that it is effectively used in this environment.		
Immunity test	Conformity	Guide – Electromagnetic environment
Electrostatic discharge (ESD) IEC/EN61000-4-2	± 6kV contact ± 8kV air	Flooring must be in wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC/EN61000-4-4	± 2kV power supply ± 1kV for input/output lines	The quality of network voltage should be equivalent to that of a typical commercial or hospital environment.
Surge IEC/EN61000-4-5	± 1kV differential mode ± 2kV common mode	The quality of network voltage should be equivalent to that of a typical commercial or hospital environment.
Voltage dips, short interruption and voltage variation IEC/EN61000-4-11	< 5% Ut for 0,5 cycle 40% Ut for 05 cycle 70% Ut for 25 cycle <5% Ut for 5 sec.	The quality of network voltage should be equivalent to that of a typical commercial or hospital environment. If the user of the ALYA lamp requires continuative use, also when electrical power is not present, it is advisable to have a UPS (uninterrupted power supply).
Power frequency magnetic field IEC/EN61000-4-8	3A/m	Level of the magnetic field at the typical power main frequency of a commercial or hospital environment.
Conducted immunity IEC/EN61000-4-6	3Vrms 150kHz to 80MHz (for non life-supporting equipment)	Portable and mobile RF communication devices should not be used near any part of the dental unit, including cables, except for when they comply with the recommended separation distances calculated by the equation applicable for the transmitter frequency. Recommended separation distances: $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P}$ from 80 MHz to 800 MHz $d = 2,3\sqrt{P}$ from 800 MHz to 2,5 GHz Where P is the maximum nominal power emitted from the transmitter in Watts (W), according to the transmitter manufacturer, and d is the recommended separation distance in metres (m). The intensity of the fixed RF transmitter field, as determined in an electromagnetic study of the site a, could be less than the level of conformity of each frequency interval. Interference could be expected in proximity to devices bea-
Conducted immunity IEC/EN61000-4-6	3Vrms 80MHz to 2.5GHz (for non life-supporting equipment)	

Note: Ut is the value of the voltage.

Note 1: At 80 MHz and 800 Mhz the higher frequency interval is applied.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is influenced by the absorption and reflection of structures, objects and people.

a) The ISM bands (industrial, scientific and medical) between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz and 40.66 MHz to 40.70 MHz.

b) The levels of conformity of the ISM bands between 150 kHz and 80 MHz and in the bands 80 MHz to 2.5 GHz are indicative. The probability of a portable transmitter causing interference if accidentally brought in the range of a patient decreases depending on environment conditions. This is why an additional factor of 10/3 has been incorporated into the formula used to calculate the separation distance from transmitters.

c) Field intensities for fixed transmitters like base stations for radio-telephones (cellular and cordless phones) and terrestrial transmitters, amateur radio stations, AM and FM radio transmitters and TV transmitters cannot be predicted theoretically and precisely. To establish an electromagnetic environment caused by fixed RF transmitters, an electromagnetic investigation of the site should be done. If the intensity of the field measured in the area of use of the dental unit is higher than the applicable level of conformity described above, it is necessary to monitor the lamp to check that it functions properly. If abnormal performance is detected, additional measures could prove necessary, such as a different orientation or position of the lamp.

d) The field intensity in an interval of frequencies from 150 kHz to 80 MHz should be less than 3 V/m.

Recommended separation distances between portable and mobile radio-communication devices and the dental unit			
The ALYA lamp is designed to function in an electromagnetic environment in which disturbances from RF emissions are kept under control. The client or operator of the unit can contribute to preventing electromagnetic interference by assuring a minimum distance between mobile and portable RF communication devices (transmitters) and the dental unit, as recommended below, in relation to the maximum output power of the radio-communication devices.			
Maximum nominal output power of the transmitter W	Separation distance at the transmitter frequency m		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters with a maximum nominal power emission not listed above, the recommended separation distance d in metres (m) can be calculated using the equation applied to the transmitter frequency, where P is the maximum nominal power emission of the transmitter in Watts (W) according to the manufacturer.			
Note: At 80 MHz and 800 MHz, the highest interval is applied. These guidelines may not apply in all situations. Electromagnetic propagation is influenced by the absorption and reflection of structures, objects and people.			

2. GENERAL FEATURES

2.1 DESCRIPTION OF THE PRODUCT

The device is used by the user to light the operative field when treating dental pathologies.

The light source on the head is composed of two LED whose light reflects against two parabolas passing through two secondary lenses.

The parabolas enable a regular and uniform spotlight to be obtained at each lighting level and to evenly distribute the light in the operative field, without creating shadows or obstruction by the operator.

Adjustment of light intensity can be carried out with a joystick or proximity sensor. The proximity sensor switches on and off the light without having direct contact, thereby eliminating the possibility of cross infection on the command.

The “automatic switch on” or “Auto-on” function allows the light to switch on automatically each time the light power supply is activated.

The remote cable allows bringing the light commands to the combined unit. Comply with the instructions provided in the installation paragraph.

On the head near the joystick and/or the proximity sensor there is a button which allows for activation of the synchronization function with the environment lamp produced by Faro.

The synchronization function allows the Alya lamp to control the level of illumination of the environment lamp in order to guarantee a more uniform level of illumination between the operating area and the surrounding zone, in order to reduce the effect of glare and to improve comfort. Visual.

In the version with light on the rear arm called “Alya with Theia Tech” the light source is composed of a series of LED whose light passes through a diffuser for distribution in the environment below.

Adjustment of the lighting level takes place in synchronous mode with that of the head for which on reducing or increasing lighting produced by the head, that on the arm also adjusts.

The “Theia Tech” light version also has a local command intended only to provide on/off on the fixed arm.

Once the local command is switched on, the light automatically synchronises with the level of head intensity. If the light on the head is switched off, the light on the rear arm switches on to maximum intensity. On switching on the light on the head, the light on the arm will synchronise automatically.


The light on the arm performs a visibility increase function in the pre-operative zone to reduce glare that generates following and/or after viewing the operative field.]

Maintenance is facilitated thanks to application of the new technologies which take into consideration the various needs for safety, ergonomics and hygiene.



The handles can be removed and sterilised. Comply with the specifications defined in the specific section.

For the electrical connections, comply with the instructions supplied in the installation paragraph and the wiring diagrams included in this manual


2.2 DESCRIPTION OF THE PARTS

- 1 – Head
- 2 – Central arm
- 3 – Rear arm without transformer with or without light (Theia Tech)
- 4 – Joints
- 5 – Pin for combined or application connection


- 1 – Joystick
- 2 – LED secondary lens
- 3 – Parabola
- 4 – Sterilisable handles
- 5 – Proximity

- 6 – Electronic board slot
- 7 – Switch on and adjustment symbols
- 8 – Synchronization button

Joystick version

Sensor version

2.3 DEVICE IDENTIFICATION

The variants on the market are differentiated by:

- Type of device (complete light, complete light with Theia Tech or head)
- Switch on and adjustment interface (**Joystick or proximity sensor, for complete lights and head**)
- Combined control mode (on-off function, remote control; for complete light and head)
- Type of mount (combined, ceiling; for complete light only)
- Arms length (for complete light only)
- Power supply (with or without transformer, for complete light only)

Codes are developed as follows:

ALYA – Complete light					
Type	3° digit Mounting and combined control	4° digit – Voltage interface	5° digit – rear arm x central arm (mm)	6° digit Available	7° 8° 9° digits Customisation
5 1	0 Standard combined	0 Joystick 17-24 V AC 22-35 V DC	0 810x550	0	000 (std FARO) JJJ
	1 Ceiling	1 Proximity sensor 17-24 V AC 22-35 V DC	1 960x550		
	2 Combined Auto-on		9 810x 855		
	4 Combined Rem cable				



ALYA – Complete light with Theia Tech					
Type	3° digit Mounting and arm dimensions	4° digit Head type	5° digit interface type	6° digit Voltage	7° 8° 9° digits Customisation
52	1 Combined 550*810	1 Standard	1 joystick	1 24Vac 50/60Hz 24 V dc	000 (std FARO) JJJ
	2 Combined 550*960	2 With extension	2 Joystick auto-on		
	3 Combined 855*810		3 sensor		
	4 Combined 855*960		4 Auto-on sensor		
	5 Ceiling 550*810		5 Joystick + remote cable		
	6 Ceiling 550*960		6 Auto-on joystick + remote cable		
	7 Ceiling 855*810		7 Sensor + remote cable		
	8 Ceiling 855*960		8 Auto-on sensor + remote cable		

ALYA HEAD					
Type	3° digit Temperature col and combined control	4° digit – Power supply and control	5° digit Available	6° digit Colours	7° 8° 9° digits Custom
3 1	5 Optical unit 5000 K	0 Joystick 17-24 V AC 22-35 V DC	0	0 White	000 Std FARO
	6 Optical unit 5000 K On/Off	1 Proximity sensor 17-24 V AC 22-35 V DC		3 Grey	
	8 L.U. 4000 K				

3. DEVICE INSTALLATION

	Warnings for electrical danger, for the light and head:
	The device must be installed by specialist technicians. On installation, the power supply must always be disconnected. Refer to the wiring diagrams in the manual. Check the plate data before installation
	Note for installation
	The power supply cable on the complete light is supplied without any connector or terminal to allow connection according to the specifications of the combination or application. The functionality and safety of the light does not depend on the polarity of the power supply current. Therefore inversion of the power supply cables will not pose a risk of malfunctioning.

3.1 DIMENSIONS

3.2 COMPLETE DENTAL LIGHT

3.2.1 Electrical requirements

The requirements for correct installation for any application (combined, wall, floor or ceiling) are as follows:

Power supply	Power cable	Type of power supply and safety requirements	Classification	Compliance with IEC 60601-1
Complete light version 17-24 Vac 50/60 Hz	2 x 0.5 mm ² 300 V 105°C PVC insulation diameter insulation 1.85 mm Only use certified terminals and connectors with resistance to flame V1 or similar.	transformer complies with IEC/EN 60601-1 third edition and IEC/EN 60601-1-2 with thermal protection or protection downstream with at least one appropriate fuse: • T1.6AL 250V Minimum requirements: • Output: 17-24 Vac; • Power: 26 VA; • Class B; • Rigidity over 4000 V. • Thermal protection	Component built-in	The medical system must be declared compliant with IEC/EN60601-1 by the installation technician or manufacturer. Note for the installation technician: ensure the combined version on which the light is installed is certified to host the complete light.
Complete light version 22-33Vdc		Power supply unit complies with IEC/EN 60601-1 third edition and IEC/EN 60601-1-2 with thermal protection or protection downstream with at least one appropriate fuse: • T630mA 250V Minimum requirements: • Output: 22-33 Vdc; • Power: 14 VA; • Class B; • Rigidity over 4000 V. Continuous protection from short circuit or overcurrent		
Complete light version with Theia Tech 24Vac 50/60Hz		transformer complies with IEC/EN 60601-1 third edition and IEC/EN 60601-1-2 with thermal protection or protection downstream with at least one appropriate fuse: • T2AL 250V Minimum requirements: • Output: 24Vac; • Power: 40VA • Class B; • Rigidity over 4000 V. Thermal protection		

Power supply	Power cable	Type of power supply and safety requirements	Classification	Compliance with IEC 60601-1
Complete light version with Theia Tech 24Vdc	2 x 0.5 mm ² 300 V 105°C PVC insulation diameter insulation 1.85 mm Only use certified terminals and connectors with resistance to flame V1 or similar.	transformer complies with IEC/EN 60601-1 third edition and IEC/EN 60601-1-2 with thermal protection or protection downstream with at least one appropriate fuse: • T2AL 250V Minimum requirements: • Output: 24 Vdc; • Power 28 VA • Class B; • Rigidity over 4000 V. Thermal protection	Component built-in	The medical system must be declared compliant with IEC/EN60601-1 by the installation technician or manufacturer. Note for the installation technician: ensure the combined version on which the light is installed is certified to host the complete light.

Tab 1 – Requirements for electrical connection and compliance with IEC 60601-1.

Check the packaging contains the following parts:

- Dental light / Head (in the required version)
- Sheet to download for site instructions www.faro.it/download

3.2.2 Safe working load

The dental light EDI Led type suits various installation positions:

DENTAL UNIT, CEILING, WALL, FLOOR.

SAFE WORKING LOAD

	Total load (SAFE WORKING LOAD)	Safety load (MINIMUM BREAKING LOAD)
Mobile arm 855 mm	29.2 N	235 N
Mobile arm 550 mm	25.6 N	205 N

3.2.3 Complete combined version light assembly

With a digital level, ensure the connection element is perfectly parallel to the ground.

Install the light by inserting the light terminal pin in the specific combined compartment.


Connect the power supply cable according to the specifications outlined in Tab. 1.

Check the light stays balanced in all positions. If necessary, use the adjustment system on the spring to balance the light.

Check switch on and adjustment and (if present) the Auto-on command and the remote cable.

3.2.4 Remote cable connection

Connect the cable to the two buttons (A and B) with contact normally open (already supplied) according to the following diagram:

3.2.5 Installation applications

- The applications are not supplied with the light
- The device has to be installed by specialized technicians only



The light has to be installed with FARO applications only.

The light is supplied with rotation limit switch between the fix and the mobile arm. THE LIMIT MUST NOT BE PASSED OVER OR FORCED.

CEILING MOUNTING INSTALLATION

1. Ceiling flange
2. Expander
3. Screw
4. Washer
5. Cable gland
6. Terminal board
7. Flange
8. Screw
9. Nut
10. Washer
11. Washer
12. Nut
13. Screw
14. Support column
15. Ring
16. Ceiling light support
17. Plug
18. Screw
19. Screw
20. Column coupling bushing
21. Lamp pin
22. Key switch
23. Fixing guide

CEILING MOUNTING

NB1. The device has to be installed by specialized technicians only

NB2. Always cut off the power supply prior to installation

NB3. Before proceeding with the assembly, ensure the ceiling can bear the application.

The authorized construction materials of the ceiling are concrete and natural stone. Use the supplied fittings or equivalent ones.


NB4. Maximum applicable load: 70 kg

NB5. Installation is allowed only in spaces provided with electrical plant compliant with the national legislation in force as regards spaces for medical use

NB6. The lamp without transformer must be powered with low-voltage (12-24Vac o 17-33Vdc) using a safety transformer (according to IEC/EN 60601-1) with thermal protection or at least one suitable fuse (T500mA/250V~).

INSTALLATION STEPS


A. After having fixed as reference point the chair centre (A), install at a distance of 650mm and 150mm, according to the directions given in the figure.

B. Unfit the flange (7) by removing the nuts (12) and washers (11).

C. By using as guide the flange (1), carry out on the ceiling 4 bores with Ø14 drill. Fit the expanders (2) in these bores.


D. Take the flange (1). Pass the power cable through the cable gland (5), then push the flange (1) against the ceiling; do not choke the cable between the flange (1) and the ceiling. Pass the screws (3) and washers (4) through the 4 bores used as reference in drilling the ceiling. Lock with the special wrench (installation accessory) the screws (3)

E. Connect the power cable to the terminal board (6) (see electrical diagrams)


CEILING MOUNTING

- F.** Match the 2 fixing guides (23) to the screws (8) and fix then by means of the nuts (9) and washers (10)

- G.** Calculate the proper length of the column (14), according to the formula $L=H-1900\text{mm}$. Cut the exceeding column (14) part on the side where NO LATERAL BORES ARE PRESENT

- H.** Insert the column (14) in the flange (7) and mark on the column (14) the position of the bores on the flange (7). Pay attention the orientation of the column in relation to the dental unit. Remove the column and carry out two $\varnothing 8$ bores at the marked points


- I.** Fit on the column (14) the ring (15) at about 300mm height (it is only temporary position to allow the assembling)

- J.** Insert the ceiling light support (16) on the column (14)


- K.** Fit the column (14) in the special bore of the column fixing flange (7)

- L.** Lock the screw (13) and the two screws (18) with hexagonal spanners (installation accessory). Tighten sturdily the screw (13) and make sure the screws (18) have passed through the bores on the column (14)


CEILING MOUNTING

- M.** Hook the unit freshly assembled (column fixing flange (7) + column (14)) to the fixing guides (23), by matching the 4 bores of the flange (7) to the screws (8) of the flanged on the ceiling (1)


- N.** Screw (without locking) the nuts (12) and the remaining washers (11) on the screws (8) of the ceiling flange (1)

- O.** Unscrew the three screws (19) of the column (14) and remove the bushing (20)


- P.** Insert the bushing (20) on the lamp pin (21)

- Q.** Insert in the pin (21) groove the key switch (22)



CEILING MOUNTING

- R. Slip - from the top - inside the column (14) a traction cord.
- S. Connect the lamp conductor to the traction cord.
- T. Fit the lamp on the column (14) and fix it with the three screws (19); make sure the bores of the bushing (20) match the screws seat on the column (14) and tighten. Simultaneously pull the traction cord to push out the conductor of the lamp from the column fixing flange (7) of about 200 mm.
- U. Connect the lamp conductor to the terminal board (6) (see electrical diagrams).
- V. Verificare la perpendicolarità della colonna agendo sui dadi (9).
- W. Stringere i dadi (12) e le rondelle (11) per fissare la flangia (7), rendendola indipendente dalle guide di fissaggio (23).
- X. Far aderire la plafoniera (16) al soffitto, spingendovi contro l'anello (15).

WIRING DIAGRAM



CEILING MOUNTING WITHOUT transformer

WALL MOUNTING

LIGHT MOUNTING WALL MOUNTING

1. Plug
2. Screw
3. Wall mounting
4. Expanders
5. Terminal board
6. Terminal board cover
7. Cover
8. Screw
8. Washer
8. Cable gland

NB1. The device has to be installed by specialized technicians only

NB2. Always cut off the power supply at the installation site prior to installation

NB3. Before proceeding with the assembly, ensure the wall can bear the application.

The authorized construction materials of the wall are concrete and natural stone. Use the supplied fittings or equivalent ones.

NB4. Maximum applicable load: 70 kg

NB5. Installation is allowed only in spaces provided with electrical plant compliant with the national legislation in force as regards spaces for medical use

NB6. The lamp without transformer must be powered with low-voltage (12-24Vac o 17-33Vdc) using a safety transformer (according to IEC/EN 60601-1) with thermal protection or at least one suitable fuse (T500mA/250V~).

The resulting medical system has to be declared by the installer as being compliant with IEC/EN 60601-1

WALL MOUNTING


fig.A



fig.B

Always cut off the power supply at the installation site prior to installation

- After having fixed as reference point the chair centre (see fig. **A-B**), carry out on the wall three d.14 bores at the bores of the wall mounting (3); pay attention the perpendicularity of the bore.
- Insert the expanders (12) in the d.14 bores previously carried out and block with the special hex spanner (installation accessory) the 2 screws; do not choke the cable between the application (3) and the wall.
- Apply three plugs (1) on the wall mounting (3) bores.
- Unscrew the screw (8) Remove the cover (7), insert the lamp arm inside the application and grease the pin. Connect the lamp cables to the terminal board (5) (see the electrical diagram below), including the earthing cable. Connect the wires which come out of the wall to the terminal board, if previously walled in. In the absence of such arrangement, the connection as to be carried out with a movable cable, to be passed in the cable gland (10).
- Fit back the cover (7) by means of screws (8).

WIRING DIAGRAM



Wall mounting

FLOOR MOUNTING

FLOOR MOUNTING

1. Pin
2. Bushing
3. Screw
4. Support column
5. plugs
6. Screw
7. Washer
8. Screw
9. Cover
10. Terminal board
11. Cable gland
12. Expanders
13. Floor support
14. Dowels
15. Plug

FLOOR MOUNTING

The device has to be installed by specialized technicians only

NB2. Always cut off the power supply prior to installation

NB3. Before proceeding with the assembly, ensure the floor can bear the application.


The authorized construction materials of the wall are concrete and natural stone. Use the supplied fitting or equivalent ones.

NB4. Maximum applicable load: 70 kg

NB5. Installation is allowed only in spaces provided with electrical plant compliant with the national legislation in force as regards spaces for medical use

NB6. The lamp without transformer must be powered with low-voltage (12-24Vac o 17-33Vdc) using a safety transformer (according to IEC/EN 60601-1) with thermal protection or at least one suitable fuse (T500mA L250V~). The resulting medical system has to be declared by the installer as being compliant with IEC/EN 60601-1


A. After having fixed as reference the centre of the chair “a”, install at a distance 650mm and 150mm according to the directions given in figure “C”

- Always cut off the power supply at the installation site prior any operation
- After having fixed as reference point (a) the chair centre (see fig.C), drill four d.14 bores in the floor, so that they match the bores on the support (13).
- Prepare the floor support (13) by passing through the washer (7) and screw (6), then screw the expanders (12) on the screws (6) a few turns, pass the power cable through the cable gland (11).
- Insert the expanders (12) in the d.14 bores previously carried out and block with the special hex spanner (installation accessory) the 6 screws; do not choke the cable between the support (13) and the floor.
- Apply four plugs (5) on the floor support (13) bores.
- Untighten the screws (8) and remove the cover plate (9) Connect the power supply cable to the terminal board (10).
- Fix the column (4) to the floor support (13); while securing it, ensure the column is perpendicular. Secure with the three screws (3) the bushing (2) to the column (4); ensure the holes of the bushing (2) match the screws on the column (4).
- Connect the light cables to the terminal board (10).
- Fit the cover plate (9) to the floor support (13) with two screws (8).

WIRING DIAGRAM

Floor mounting

3.3 HEAD


3.3.1 Mechanical requirement

For mechanical connection you must have an adequate space for the pin in the head and the G. fastening components.

The support system must be designed to support the following loads, multiplied by the safety factors outlined by IEC 60601-1 or IEC 80601-2-60

Alya Head	Alya Screen
1,80 kg	0,35 kg
2,15 kg	

For mechanical connection, follow the procedure below:

1 - Support the head and insert the washers in the threaded pin in compliance with the sequence in the figure. 2 - Then insert ring nut G in compliance with the sequence indicated in the figure and screw in with adequate equipment. The ring nut must be screwed in to give the right rotational force to the head.	
3 – Screw in the 2 safety screws F.	 UN42-0015BF
	Attention The central arm without the head load tends to rise in a sudden manner with the risk of knocking against parts of the body. During the entire installation, keep the central arm in position and do not release it until head installation is complete.
	Warnings for danger of suspended loads falling Attention - Risk of head falling after installation: - only use screws supplied by FARO. - screw in the safety screws together.

Once mechanical connection is complete, complete electrical wiring.

3.3.2 Electrical requirements

The requirements for correct installation of **the head** are as follows:

Power supply	Power cable	Type of power supply and safety requirements	Classification	Compliance with IEC 60601-1
17-24 Vac 50/60 Hz	Power cable: 2 red unipolar cables: UL Style 1061 300 V T 80°C 1x26 AWG VW 1 Ø max 1.02mm Std connector: molex 51021-0300 with 3 poles	Transformer complies with IEC/EN 60601-1 third edition and IEC/EN 60601-1-2 with thermal protection or protection downstream with at least one appropriate fuse: • T1.6AL 250V Minimum requirements: • Output: 17 - 24 Vac; • Power: 26 VA; • Class B; • Rigidity over 4000 V. • Thermal protection	Component built-in	The medical system must be declared compliant with IEC/EN60601-1 by the installation technician or manufacturer. Note for the installation technician: ensure the combined version on which the light is installed is certified to host the complete light.
22-33Vdc		Power supply unit complies with IEC/EN 60601-1 third edition and IEC/EN 60601-1-2 with thermal protection or protection downstream with at least one appropriate fuse: • T630mAL 250V Minimum requirements: • Output: 22-33 Vdc • Power: 14 VA; • Class B; • Rigidity over 4000 V. • Continuous protection from short circuit or overcurrent		

4. INSTRUCTION FOR USE

Read paragraph 1 carefully for safe use of the device.

The device must be clean before use (see Device Cleaning paragraph).

	Attention Simultaneous use of the light with electro-surgical scalpels can cause its malfunctioning.
	Attention The control joystick must be handled with care to avoid breakages. Never move the light using the joystick to grip.
	Note Each time you switch on the light, the light intensity will be saved on previous switch off.
	Warning - danger of contact with powered parts do not use the device if parts or casing are damaged.

4.1 SWITCH ON AND OFF

Refer to §1.1 for the switch on and adjustment symbols.

Complete Light

To switch on and off the light, press and release the joystick lever using the right or left side. The light intensity on switch on will always be the last used before switch off.

Complete light with Theia Tech

Same operations of the complete light plus those on the fixed arm switch on and/or off in sync with that of the head.

The light on the fixed arm can be switched on/off even using the button on the arm. If the light is switch on, the light will be automatically synchronised, if the light on the head is off, it will be regulated to maximum intensity.

4.1.1 Adjustment:

a) to reduce light intensity keep the lever pressed on the joystick on the left side (light rear view) until the desired intensity is reached.

On reaching the minimum intensity, you will hear an acoustic signal (1 beep).

a) to increase light intensity keep the lever pressed on the joystick on the right side (light rear view) until the desired intensity is reached.

On reaching the maximum intensity, you will hear an acoustic signal (1 beep).

c) to skip to maximum intensity, press and release the joystick lever using the front or rear side. On subsequent pressing on the front or rear side, the light intensity will return to that previously saved.

The light on the fixed arm is adjusted in sync with that of the head, it cannot be adjusted independently.

4.1.2 Light / Head WITH PROXIMITY

Switch on/off

Complete Light: To switch on or off, approach the sensor once up to a maximum distance of 3 cm.

Complete light with Theia Tech: same operations of the complete light plus the light on the fixed arm switch on and/or off in sync with that of the head.

Adjustment

Complete light: To adjust light intensity, you need to stay still near the sensor until the desired intensity is obtained. Adjustment enables passage from the maximum value to the minimum value again to the maximum. On reaching the maximum intensity, you will hear an acoustic signal (2 beeps). On reaching the minimum intensity, you will hear an acoustic signal (1 beep).

	Note Each time you switch on the light, the light intensity will be saved on previous switch off.
---	--

Complete light with Theia Tech same operations of the complete light plus the light on the fixed arm will adjust in sync with that of the head.

4.1.3 Light / Complete light with Theia Tech / “ALYA” head WITH REMOTE CONTROL

Switch on / Switch off / Adjustment

- To switch on and/off press and release button “A”.

- Adjustment:

a) to reduce the light intensity keep button “A” pressed until reaching the desired intensity.

On reaching the minimum intensity, you will hear an acoustic signal (1 beep).

b) to increase the light intensity keep button “A” pressed until reaching the desired intensity.

On reaching the maximum intensity, you will hear an acoustic signal (2 beeps).

c) to immediately reach the minimum light intensity, press button “B”.

On reaching the minimum intensity, you will hear an acoustic signal (1 beep).

Subsequent button pressing will bring the light back to the previous light intensity selected.

	Note Each time you switch on the light, the light intensity will be saved on previous switch off.
---	--

4.1.4 LAMP/ COMPLETE LAMP WITH THEIA TECH / “ALYA” HEAD WITH SYNCHRONIZATION COMMAND

Where foreseen, it is possible to connect the Alya lamp to the Faro environment lamp in wireless mode in order to turn them into a synchronized lighting system named “Synchro”.

The “Synchro” mode has been purposefully designed to improve the comfort of the medical dentist/orthodontist, in order to reduce the effect of glare generated when passing from observation of a highly illuminated surface (e.g. oral cavity with a dental lamp) to a poorly lit surface (e.g. dental quiver).

With the method named “Synchro”, activatable using the button situated at the head of the Alya lamp, it is possible to change the illumination produced by the Faro environment lamp automatically on the basis of the value of illumination produced by Alya.

Note: Between the dental lamp and the environmental one, a small delay in synchronization may be experienced. This is due to the communications protocol. This effect is normal and does not signify a defect.

The “Synchro” function, in order to be enabled, requires a “Pairing” process (which should only be carried out once) in order to create a link between the two lamps. Subsequently, the “Synchro” function can be enabled and/or disabled as required by the user, from the button situated on the dental lamp.

“PAIRING” PROCEDURE

NOTE:

- The “Pairing” procedure is only necessary at the first connection, although it can be repeated in the event of changing the Alya lamp or the electronics of one of the two inter-connected lamps.

- If there should be more environment lamps in the surgery, ensure that the other lamps are turned off or have been on for more than 60 seconds.

To carry out the **“Pairing”**, proceed as follows:

1. Apply current to the Faro environment lamp you wish to couple.
The environment lamp prepares itself for a Pairing connection, for a maximum time of 60 seconds.
2. Within the 60 seconds, press the “synchro” button on the dental lamp for at least 3 seconds, but for no more than 6 seconds, otherwise the procedure will be cancelled.
On receipt of the **“Pairing”** request by the dental lamp, the environment lamp activates the blue LED on the aluminium body. If the blue LED should not activate, it is possible (within 60 seconds from switching on) to carry out other attempts. Once the time has been exceeded, it is necessary to repeat the procedure from point 1.
3. From the blue LED illuminating on the environment lamp, there are 60 seconds to confirm the **“Pairing”** by pressing the programming button  situated on the remote control of the environment lamp. At this point, the blue LED (of the environment lamp) flashes at double speed and then goes out. If the button  on the remote controls is not pressed within 60 seconds, the blue LED goes out and the procedure must be repeated from point 1.

Having carried out the **“Pairing”**, synchronization between the 2 lamps is now enabled.

To **ACTIVATE THE SYNCHRONIZATION FUNCTION**, proceed as follows:
Press the Synchro button for 2 seconds then release it.
On release, a sound signal (Beep) will be heard and the blue LED on the environment lamp will light to indicate that synchronization has been activated.


To **DE-ACTIVATE THE SYNCHRONIZATION FUNCTION**, proceed as follows:
Press the Synchro button for 2 seconds then release it.
On release, a sound signal (Beep) will be heard and the blue LED on the environment lamp will go out to indicate that synchronization has been disabled.

Note connected with synchronization:

- When the Faro environment lamp is synchronized with the Alya lamp (or it adjusts itself automatically), the blue LED on the body is lit and steady. When the LED is off, the synchronization is not active.
- The remote control is always enabled, so it is therefore possible to change the illumination value. Nevertheless, if the environment lamp should find itself in a synchronized state (with the blue LED lit), as soon as a new adjustment is carried out on the Alya lamp, the illumination value will immediately be updated.
- If the Alya lamp should be turned off, the environment lamp will remain lit at the illumination value in use.
- If the Alya lamp should be turned on, the environment lamp will turn on automatically.

4.2 JOYSTICK ALYA JOYSTICK LEVER ASSEMBLY

- Fully insert lever "A" on the pin of the joystick.
- The hole in lever "A" must be positioned corresponding to surface "B".
- Fully screw in grub screw "C" with the wrench supplied.

Attention

The control joystick must be handled with care to avoid breakages. Never move the light using the joystick to grip.

5. ORDINARY MAINTENANCE

There are no ordinary maintenance operations.

6. CLEANING

Warning against danger of wear and corrosion and falling suspended weights

For all metal or plastic parts of the light, it is strictly forbidden to use substances that are abrasive, acid, substances containing chlorine or chloride ions, detergents with Trilene base, petrol, white spirit or similar.

It is forbidden to directly spray any chemical substance on the device.

6.1 PARABOLA CLEANING

Cleaning must be carried out using a soft cloth in cotton or absorbent cotton with ethyl alcohol or the specific PERFLEX detergent. Water-alcohol based disinfectants are suitable with 70% isopropyl alcohol or ethanol.

Attention - potential damage or wear on the parabolas

Never spray detergent directly on the parabolas.

Cleaning operations on the parabolas must be carried out wearing gloves, to avoid leaving fingerprints on the surfaces.

Never use detergents containing surfactants or water-repellents that depositing can leave streaks. Slight streaking will not prejudice the quality of the light.

Products differing from those suggested could damage the parabolas.

If in doubt, contact FARO customer care.

6.2 CLEANING THE HEAD

Cleaning must be carried out using a soft cloth in cotton or absorbent cotton with ethyl alcohol or the specific PERFLEX detergent. Water-alcohol based disinfectants are suitable with 70% isopropyl alcohol or ethanol.

 	Warning against danger of wear of the plastic and falling suspended weights
	Never spray detergent directly on the head. Do not use detergents-disinfectants containing the following to clean plastic parts: • AMMONIUM HYDROXIDE • SODIUM HYDROXIDE • METHYLENE CLORIDE • ME THYL ALCOHOL. • ALL KINDS OF ACIDS Faro has tested and suggests the following disinfectants: Eco Jet-1 (Cattani Group) / Sporekin Plus DS (Ims srl) / Zefirol Quick (Molteni Dental) / Durr FD366 Sensitive

6.3 ARM CLEANING

Always use a cloth soaked in disinfectant approved to disinfect the surfaces and pass it over.


 	Warning against danger of corrosion and mechanical collapse with falling suspended weights
	Never spray chemical substances directly on the arms or joints and their openings.
 	Warning against danger of wear of the plastic with falling suspended weights
	Do not use detergents-disinfectants containing the following to clean plastic parts: • AMMONIUM HYDROXIDE • SODIUM HYDROXIDE • METHYLENE CLORIDE • ME THYL ALCOHOL. • ALL KINDS OF ACIDS Faro has tested and suggests the following disinfectants: Eco Jet-1 (Cattani Group) / Sporekin Plus DS (Ims srl) / Zefirol Quick (Molteni Dental) / Durr FD366 Sensitive

7. STERILISING THE HANDLES

	Warning - danger of cross contamination
	The handles are not supplied sterile, they must therefore be sterilised before use. The handles must be sterilised before each patient.

7.1 Removing the handles

To remove the handle, unscrew knob "A" and remove it from the support.

7.2 Decontamination and disinfectant

Before sterilising the handles, they must be decontaminated and disinfected.

To disinfect, Faro has tested the following products for disinfection: Eco Jet-1 (Cattani Group) / Sporekin Plus DS (Ims srl) / Zefirol Quick (Molteni Dental) / Durr FD366 Sensitive

 	Attention - danger of plastic breaking
	The handles cannot be disinfected by thermo-disinfection.

7.3 Sterilisation

The handles must be packaged in compliance with EN 868-5.

The handles can be sterilised with standard cycles 121°/134° C up to two hundred (200) cycles or however up to loss of the mechanical performance.

The parameters of the sterilisation cycle are as follows:

Ciclo EN 13060	Temperature	Pressure	Minimum Holding Time
B	121°C	207 KPa	15 min
B	134°C	308 KPa	3 min

8. PERIODIC CHECKS

Operation	Frequency	Applicability		Procedure	Enabling
		LD	TE		
Check the absence of range between the arm joints	Annual	x	N/A	Check the gap between the joints 5 and the arm has not changed since first use	User
Check data plate is legible	Annual	x	x		User
Check casing is intact	Twice-yearly				Service Engineer
Check electrical safety EN 62353 1. Rigidity 2. Dispersion	Twice-yearly	x	x	Measure dielectric rigidity and dispersion of the casing. Limits defined in IEC 60601-1	Service Engineer
Light parameter checks.	Twice-yearly	x	x	With a spectroradiometer, measure the values of: • Maximum lighting: >35000 lux • CRI decay: <20%. • Underlying value of the blue light on the spectrum emitted measured in: <100 W/m2	Service Engineer

Service Engineer: competent person in electro-mechanical equipment maintenance.

9. ACOUSTIC ALARM SIGNALS

9.1 Acoustic Signals

OpL** = Beep 30 seconds

OTP* = Beep 30 seconds

* OTP: LED overtemperature protection.

** OpL: Disconnected led load

9.2 GUIDE TO PROBLEMS

The table below represents a guide to the potential defects of the light.

If you cannot resolve the problems, call technical support

Effect	Cause	Action (Service Engineer - SE)	Resp
The light does not switch on	Power supply not inserted or inserted incorrectly.	Check the power supply is inserted and the combined unit is on.	User
	Interference with electro-surgical scalpels or high energy tools	Switch off the electro-surgical scalpels and check the permanence of the effect.	User
	Command on joystick applied incorrectly.	To switch on and off the light, press and release the joystick lever using the right or left side.	User
The light will flicker	Interference with electro-surgical scalpels or high energy tools.	Switch off the electro-surgical scalpels and check the permanence of the effect.	User
The light does not adjust light intensity	Command on joystick applied incorrectly.	Use the command correctly as described in this manual.	User
	Interference with electro-surgical scalpels or high energy tools.	Switch off the electro-surgical scalpels and check the permanence of the effect.	User
The light intensity is considerably reduced	Parabolas or secondary lens dirty.	Clean the parabolas and the secondary lenses.	User
	Use of wrong procedures.	Check you have maximum adjustment with the command.	User
Stains have appeared on the parabolas or the reflective layer has come away.	Use of non-approved products.	Clean the surfaces with specific "Faro Perflex" product. Clean the surfaces with isopropyl alcohol. To restore the surfaces, you need to get service to replace the parabolas.	User
The light does not stay balanced and tends to lower	Excessive load on head (small mirrors, cameras, etc..).	Remove the excess loads.	User
The lamp does not control	Synchro function off	Activate the function; see 4.1.4	User

10. TECHNICAL FEATURES

Complete light:

Power supply voltage (without transformer):

- $17 \div 24 \text{Vac} \pm 10\%$ - 50/ 60Hz;
- $22 \div 35 \text{Vdc} \pm 10\%$

Absorbed power:

- 26VA (version $17 \div 24 \text{Vac}$);
- 14VA (version $22 \div 35 \text{Vdc}$)

Recommended fuses:

- Version $17 \div 24 \text{Vac}$: T1.6AL 250V
- Version $22 \div 35 \text{Vdc}$: T630mAL 250V

Protection against electrical danger:

- Class II device

Classification EN 62471:

- Exempt class

Complete light with Theia Tech

Power supply voltage (without transformer):

- $24 \text{Vac} \pm 10\%$ - 50/ 60 Hz
- $24 \text{Vdc} \pm 10\%$

Absorbed power:

- 40VA (version 24Vac)
- 28VA (version 24Vdc)

Recommended fuses:

- T2AL 250V

Protection against electrical danger:

- Class II device

Classification EN 62471:

- Exempt class

Optical features of the light produced by the head in compliance with ISO 9680

Spot light dimensions: 180 mm x 90 mm

Lux: 3.000*-50.000* lux @700mm

Colour temperature: 5000 K*


CRI (colour rendering index): >95*

* Typical values subject to tolerances.


Labelling in compliance with EN 62471: not necessary

10.1 WIRING DIAGRAMS


- **Complete light:**
Electrical circuit in light head

Wiring Diagram – Alya with transformer


- Complete light with Theia Tech:
Electrical circuit in light head

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
GUARANTEE CERTIFICATE

FARO offers the final customer a **24 month** guarantee starting from the date of purchase.

Repairs under guarantee must be performed at FARO; expenses and transport risks are at the risk of the purchaser. **Repair under guarantee is considered valid only when:**

- **all sections of the certificate have been filled in and sent in advance to FARO by Fax (039.6010540).**

The guarantee covers faults due to the bad quality of the material or manufacturing defects; in the case of valid claims, the guarantee covers free repair or replacement. **Claims for damages and/or interest are excluded.** The guarantee is not considered valid, at the sole discretion of FARO, if the fault is due to tampering, damage, incorrect use, improper maintenance and normal wear and tear.

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FARO S.p.A.
FARO SpA si riserva il diritto di modificare, senza preavviso, le caratteristiche indicate nel presente manuale.
FARO SpA reserves the right to change the specifications of this equipment without notice.
FARO SpA se reserve le droit de modifier, sans préavis, les caractéristiques dans ce manuel.
FARO SpA behält sich recht vor, jederzeit stillschweigend technische oder bauliche Änderung vorzunehmen.
FARO SpA se reserva el derecho de modificar sin aviso previo la características incluidas en el presente manual de uso.