Benutzerhandbuch *User Manual*

MELAcontrol® Pro

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Durable test body system in accordance with EN 867-5 and EN ISO 11140-1

Intended use


The MELAcontrol Pro is a durable test body system in accordance with EN 867-5, consisting of a test body and an indicator strip. In accordance with EN ISO 11140-1, MELAcontrol Pro is to be classified as a type 2 indicator. MELAcontrol Pro can be used as a re-usable test body system in accordance with EN 13060 for products with narrow lumen and in accordance with EN 285 for hollow loads.

MELAcontrol Pro supplements the physical measurements and is only suitable for testing the following steam sterilizers:


- Steam sterilizers with type "B" cycles in accordance with EN 13060
- Large steam sterilizers in accordance with EN 285

Application

1. Unscrew the test body sealing cap (pos. a).

2. Insert an indicator strip (pos. b) in the slot of the plastic bracket (pos. c) of the sealing cap.

3. Return the sealing cap to the test body and screw it hand tight.

CAUTION

Danger of burns from hot condensate!

Position the test body for sterilization in a level position, otherwise hot condensate can run out of the test body after removal.

- Let the test body cool before removal.

PLEASE NOTE

The test body must be dry before re-use; if possible blow it out with compressed air. Residue condensate in the test body has a negative effect on air removal.

- Insert an unused indicator strip in the plastic bracket of the sealing cap before every test.

Use for daily routine check

1. Place a new indicator strip in the test body.
2. Slide the test body with the flat side facing downwards onto a tray in the empty sterilization chamber.
3. Start the steam sterilizer in a validated sterilization program according to the intended use.
4. Remove the test body from the steam sterilizer when the program run has ended.
5. Remove the indicator strip from the sealing cap. Evaluate the indicator strip.

Use as batch control

1. Place a new indicator strip in the test body.
2. Insert the test body with the flat side facing downwards in the sterilization chamber, together with the sterilization material.

3. Start the steam sterilizer in a validated sterilization program according to the intended use.
4. Remove the test body from the steam sterilizer when the program run has ended.
5. Remove the indicator strip from the sealing cap. Evaluate the indicator strip.

Evaluation

	Completely black discolouration	Sterilization successful sufficient air removal and steam penetration
	Red-brown discolouration	Sterilization not successful only temperature impact, no steam penetration
	Partially black discolouration	Sterilization not successful insufficient holding time and steam penetration

Documentation

The evaluation of the results must be performed by a sufficiently trained person holding responsibility for this task. The colour change is caused by an irreversible chemical reaction. If stored correctly, the indicator strips can be archived for over five years. EN ISO 15882 only requires a description of the results. This can be performed in written or digital form. Comply with the valid national and regional regulations.

Servicing

Change the seal (pos. d) in the sealing cap after every 250 tests. The seal is included in every refill pack (art. no. ME01076).

To change the seal, unscrew the white plastic bracket from the sealing cap. Using tweezers, remove the black seal from the sealing cap and insert a new seal. Screw the white plastic bracket into the sealing cap hand-tight.

Storage

Do not store the indicator strips together with chemical substances and gases (e.g. ethylene oxide [EO]), hydrogen peroxide [H_2O_2] or instruments sterilized with H_2O_2 . Incorrect storage of the indicator strips can cause unwanted colour changes. Further information regarding storage is provided on the packaging. A secure determination of malfunctions in the process can only be guaranteed up to the specified date of expiry.

Important information regarding routine operation

Please comply with the current recommendations issued by the Robert Koch Institute (RKI) and the information contained in DIN 58946-7. Further information regarding routine operation is to be found in user manual of the steam sterilizer.