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

Universal-Program

Universal-Program

Holding times during sterilization:

Universal-Program: 5:00 min (3:30min with FW 2.0)

Quick-Program B: 3:30 min

Quick-Program S: 3:30 min

Gentle-Program: 20 min (16 min with FW 2.0)

Prion-Program: 20 min

<pre> ----- !0 01100ED0E0 !1 A03350MA.PRO ----- 10 MELAG Vacuklav 40-B ----- 15 Program: Universal-Program 20 Program type: 134 °C wrapped 25 Date: 13.09.2010 30 Daily batch: 01 Whole: 00802 35 User: 1002 36 Indicators changed: Yes 37 Batch released: Yes ===== 40 Universal-Program successfully ended 42 = = ===== 45 Temperature: 135.5 +0.06/-0.14 °C 50 Pressure: 2.17 +0.00/-0.01 bar 55 Plateau time: 05 min 30 s 60 Conductivity: 6 µS/cm (698:21737.0) 65 Start time: 08:28:57 70 End time: 09:00:17 (31:20 min) ===== 80 SN:201040-B4001 ===== 81 RC V3.022 07.10.2010 82 Para V3.030 06.10.2010 83 BO V3.142 07.10.2010 ----- Step Time t[m:s] P[mbar] T[°C] Step Time t[m:s] P[mbar] T[°C] SP-S 0:01 0:01 1014 49.4 SK11 0:13 0:12 1631 75.2 . . SF22 4:34 0:39 180 104.6 SF23 5:16 0:42 1814 114.8 SF31 5:29 0:13 1295 111.6 SF32 6:10 0:41 200 106.1 SF33 6:55 0:45 1938 117.1 SF41 7:10 0:15 1294 113.0 SF42 7:36 0:26 400 108.3 SF43 8:00 0:24 1710 112.8 SH01 8:43 0:43 2695 129.8 SH02 9:04 0:21 2845 131.9 SS01 9:37 0:33 3061 134.1 SS02 15:07 5:30 3171 135.6 SA00 15:37 0:30 1289 115.7 ST01 19:40 4:03 68 83.4 ST02 21:40 2:00 698 99.0 ST03 25:40 4:00 64 87.2 ST04 26:40 1:00 657 94.1 ST05 30:40 4:00 63 83.3 SB10 30:55 0:15 800 87.0 SB20 31:15 0:20 1014 89.7 SP-E 31:15 0:00 1014 89.7 >> Never change code on follow. line << 6C000006009C028A43A2050A0FEF03273194348C >> Authentication of batch log << ----- 0.00 0.0 0.0 0.0 0.0 0.0 -edk---etm---etd---etp---etv---ett-END- </pre>	Header !0 Identification number !1 File name 10 Type of the autoclave 15 Program name 20 Sterilization parameters of the program 25 Date 30 Daily and overall batch number 35 User-ID 36 Batch indication 37 Batch release } Release process 40 Control message 42 Warning or error message for program interruption 45 Sterilization temperature with max. deviations 50 Sterilization pressure with max. deviations 55 Sterilization time 60 →Conductivity of the →feed water 65 Time at start of the program 70 Time at finish of the program 80 Serial number of the autoclave 81 Current version of the device firmware 82 Current version of the device parameter 83 Current version of the user interface Values in the program steps and legend Time Time (minutes:seconds) elapsed since the start of the program t Duration (minutes:seconds) of a program step [m:s] P Pressure in the chamber in millibar [mbar] T[°C] Temperature in the chamber in degrees centigrade At the beginning of the individual lines are initials which indicate the type of the respective program step. You will receive a listing of all step initials when you output a legend log report with the Logs output menu. Programmschritte: SK Conditioning SF Fractionation SH Hold SS Sterilization SA Pressure release ST Drying SB Ventilation SP-E End Authenticity proof (electronic signature) may never be changed; deciphering of the code by MELAG permits inference whether the data was created and changed on an autoclave from MELAG. Sensor measured values are shown here in case of a malfunction. The values are helpful for the technician.
--	--

```

-----
!0 01100ED0E0
!1 A03361GS.STR
-----
10 MELAG Vacuklav 40-B
-----
15 Program:      Universal-Program
20 Program type: 134 °C wrapped
25 Date:         12.05.2011
30 Daily batch: 06      Total: 01900
35 User:         1233
36 Indicators changed: deactivated
37 Batch released:      deactivated
=====
40 Universal-Program NOT successful
42 E112=FRM4: Maximum permissible
temperature exceeded in double jacket.
=====
45 Temperature:  --  --/-- °C
50 Pressure:     --  --/-- bar
55 Plateau time: -- min -- s
60 Conductivity: 9 µS/cm (52:7450.0)
65 Start time:   10:55:04
70 End time:     11:00:14 (05:10 min)
=====
80 SN:201040-B4002
=====
81 RC   V3.030 04.05.2011
82 Para V3.034 01.12.2010
83 BO   V3.150 04.05.2011
-----
Step  Time t[m:s] P[mbar] T[°C]
Step  Time t[m:s] P[mbar] T[°C]
SP-S  0:01   0:01   1002   98.3
SK11  0:14   0:13   1644   104.8
SK12  0:40   0:26   1292   110.0
SK11  0:48   0:08   1633   109.6
SK12  1:14   0:26   1293   110.5
SK21  1:22   0:08   1636   110.1
SU00  1:45   0:23   1515   113.7
SU00  1:46   0:01   1509   113.6
SUDP  2:28   0:42   2347   124.9
SU-A  3:43   1:15   1316   113.7
SU-V  5:00   1:17   596    111.4
SU-B  5:10   0:10   963    110.1

>> Never change code on follow. line <<
62000006DEE2FF8B56A3090A0F010371FFFE00FC
>> Authentication of batch log <<
-----
0.00 147.0 113.0 113.0   0.0 111.0
-edk----etm---etd---etp---etv---ett--END-

```

10 MELAG Vacuklav 40 B+

Legend for step token in the log file

Sterilization programs

S... Step
.K... Conditioning
SK11 Conditioning Type1 (steam entry)
SK12 Conditioning Type1 (pressure release)
SK21 Conditioning Type2 (steam entry)
SK22 Conditioning Type2 (pressure release)
.F.. Fractionation
SF12 1st Fractionation Evacuate
SF13 1st Fractionation Steam entry
SF21 2nd Fractionation Pressure release
SF22 2nd Fractionation Evacuate
SF23 2nd Fractionation Steam entry
SF31 3rd Fractionation Pressure release
SF32 3rd Fractionation Evacuate
:
SH01 Maintain steam entry
SH02 Maintain Control
SS01 Sterilization Start
SS02 Sterilization
SA00 Pressure release
.T.. Drying
ST01 Drying 1st Vacuum drying
ST02 Drying Ventilation drying
ST03 Drying 2nd Vacuum drying
ST04 Drying Steam drying
ST05 Drying 3rd Vacuum drying
SB10 Ventilation
Vacuum test
.X.. Extra programs
SX1V Evacuation phase
SX1W Equilibration phase: Temperature
before measurement
SX1M Measuring phase, measuring pressure increase
during 10min monitoring time
SX1E End of Vacuum test
.U.. Interruption, abort
SUDP Abort pressure equilibration
SU-A Pressure release after abort
SU-V Evacuate after abort
SU-B Ventilate after abort
-----END--

	Value	Sensor for	Input	in	comment
ett	theoretical temperature	(Press.)	(AIN07)	°C	computed from pressure
etv	temperature of chamber	--	AIN01	°C	for process control 10bit
etp	temperature of chamber	Temp.	AIN01	°C	for process control 14bit
etd	temperature of chamber	Temp.	AIN02	°C	for protocol 14bit
etm	temperature of double jacket	Temp.	AIN03	°C	monitoring of the double jacket 14bit
edk	pressure of the chamber	Press.	AIN07	bar	

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
User administration (setting, editing and deleting the admin / user PIN) -----

Settings made in the **User administration** can only be performed after entering the admin PIN (set as "1000" upon delivery). The pre-set admin PIN should be changed during the autoclave first installation.

1. Select the menu **Settings** (key tool) → **User administration**.
2. Enter the admin PIN and confirm with **LOGIN**. The display switches to the menu **User administration**.
3. Select the pushbutton "Continue to user list" to display the complete user list.
4. The ID 1001 is reserved for the admin PIN.
Change the pre-set admin PIN by pressing the button **EDIT** and confirm with **SAVE**.
5. All further users are added after ID 1002.
The user PINs can be added / changed in the same way as the admin PIN.

Electronic release process on the display-----

Please ensure that the release process is activated on the display before the autoclave is delivered. If the release process is not required, you can deactivate it in the following manner:


1. Select the menu **Settings** (key tool) → **User administration**.
2. Working in the window in which the options for graphic recording, the authorization options and immediate issue are determined, de-select the options "batch indication" and "batch release."

PLEASE NOTE

The user administration menu can be used to decide whether it is necessary to enter a user PIN following the end of a program run.


The **Logging** menu is used to determine all further settings regarding logging and log issue, e.g. issue media, IP addresses, transmission via FTP/TCP etc.

3. Work through the logging menu (using **CONTINUE**) and accept the settings in the last window with **SAVE**.

Entering contact addresses in the menu "Settings" -----

After the autoclave has been set up, the technician contact data of the responsible stockist / MELAG customer services can be entered in the menu **Settings** → **Service**.

Please ensure that you read the user manual before commissioning.

Please fill in:

Operator (name, address, stamp, date)
Signature operator

Few organisational measures must be taken for the documentation and the release process with the Premium-Plus-Class. We recommend to make some copies of the form, to fill in and to display it publicly.

1. We use the following MELAG autoclave

☐ Vacuklav® 40 B+
 ☐ Vacuklav® 41 B+
 ☐ Vacuklav® 43 B+
 ☐ Vacuklav® 44 B+

with the serial number _____

2. Persons who are entitled to release

The user PIN must not be documented publicly. The administrator can change the user PIN or add a new user if necessary.

User	Surname, name	Signature
1001 (Admin)		
1002		
1003		
1004		
1005		
1006		
1007		
1008		

3. Storage life

The sterilized items can stored in our practice up to ____ months after preparation and can still be used if the following requirements are fulfilled:

- Stored dust-free / free of contaminations in a cabinet or drawer, or in a separate room for sterile material.
- Protected from moisture during the entire storage period
- Protected from damage

Upgrading from a b/w display to a Colour Touch-Display (Art.-No. 67205)

Item	Article
	The conversion set (Art.-No. 67205) contains:
1	1x upper display cover (complete with screws and cover clips)
2	1x Colour Touch-Display
3	2x CF cards
4	1x CF card reader
5	1x type plate with new serial number
6	1x label with device logo
7	1x CD-Rom with documentation

Tools required: TORX spanner TX20, open-end wrench, 1 side cutting pliers

CAUTION!

Before opening the device switch off at mains and remove the plug!

Performing the electronics update-----

NOTICE

Before starting the update make a backup of all program and malfunction logs, including the status log and the system log on a CF card or a computer. Note the network settings, logging settings as well as the maintenance and batch counter values.

1. Open the autoclave door and turn off the autoclave at the mains.
2. Place the "Mainboard Update" CF card (CF card ①) loosely in the card slot of the control electronics **but not to its fullest extent!**
3. Connect the autoclave to the mains and switch on at mains.
4. As soon as the red LED next to the card slot illuminates, push the CF card (CF card ①) **immediately** in the card slot to its fullest extent. The LED illuminates as long as the update process is active (approx. 10 seconds).
5. As soon as the LED flashes, remove the CF card from the internal card slot.
6. Switch off autoclave at the mains and remove the plug.


Fig. 1

Dismantle the b/w display -----

1. Remove the Torx screws in the unit cover using the Torx TX20 spanner. Remove the unit cover.

NOTICE

To obtain free access to the rear screws of the Vacuklav® 41-B, 43-B it is necessary first to remove the rear panel and to lift the storage tank a little at the front.

2. Remove the connection cable to the b/w display (Fig. 1, pos. a) from the control electronics and slide the cable through the conduit in the holding plate of the control electronics (Fig. 1, pos. b).
3. Push the snap tape (Fig. 2, pos. c) a little outwards and fold the b/w display forwards.
4. Push the hinge guide (Fig. 2, pos. d) outwards, remove the display from the plastic bracket and remove the small pins. Guide the display cable forwards through the display cover and remove the b/w display and the cable.


Fig. 2

5. To be able to dismantle the display cover the small hinge hook (Fig. 3, detail) must be removed.
6. Remove the 6 TORX screws from the rear side (above) in order to be able to remove the displays cover.


Fig. 3

Mounting the new display cover -----

NOTICE

Open the cavity for the display fan (see Fig. 4, detail).

1. Click the display cover in the front plate and screw tight with the 5 TORX screws (one screw is not needed).
2. Re-hang the hinge hook.
3. Guide the display cable from the front through the elongated hole in the front panel (see Fig. 3, pos. e) and through the conduit in the holding plate of the control electronics (see Fig. 1, pos. b). The cable laying through the ferrit core is not applicable.

NOTICE

The storage tank must be raised forward a little in order to pull the Vacuklav® 41-B, 43-B (following the update 41 B+, 43 B+) display cable through the conduit from inside.

4. Connect the display cable connector to the connection (Fig. 1, pos. a) on the control electronics.
7. Connect the RJ45 connector in the display cover to the rear of the Colour Touch-Display (see Fig. 4).
8. Push the Colour Touch-Display from above into the rear channel of the display guide until it clicks.
9. Attach the cover clips on both sides of the display cover.
10. Re-mount the housing.



Fig. 4

Performing the software update-----

11. Slide the "**Display Update**" CF card (CF card ②) in the lateral card (see Fig. 5) slot of the Colour Touch-Display to its fullest extent! NOTE: Check the device type on the CF card!
12. Connect the autoclave to the mains and switch on at mains.
The second update now begins from the "Display Update" CF card. The Colour Touch-Display shows a progress bar. The display switches on and off many times during the update procedure.

WARNING

Never switch the autoclave off during the update procedure.

As soon as the user menu with the preinstalled language English is displayed on the screen permanently the update process is completed (see Fig. 6).

5. Remove the CF card ② from the lateral card slot of the Colour Touch-Display.

NOTICE

If a malfunction message is shown on the display after the update switch off the autoclave at the mains and then on again. The system should restart successfully.

Fig. 5


Fig. 6

Commissioning -----

1. Glue the new name plate on the rear of the autoclave next to the previous plate and place the label with the serial number next to the present plate on the left-hand side of interior next to the chamber opening.
6. Affix the label with the device logo in the centre underneath the door cover (see Fig. 7).
13. The following controls/settings must be performed in accordance with the current status log before regular operation can recommence.
 - Working in the **SERVICE menu** (castle symbol) enter the new serial number, the previously set values of the maintenance counter, the total batch counter and the number of successful batches and daily batches. (Access menu, see document "Resetting the maintenance counter").
 - Working in the **SETUP menu** (spanner symbol) set the language by entering the activation code included in the scope of delivery (see conversion and training protocol).
 - When using a water treatment unit with the "stand alone devices" the water supply must be set to "EXTERNAL" in the **SETUP menu**.
 - Working in the **SETUP menu** under the item "Logging" re-enter the IP address and FTP/TCP settings if necessary.
 - If the other parameters were changed, also re-set these.
 - Check whether the correct settings (type, serial number, logging etc) are displayed in the "**Info & Status**" menu.

Fig. 7

As documentation of proper installation and setting up as well as for your warranty claim, the conversion and instruction protocol must be filled out by the person responsible and a copy sent to MELAG.

Conversion and instruction protocol

Please send to MELAG

MELAG Medical Technology
Geneststraße 6 – 10
D10829 Berlin

Dear Madam/Sir,

within the scope of Quality Assurance we are obligated, in cooperation with you, the responsible specialist MELAG dealer, to upgrade this MELAG autoclave at the operator's in accordance with good engineering practice and provide instructions in its use.

Please copy, fill out and sign this form and send it to us after the autoclave has been successfully upgraded.

Device and contact data

We, hereinafter referred to as specialist MELAG dealer, have today upgraded the autoclave of the Premium-Class to the Premium-Plus-Class, instructed the personnel and performed the start-up as described below:

Specialist dealer (name, address/stamp) 	Operator (name, address/stamp)
E-Mail address for software update 	E-Mail address for software update
New autoclave serial number 	Date
For the 1st language.....the activation code is:..... For the 2nd language.....the activation code is:..... For the 3rd language.....the activation code is:.....	

The following persons were present during the instructions for use:

	Name, first name (please print in block letters)	Signature
From practice/clinic:		
		
		
		

	Name, first name (please print in block letters)	Signature
Performing technician:		
From the company:		

Executed work

Please tick the following points after completing the respective tasks.

Upgrading from Premium-Class to Premium-Plus-Class

☐ Upgrade performed according to separate instructions (AS_077-10_GB_40B_41B_43B_44B.pdf).

The device type has changed as follows:

	Device type before	Device type after
☐	Vacuklav [®] 40-B	Vacuklav [®] 40 B+
☐	Vacuklav [®] 41-B	Vacuklav [®] 41 B+
☐	Vacuklav [®] 43-B	Vacuklav [®] 43 B+
☐	Vacuklav [®] 44-B	Vacuklav [®] 44 B+

Check time and date setting. If necessary, set according to user manual

☐ Time and date correct/set correctly

The following address (preferable of the attending specialist dealer/customer service) was set in the Settings menu → Service:

Company _____ Contact _____
 Street _____ Phone _____
 ZIP _____ Fax _____
 City _____ E-mail: _____

Settings in the user administration/logging

Admin PIN*: _____
 *Please keep the PIN in a secure and not accessible to public place.

Release process with PIN? ☐ YES ☐ NO
 Batch indication ☐ YES ☐ NO
 Batch release ☐ YES ☐ NO

The documentation software MELAview is used as full version ☐ YES ☐ NO

If yes:

- ☐ A new license code has been ordered and MELAview has been activated.
- ☐ The practice requests a new license code.

Batch documentation

☐ The function of the batch documentation was checked and works beyond

User-defined settings restored

☐ All settings lost as a result of the conversion have been restored (see commissioning, item 3 in the conversion instructions).

Label-Printer MELAprint[®] 60

Furthermore the label printer was connected. ☐ YES ☐ NO

If yes:

- ☐ The functionality of the label printer has been checked.
- ☐ The operating personnel were instructed

Test run with function test

- ☐ Test run with following log output performed.
- ☐ Function test also performed for connected accessories, e.g. log printer.

Electrical check in accordance with VDE 0701 carried out

☐ Electrical check in accordance with VDE 0701 carried out

Instructing the operating personnel and handing over the technical documentation

- ☐ Instructions to the batch documentation with CF card, for archiving and printing the files successfully.
- ☐ Operating personnel were instructed in new user guidance
- ☐ New User Manual was handed over

☐ New Technical Manual was handed over

Print a system and status log and attach it with the conversion protocol

- ☐ System and status log printed out

The device was properly upgraded by the above named technician

- ☐ without complaint
- ☐ with complaint

.....

Logs

Affix logs here and fold upwards. If no log printer is available, enter the successful completion of the program sequence.

Test run	Test label of MELAprint [®] 60 where applicable
System log	Status log

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

MELAG order-no. BS042+BS043

tools needed:	Torx TX20, Allen key 5 and 6 mm, clamping ring pliers A2, engineers pliers, hammer, arbor, sensor gauge 1,5 and 2 mm, sensor gauge 3 and 10mm
basic device state:	Door opened, power plug disconnected


Overview assembly places (device without back plane and casing):


Notice!

Mass spring water tank and hose to the elbow fitting „waste water“ are located on the right side of the device back plane.

	figures	working steps/comments
1st		Remove casing and back plane.
		Information for autoclave types: Vacuklav®41-B and 43-B 1. Remove the cover of the water tank
		2. Remove the wire frame out of the rubber fringe of the cover.
		3. Pull the fringe off the casing plate and push it inside the device. Remove the casing.


2nd		Loosen the door panel from the door bar and put it to the side. Remove the 2 fixing screws of the door plate (Torx TX20). Caution!  Notice the spring force of the parts which are located behind it. Notice!  Notice the assembling position of the pertinax-parts for the reassembling.
3rd		Remove the door bushing. Loosen the clamping ring (pliers A2) and remove all related parts. Caution!  Notice the spring force of the parts which are located behind it.
4th		If the safety pin of the door bushing is bend, you have to exchange it. Otherwise push the pin with the mandrel and hammer as far as there is no supernatant inside the door bar (figure from inside, view on the door bar).
5th		Screw the under item 3 removed door bushing clockwise on the assembled spindle as far as the axis split-pin is accessible (connection motorshaft/ door spindle).

6th		Remove the split-pin (engineers pliers). The split-pin is located behind the bearing block and connects the door spindle with the motor axis.
7th		Pull off the door motor in direction to back plane.
8th		Loosen the bearing block. You may have to loosen the cables which lay in front of it for the purpose of easier dismounting. It is illustrated the lower screw fitting of the bearing block. (Allen-key 5 mm, see arrow).

Information for autoclave types: Vacuklav®41-B and 43-B

1. For the purpose of reaching the upper bearing block easier you have to set up or remove the water tank. Therefore remove the rubber cover of the tank and drain it completely.


2. Unhinge the mass spring at the rear area of the water tank.




3. Loosen the hose connection from the tank to the elbow fitting „waste water“.




4. Lift the water tank at the front of the device front a bit until the upper bearing block screw fitting is accessible. Fix the tank with proper tools.



The upper bearing block screw fitting is accessible from the upper side of the device (view from the top inside the device, left side in direction to the display)


9th		Removal of the old door spindle. Loosen the door bushing that was screwed before and take it off with the bearing block, spring, ball bearing and bearing shell.  Caution! _____ Notice the spring force of the parts which are behind it.
10th		 IMPORTANT! The ball bearing has to be lubricated with the supplied grease Molykote Multilab (MELAG order no. 88925) before the re-assembling of the door spindle. Mount the new door spindle with the lubricated ball bearing, spring and bearing shells and push them towards the spring force to the front. As done with the removal screw the door bushing towards it and pull the spindle to the front as far as possible.  Caution! _____ Note the assembling order of the bearing shells: The big outer ring has to show to the bearing block, small outer ring to the door bushing.
11th	 	Fix the new bearing block with the two hex key screws and the lock washers. It is illustrated the lower screw fitting of the bearing block (Allen-key 5 mm, see arrow).  Caution! _____ Adjust the bearing block so that it is collateral to the motor shaft and does not cant.  Note! _____ Put same grease on the upper screw and then put it on the allen-key so that the screw cannot fall down. (view from the top inside the device, on the left side in direction to the display)

		Information for autoclave types: Vacuklav®41-B and 43-B 1. Drop the water tank again respectively mount it again, remount the mass spring again. 2. Fix the hose to the elbow fitting „waste water“ 3. Plug the rubber cover on the water tank again. Fill the right side of the tank with fresh water up to where it is marked.
12th		Guide the door motor with the axis into the door spindle. Caution! Before the assembly check the 3 motor guidings in the gearbox if they are fixed. Note! Adjust the spindle and motor axis by rotating the door bushing for the purpose of easier assembly.

		Put the split-pin in again. Unscrew the door bushing off the spindle.
13th		Adjust the door contact switch at the motor. The adjustment is done while the door is open (motor in rear position) by loosen the both fastening screws (see arrow, view from the rear side on the left side of the device) and by moving the switch in its slotted holes. While the door is open adjust the distance x between the rocker switch of the door contact switch and the motor at 3 mm (from SN 0540-B1337 and SN 0544-B1169 at 10mm) by using a sensor gauge. Move the rocker switch of the micro switch in direction to the chamber until the contact closes (you can hear a “click”). When moving the rocker switch back, the micro switch of the rocker switch has to close briefly before it hits the sensor gauge, if necessary adjust the position of the micro switch.
14th		Mount the door bushing Guide the bushing from outside through the door bar and put washer, spring, washer on it. Push it towards the spring force and place it at the snap ring. Push the snap ring in the key-way. Put same supplied oil on the door bushing.

15th		Connect the device with the power supply system and switch it on without closing the door panel. Start the vacuum test and find out the distance between door bushing and device front plate (figure shown without door for the purpose of better understanding).
		When required adjust the distance by loosen the fastening screw (Allen-key 5 mm) and by rotating the middle regulating screw* with the allen-key wrench 6 mm (but not the 4 mm screw!) on gauge dimension. from Software-Version 2.2 gauge dimension 2,0 mm before Software-Version 2.2 gauge dimension 1.5 mm *middle drill, door bar, see following page
		The middle regulating screw is located behind the middle drill of the door bar (see arrow). Start the vacuum test and check the distance again. If necessary repeat the procedure. Fix the fastening screw.
16th		Fix the door panel. Set the Pertinax-parts with the round side to the panel in and provide them with the springs. Mount the door plate with both screws.  Caution! _____ Note the correct applying of the screws, otherwise the door panel could be damaged.  Note! _____ Grease the rounded side of the pertinax- parts lightly.
17th	test run, visual inspection	

18th	Mount the back plane and casing
	Information for for autoclave types: Vacuklav[®]41-B and 43-B: 1. Plug the rubber fringe for the tank cover on the casing plate again.
	2. Put the wire frame in the fringe.
	3. Put the cover of the water tank on the device again.

Instructions for replacing the O-ring at the tube-heater

Included in delivery:

items	order-no.
1x O-ring for tube heater and Molykote grease M33	MELAG order-no. 45895
1x threaded pin 4x16	MELAG order-no. 44600
1x toothed lock washer	MELAG order-no. 22470
1x nut M4	MELAG order-no. 20661

Tools needed:

Phillips-Screw driver PH1, 7 mm socket key,
2mm allen-key, tweezers

Warning!

Before beginning with the works it is necessary to drain the double-jacket completely and to carry out a pressure release! Also see separate instructions in the appendix of the maintenance instructions "Pressure release and drain of the double-jacket". Check the pressure reading at the pressure gauge!

In case of non-observance you can suffer burns/serious injuries.


For the autoclaves Vacuquick®13-B/14-B please notice, that you have to dismount the vacuum pump before beginning, for it is necessary to reach the tube heater.

The already mounted Phillips-screw for fixing the tab washer is to be replaced by the threaded pin with nut when reassembling it.


Working steps:

Workings are executed while casing and back plane are removed.


1. Remove the four blue connections from the tube-heater.

2. Loosen the lower screw from the tab washer with a Phillips-screw driver and remove it.
3. Take off the tab washer carefully, please note the conductor of the temperature sensor.

4. Pull the temperature sensor out of the chamber. For a free access to the threaded hole of the tube-heater bend the temperature sensor aside.


5. Unscrew the tube-heater out of the chamber.


Warning!

Tube-heater may be hot.
Risk of burning oneself!


Before reassembling it grease the windings and the sealing cone with the supplied silicon-grease so you can see a light film of grease on it.

6. Remove the red old O-ring with the tweezers out of the nut of the inner thread of the chamber.

7. Grease the new black O-ring a little completely, so that there is a light glossy film on it, and put it in the notch instead of the old O-ring.


8. Screw the greased tube-heater with your hands in the hole of the chamber again as far as it will go.


9. Put the temperature sensor in again.

i Important!


Please note, that the threaded hole for the Phillips-screw has to lay below, because the drill also lays at the bottom of the tab washer! If necessary screw the tube-heater back a little.

10. Put the threaded pin on it and screw it in with the 2 mm allen-key.


11. Place the toothed lock washer and the nut on the threaded pin and tighten it with the socket key.

12. Plug the contacts on it again.

i Important!

Please notice that the wire of the temperature sensor has to be layed behind the notch of the tab washer as you can see on the left figure.

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Connecting the autoclave to a network

The best way to connect the autoclave to a network is with the aid of the computer support company with which your practice works.

You can directly connect a computer to the autoclave or over a network. A prerequisite is that the computer has a network card with a RJ45 connector jack (LAN). An FTP server (communications program) or an FTP service is required for connection over FTP. For log output via TCP there is no FTP program, but an applicable program, for example MELAview2 or MELAssoft is required.

Connect the Ethernet cable

Determining the IP address of the computer

Determine the network situation

Entering a fixed IP address

Specifying the computer for the FTP server program


Setting up the server (MELAG FTP server)

1. a) Autoclave – computer with 1:1 cable
b) Autoclave – network with crossover cable
If a so-called “intelligent switch” is used in the network, the cable connection is the same. Every interface can be employed at the autoclave.
2. Find out the IP address of the computer from the computer support administrator or determine this yourself (see question 6, page 28).
3. Fundamentally important: The chosen computer must always have a fixed IP address that is not dependent on automatically or manually configured networks.
For automatically configured networks, the DHCP server in the area must be informed of the range of the number or the number itself as a static IP address (es). The computer can also be assigned several IP addresses, in case the one already existing in the computer cannot be used.
The autoclave and the computer must belong to a subnetwork. In most cases this means that the IP address of the subnetwork must be the same for the first three numbers sets (e.g. 192.168.40.xx).
The IP address of the autoclave and the computer must differ in the fourth number set (e.g.: IP of the autoclave: 192.168.40.40 and IP of the computer: 192.168.40.140).
The default setting of the IP address in a MELAG autoclave meets these requirements.
4. Enter the fixed IP addresses for the autoclave and the connected devices into the autoclave (see **Verifying and changing IP addresses**, page 15).
5. Because the autoclave itself is equipped with a →hub, direct connection to a computer creates a network. For this reason the procedures for direct connection only minimally differ from connection over a network. In a network a computer must be specified on which the small program →FTP server can run. This program receives the logs via data transfer. The autoclave seeks this computer using the IP address that is set in the autoclave and logs on via the FTP server. The completed program logs are also later stored on this computer. In the selection of this computer, note that it is advantageous to have the stored logs included with the automatic data backup system of the practice.
6. MELAG offers a free, specifically developed FTP program for this. With the MELAG FTP server, it is possible to log on as a user at several devices simultaneously and to receive data from the autoclave and other devices, for example a washer-disinfector, in parallel. The FTP server supports the so-called multithread capability. You can determine for yourself the destination folders in which the device directory along with log data are stored in the FTP server program.

Connection to a practice network

Manually configured network

Automatically configured network (DHCP)

In the standard ex factory delivery, the autoclave contains IP addresses that all belong to the network with the subnetwork mask listed in Table 2. To integrate the autoclave into a medical practice network (and only then!), two possibilities are to be distinguished:

If the network is manually configured, then the subnetwork mask of the medical practice network must be known. If a subnetwork mask is used other than the one preset in the autoclave, a software expert should adapt the IP addresses in the autoclave. If the masks conform and the above-mentioned IP addresses are not assigned, then a log printer or label printer can be connected directly to the autoclave.

The network is automatically configured in a dynamic network.

- In this case, in the **Settings** → **Logging** menu set the IP address of the autoclave on DHCP.

If logging takes place over a computer via FTP, then a fixed IP address must be determined for the computer, which is then registered in the autoclave.

The IP addresses for the following devices are saved in the autoclave as parameters:

Device	IP address	Comment
Autoclave itself	192.168.40.40	Default setting ex factory
Computer	192.168.40.140	Default setting ex factory
Log printer	192.168.40.240	Default setting ex factory
Label printer	192.168.40.160	Default setting ex factory
Gateway	192.168.40.244	Not relevant within a network
Subnetwork mask	255.255.255.0	Possibly take over from the customer network

Table 2: IP addresses saved in the autoclave


Warning!

- The input of false IP addresses can lead to disturbance of the practice network.

Verifying and changing IP addresses

The IP addresses are changed in the **Settings** → **Logging** menu.


- Follow the logging assistance until you reach the window in which the IP addresses for the specific devices are listed.
- Select, for example, the autoclave. It opens the next window.
- Select the numeric keypad which you wish to change directly by choosing each respective number.
- Delete the digits with the **C** button, then enter a new number and confirm with **SAVE**.
- Proceed in the same way with the other devices that are to be connected to the network.


Connection via FTP or TCP

If you have selected the computer as the output medium in logging help (**Settings menu** → **Logging**), determine in the following window whether the connection will be via FTP or TCP.


For connection via FTP you must first set up a user name with password in the FTP program of the computer. The user name and password must then be registered in the autoclave in the **Settings menu** → **Logging**:

- Select Connection **via FTP**.
- In the lower button you will see the currently set user data (**standard username**: year of manufacture + serial number; **password**: MELAG12345). If you wish to change this, proceed to the settings window by pushing the buttons.

- Enter the user name and password and confirm with SAVE.

For connection via TCP you must register the TCP port number in the autoclave:

- Select Connection **via TCP**.
- In the lower button you will see the currently set TCP port (**standard TCP**: 65001).
- To change the TCP port indicated at the bottom, press the button. If you wish to change it, proceed to the settings window by pushing the buttons.
- Use the numeric keypad to select the corresponding TCP port. The current TCP port must first be deleted with the **C** button.
- Confirm with SAVE.

To output logs via TCP, an additional special program, e.g. MELAsot/MELAvi2 is required.

Frequently asked questions (FAQ) about the software

Warning

1 How can I format a CF card at the autoclave itself?

- All saved data on the CF card is deleted during formatting! Therefore you should check the CF card contents and save possibly existing backup records or other data on the PC or on some other storage medium.

Procedure for formatting

1. Insert the CF card correctly into the card slot (the noticeable raised rib at the edge points back and to the right). Never apply force.
2. Select Settings → **Format CF card** → in the window **Format CF card now?** Press **OK**
3. Confirm the confirmation message with **YES**.
4. If the message **Formatting CF card... completed** appears, you may remove the CF card.


2 How is a CF card formatted at the PC so that it is recognized by the autoclave?

The →CF card must be formatted with the file system FAT16 or FAT32. The autoclave can only store or read data on cards formatted in such a manner. Cards delivered by MELAG are already correctly formatted.

The formatting on the PC can be done without problems for example with Windows XP:

1. Plug the CF card into the MELAflash card reader which is installed and connected at the PC
2. Click **START**
3. Click **Execute...**
4. Enter the format command as specified below in the input field of the **Execute** window. X stands for the letter of the hard drive of the card reader, A determines the cluster size on the CF card, which varies depending on the memory capacity of the CF card.
5. Confirm the entry of the format command and if desired enter a name (max. 8 characters) for the CF card.


capacity of the CF card	format instruction to be entered
2 GB	format X: /A:32K
1 GB	format X: /A:16K

By adjusting parameter A you for instance prevent Windows XP from otherwise automatically formatting with FAT16 or FAT32 depending on the card size. We recommend using MELAflash CF cards with a minimum of 1GB memory capacity in accordance with the latest technical standards.

3 Where do I find the outputted logs?

After log output in the device directory

After a log output, a directory is automatically created with the encoded serial number of the respective autoclave. The name of the directory consists of five characters, which are identical to the first five characters of every log, e.g. 810RS. Under this directory there are subdirectories with the months of the log creation, e.g. 07_2010 for July 2010. The subdirectories contain all the logs created by the autoclave in the respective month. The device directory in the main directory is created on the CF card.

When using the MELAG FTP server

Direct transfer of the log on a computer via the network and using the MELAG FTP server allows you to determine directly in the FTP server where on your computer the device directory with log files is to be saved.

Do not rename directories

For every type of log output (immediate output after a completed cycle or transfer of multiple log records at once), the autoclave also verifies the storage medium and automatically creates a directory of the device and the month if none exists. Do not rename the directories, since then the logs would be stored both in the renamed directory as well as in the device directory automatically recreated by the autoclave.

Directory "doubled"

If multiple log records are outputted on one and the same storage medium, then in the device directory a directory "double" is created.

Position	1	2	3	4	5	6	7	8		X	X	X
Example	A	4	D	0	J	0	B	M	.	P	R	O
Meaning	Serial number					→Batch				Log file extension		
	Year of manuf.	Type	Manuf. No.									
Explanation	A...2010 B...2011 C...2012	0...40 B+ 1...41 B+ 3...43 B+ 4...44 B+								.PRO .STR .STB .LOG .STA .LEG .GPD (.DEM) (.DES)	successful program malfunction malfunction in standby system log status log abbreviation legend graphic log (demo log) (demo malfunction)	

Log files must be allocated to a text editor in order to open them easily

A computer is initially not aware of the extensions of the log files. If you permanently allocate a log file with the corresponding typical file extension to the text editor program, from then on you can always open this file type with a double-click. Alternatively, you can open the log files with the program MELAview2 (see question 4)

Do not rename log file names

A complete encryption of both the serial number as well as the total batch number is carried out directly in the eight-place log file name. A manual renaming of a file is therefore discernible, but not recommended. A file name is not allocated twice. The log file name ensures a good sortability of records.

Date and time of the log files

Date and time of the log files in Windows Explorer are identical with the date and time of the program start, provided that these were saved on the respective medium by immediate output. This information is lost on output collected later on a medium or when sent by e-mail.

Serial no. and total batch number are found in the log

It is not necessary to know the encryption within the name of the log file, since the serial number and total batch number can be seen by double-clicking on the file the contents. The prerequisite is the above mentioned allocation of the log file to a text editor (see question 4).

Total batch number must be entered when the circuit board is exchanged

In the rare event that a circuit board must be exchanged, the total batch number must be entered. This is possible for →authorized experts in the **Service** menu.

4 How can I open and print a log file? (MELAviiew2)

The program MELAviiew2 for displaying, editing and printing log files

MELAviiew2 is a program developed by MELAG for displaying and sorting sterilization logs. Furthermore, you have the option of logging graphic logs during a program run. Text and graphic logs can be printed from MELAviiew2 or exported in various data formats.

For test purposes you can display DUMMY logs without a license code. The complete autoclave log function of MELAviiew2 can only be used with a license which has to be purchased.

An installation assistant leads you through the installation of the program.

Regular transfer of log data from CF card to the PC

The MELAviiew2 program has no administrative functions. This means that the log files are not automatically synchronized from a CF card with the folder for log archiving on your computer.

We therefore recommend that you regularly (weekly or monthly) transfer the data from the CF card to your computer. To do this, create a folder on the hard drive of your computer named, for example, MELAG Autoclave Logs, where you can save and archive the log files. You can copy the log files from the CF card to the computer using Windows Explorer.

Functions of MELAviiew2

MELAviiew2 is able to display complete log lists. You can sort by double clicking on the headers within the list and filter the files via the selection field below.

A statistical analysis of the logs is carried out automatically (number of successful programs, program type etc).

You can save preferred paths by selecting the path and right-clicking the folder symbol.

Graphic logs

You can record graphic logs during a program run using MELAviiew2 and display them as curve diagrams in MELAviiew2. You can export the graphic logs from the program in various data formats, e.g. for archiving.

Log files can also be read with other text editors

All log files are text files in ASCII format, which you can read, edit, and print with the simplest text editors of every operating system or with text processing or table calculation programs.

Graphic logs can only be displayed with MELAviiew2.

Link with the Windows Editor

An operating system uses the file extension to recognize which program it should use to open or edit the file. Using the following example of the Windows editor, you will learn below how you link other Windows programs (here the Text Editor) with a certain file extension.

Your computer normally does not yet recognize the extensions of the log files explained under question 3. Therefore, you have to "tell" the operating system that a file with the extension PRO, STR or LOG, for example, is always to be opened with the Text Editor.

Example: file.END represents a log file whose extension "END" is not yet recognized by the operating system.

1. Double-click on the log file in the Windows Explorer.
2. If the extension is not known, Windows XP issues the following message. Choose **Select program from a list** and confirm with **OK**.
3. Select the Editor from the list of programs in the window that opens. It is important that you check the box for **Always open this file type with the selected program** so that the operating system will recognize the extension in the future as well. Confirm with **OK**.
(From now on you can always open files with this extension by double-clicking the Windows Editor.)

5 What do the terms DHCP, IP address and subnetwork mean?

IP address

The IP address is the identification of the respective computer in a network expressed in numbers and unambiguously identifies the respective computer or autoclave with four number blocks (e.g. 192.168.88.8).

Subnetwork and IP address

Every IP address is divided into a network and an appliance part (host part). The division is carried out by the network mask (subnetwork/Sub Net Mask). The network part of the IP address must be the same for devices to be able to communicate with each other on the network. With a network mask = 255.255.255.0 (the most frequent), the first three number sets (e.g. 192.168.88.x) must match. The device part of the IP address is assigned individually and only once. The first (network itself) and the highest (broadcast) device address may not be assigned here.

6 How do I ascertain the IP address and network settings of a PC (Windows XP)?

Determination of the IP address of a PC under Windows XP:

4. **Start** → **Settings** → **Network connections** → **LAN Connection** → **Register Network support** → **IP address**
5. For instance you can discern whether a PC is integrated into a dynamic network: **Start** → **Settings** → **Network connections** → right mouse button on **LAN Connection** → **Properties**: Scroll in the window **This connection uses the following elements** to **Internet protocol (TCP/IP)** and click with left mouse key, then press on button **Properties**. If in the window **Properties of Internet protocol (TCP/IP)** the option **obtain IP address automatically** is selected, the PC in the network is dynamically addressed.

7 How is the MELAprint[®]42 log printer directly connected to the autoclave?

!
Warning!

A direct connection of the log printer to the autoclave without additional adjustment of the IP address of the log printer is only possible under one of the following conditions:

- The factory-set IP address of the autoclave, or at least the first three number sets of the IP address, are retained (see p. 15, Table 2).
- If the autoclave is to be connected to a (practice) network the IP address must (the first three number sets) of the practice and the factory set IP address of the autoclave must be the same (see p. 15, Table 2).

Procedure when connecting the log printer

1. The actual network function is carried out by a MELAG network adapter (Art.-No. 40295).
2. Prerequisite is therefore the network adapter connected with the MELAprint[®]42. The MELAprint[®]42 is exchangeable at any time; re-install after an exchange of the network adapter.
3. Connect the network adapter with a network cable (crossover, standard RJ45 Ethernet cable) at the autoclave: MELAprint[®]42 and network adapter ↔ autoclave.
4. Switch on the autoclave, if it is not already switched on.
5. Set up the log printer in the autoclave (see User Manual).

If not successful

If the aforementioned steps do not lead to success, the IP addresses set as parameters in the autoclaves and the MELAprint[®] probably no longer belong to a subnetwork (see **Verifying and changing IP addresses**), because the original factory settings were changed.

8 How is the MELAprint[®]60 label printer directly connected to the autoclave?

!
Warning!

A direct connection of the log printer to the autoclave without additional adjustment of the IP address of the log printer is only possible under one of the following conditions:

- The factory-set IP address of the autoclave, or at least the first three number sets of the IP address, are retained (see p. 15, Table 2).
- If the autoclave is to be connected to a (practice) network the IP address must (the first three number sets) of the practice and the factory set IP address of the autoclave must be the same (see p. 15, Table 2).

IP address space of the (practice) network and the autoclave are different

If the factory set IP address space of the autoclave and the (practice) network do not correspond, both the IP address of the autoclave and the label printer must be changed in the (practice) network) as follows:

- Change the IP address of the autoclave according to the instructions on page 15, **Verifying and changing IP addresses**.
- Change the IP address of the label printer according to the instructions in the operating instructions for the label printer.

9 How do I register the autoclave as a user with MELAG FTP server?

Server installation

In the following example, the registration of an autoclave is described with the FTP server of MELAG as an example:

1. Copy the installation file of the FTP server to your computer and start the installation.

NOTICE!

In order to be able to install the MELAG FTP server, the free media library "Microsoft .NET Framework 2.0" or a higher version must be installed on your computer.


2. Start the program by left double clicking on the FTP symbol on your desktop.

In the register **user**, create the autoclave as a new user. User name = serial number of the autoclave, e. g. 201040-B1000; user password = MELAG12345. Click on the folder and choose the desired directory.

NOTICE!

Enter the four digits of the year of construction as the user name. Device labels show only the last two digits of the year of construction. Example: If the label contains "1040-B1000", enter exactly "201040-B1000" as the user name.

NOTICE!

Moving the FTP server to another folder on the computer results in the loss of the user settings.

Testing

3. The FTP server is started by double-clicking on the link symbol on the desktop or the FTP symbol in the Windows start menu. Depending on the computer's configuration, it might be necessary to inform the computer's firewall that the FTP server is an exception. In the **Logs** menu, you can create system and status logs, since these are also always available without a running program. The logs document all current settings of the autoclave and can be employed for testing. Furthermore make a test run with the immediate output of a log. In the event of problems with the FTP server it is advisable to send a ping to the IP address of the autoclave.


Figure 12: Sending a ping to the IP address of the autoclave on the Windows XP input prompt

Start the input prompt at the computer with the FTP server (**Start** → **programs** → **Accessories** (input prompt)). Enter the following lines: xxx.xxx.xxx.xxx (X = IP address of the autoclave)
If this ping is answered successfully as displayed here, then the data transfer between computer and autoclave has been established. Cable, network or the IP address of the autoclave can be excluded as causes for errors that can occur.


Messages which can enter in the window of the FTP server	Possible cause/significance
The following line remains on the screen for a longer time: ! Server started on Port 21	▪ Cable connection defective ▪ IP address of the computer not set correctly in the autoclave (see Verifying and changing IP addresses, page 15).
Over a longer time period, the following lines are repeated: < 14:12:22 192.168.40.100 USER: 201040-B1000 < 14:12:22 192.168.40.100 USER: 201040-B1000 << Incorrect login or password >>	▪ The password was not entered correctly
The following rows enter c. every eight seconds: < 14:03:22 192.168.40.100 USER: 201040-B1000 > 14:03:22 192.168.40.100 USER: 201040-B1000 << is logged in.>> > 14:03:52 192.168.40.100 USER: 201040-B1000 << C:\Documents and Settings\Admin\own files \protocols >> < 14:12:22 192.168.40.100 QUIT	▪ Successful →handshake between autoclave and FTP server in which however no autoclave protocol was imported


Table 17: Possible messages in the window of the FTP server and their significance

10 How can I verify the version of the autoclave software?

Version check in the Info and Status menu

You check the version of the software of the autoclave in the **Info & Status** menu on the running autoclave.

Legend DC outputs

OFF	0/3	3 =o.k.
ON	1/3	
ON	1/2	2=short-circuit
OFF	0/1	1=electr.
ON	1/1	interruption

*After having made changes in the SERVICE menu, the SETUP menu must be left via „MENU exit and SAVE“.

Hot keys	Current operating situation	Function/explanation
1 + 3	in the second display „We1 come“: picture with waterdrop, key 3 [DOOR] and key 4 [CONTINUE]	Skip start process: display immediately changes to the SETUP menu (only for 5 sec. available)
2+3	in all menus and displays	Status display with information about the autoclave including display of AIN, DIN, DINZ
1 + 2	SPECIAL menu → Conductivity selected	Quick-Test: shortened program run includes sterilization of 10 Sec. (only for testing!)
3 + 4	DOCU menu → 06 All program logs selected	All program logs are deleted from the internal log memory, a security query appears first
3 + 4	DOCU menu → 11 All error logs selected	All malfunction logs are deleted from the internal log memory, a security query appears first
2 + 4, then 4h	SETUP menu → DIAGNOSIS+SERVICE selected	Access to the DIAGNOSIS menu
3h, then 3 + 4	DIAGNOSIS menu → SERVICE menu (Parameters) selected	Access to the SERVICE menu, in which parameters can be set among others
4h	all menus, in which values can be set, e.g. serial number, batch counter, FTP-Server Password	Mark the value to be changed (text above key 4 changes from EDIT to OK) then hold key 4 [OK]. By pressing key 4 again shortly the setting is reset.
4h	DIAGNOSIS menu → AC outputs (AC-OUT)	Switch on the components; switch off by pressing key 4 again shortly.

Legend

123:456 Above right in the display is the index (digit) for display status to identify the operating situation

1 + 3 Simultaneously press key 1 + 3 shortly (< 1 Sec.)

1h Hold key 1 for at least 1 second

The device parameters can be changed in the **SERVICE** menu:

1. **MAIN** menu → **SETUP** menu → **DIAGNOSIS+SERVICE** menu

In order to get to the protected area, proceed as follows:

2. Press keys 2 + 4 simultaneously (short). The digit sequence in the top right corner of the display changes.
3. Hold key 4 until the display changes to the **DIAGNOSIS** menu.
4. In order to change the parameters in the **DIAGNOSIS** menu navigate to the **SERVICE** menu.
5. Hold key 3 until the digit sequence in the top right corner of the display has changed.
6. Select the **SERVICE** menu by short simultaneous pressing of the keys 3 and 4.
7. Select **03 Filter Parameters** in the SERVICE menu to change the parameters.

Type	Category	Number (v2.3x)	Number (v2.4x)	Name	Value	Description
Vacuum test						
07	08	016	017	pexvt	80	discharge pressure
07	08	017	018	pyxvt	300	evacuation pressure that has to be reached at least
07	08	027	028	pvtx	300	maximum pressure after evacuation
Conductivity						
00	05	046	046	km1	40	display message "Feed water quality declining"
00	05	047	047	km2	60	display message "conductivity insufficient"
Maintenance						
00	06	054	054	nwrt	3000/2000 (as of v2.3) 1000 (before v2.3)	batches until next maintenance
00	06	268	268	dwrt	730	days between maintenance
Temperature of cooling water vacuum pump						
00	11	080	080	tkw	30 (until v2.3) 27 (as of v2.4)	temperature from which SV cooling water opens
Lock programs						
01 - 11	02	004	004	fpren	00001	Lock program: set to 0
Run times						
01 - 06	15	005	005	zsr	210 - 1200	plateau times sterilization
Advanced vacuum tests as of v2.4x						
07	02	---	034	vtmod	1	standard vacuum test as well as in cold and in hot state
07	02	---	034	vtmod	2	vacuum test incl. double chamber (SV steam open)
Advanced vacuum tests for Vacuquick® 13-B, Vacuklav® 41-B, 43-B as of v2.4:						
07	02	---	034	vtmod	3	vacuum test incl. cooler (SV* vacuum 1 open)
07	02	---	034	vtmod	4	vacuum test incl. vacuum pump (SV vacuum 1+2 open)

*SV = solenoid valve

Control electronics-Main Board Version E (Art. No. 24196/24197)

For newly supplied control electronics with the Art. No. 24196 and 24197, there is no "onboard" firmware and user interface (UI). Only the current boot loader is present, which enables programming of control electronics on the Premium-Plus-Class, the Premium-Class and MELAtherm®10. An update CF card with the appropriate software is therefore required. The update files are available on the MELAG Internet page.

Premium-Class – Changing the control electronics-----

Update through the boot loader:

1. Download the relevant files from the download area of the MELAG Internet page.
2. Copy all necessary files to the main directory of the CF card.

<i>Firmware</i>	<i>FW_2431.mot</i>
<i>User interface</i>	<i>BUS_2431.mot</i>
<i>Display symbols</i>	<i>DSY_2420</i>
<i>Language 1</i>	<i>B1L_2431D.mot</i>
<i>Language 2</i>	<i>B2L_2431GB.mot</i>
<i>Parameter (Type-specific)</i>	<i>P40B_2432.mot</i>
<i>for the update using the boot loader</i>	<i>11030.mot</i>
3. Insert the CF card, but do not push it right in. Switch on the autoclave at the mains.
4. When the LED lights up, insert the CF card right in to the end.
5. Wait for the update. On no account should you pull out the CF card while an update is running!
6. Remove the CF card when "Welcome" appears on the screen.

Premium-Plus-Class – Changing the control electronics-----

Autoclave update on the colour touch display:

1. Look for the update files for the country concerned on the download area of the Internet page.
2. Save the 3 files concerned on the CF card.

<i>Firmware</i>	<i>CFW_3036.Mot</i>
<i>Device update file</i>	<i>UP40B_3032_XSa.MTK</i>
<i>Dongle file (Type specific)</i>	<i>xTyp_40B.Mtd</i>
3. Copy to the main directory of the CF card.
4. Insert the CF card into the card slot on the side of the colour touch display.
5. Switch on the autoclave at the mains, wait for the update.
6. Remove the CF card.

Premium-Plus-Class – Changing the colour-touch display -----

The user related settings (e.g. user account, contrast settings, etc.) are lost after this update. Custom set parameters are not lost.

1. Look for the update file for the country concerned on the download area of the Internet page.
2. Save the file concerned on the CF card.

<i>File</i>	<i>MR_3032_XSa.MTK</i>
-------------	------------------------
3. Copy to the main directory of the CF card.
4. Insert the CF card into the card slot on the side of the colour touch display.
5. Switch on the autoclave at the mains, wait for the update.
6. Remove the CF card.

Premium-Plus-Class – Complete Update-----

This update is required if a new software version was released. The user specific settings are adopted in this update.

1. Look for the update files for the country concerned on the download area of the Internet page.

<i>File</i>	<i>MR_3032_XSa.MTK</i>
-------------	------------------------
2. Copy to the main directory of the CF card.
3. Insert the CF card into the card slot on the side of the colour touch display.
4. Switch on the autoclave at the mains, wait for the update.
5. Remove the CF card.

Premium-Plus-Class – Changing the parameter file-----

Example: The device type has to be changed for a Premium-Plus-Class autoclave, e.g. from 44 B+ to 40 B+. In this case, a complete software update is required.

Autoclave update on the colour touch display:

1. Look for the update files for the country concerned on the download area of the Internet page.
2. Save the 3 files concerned on the CF card.

<i>Firmware</i>	<i>CFW_3036.Mot</i>
<i>Device update file</i>	<i>UP40B_3032_XSa.MTK</i>
<i>Dongle file (Type specific)</i>	<i>xTyp_40B.Mtd</i>
3. Copy to the main directory of the CF card.
4. Insert the CF card into the card slot on the side of the colour touch display.
5. Switch on the autoclave at the mains, wait for the update.
6. Wait for 7 seconds.
7. Switch off the autoclave at the mains.
8. Switch the autoclave on at the mains again.
9. Wait for the update.
10. Remove the CF card.

Trying a "simple" update the loading of the parameter file would be prevented through the dongle file. This update using the boot loader replaces the 44 B+ dongle file with the 40 B+ dongle file and thus loading the parameter file is made possible.

Premium-Plus-Class – Update or add a language-----

New languages to be added have to be activated after the update. The required activation code is available from MELAG.

Please request the desired language by E-mail, indicating the serial number.

1. Look for the update files for the country concerned on the download area of the Internet page.

<i>Language file</i>	<i>Lng_3032_XSa.MTK</i>
----------------------	-------------------------
2. Copy to the main directory of the CF card.
3. Insert the CF card into the card slot on the side of the colour touch display.
4. Switch on the autoclave at the mains, wait for the update.
5. Enter the relevant activation code, to activate the language.

Premium-Class – Conversion to Colour Touch Display-----

Currently, conversion to a colour touch display has only been approved for control electronics with the main board version E. If the control electronics previously had software for a black and white display, an update of the boot loader is required.

→ Update of the boot loader:

1. Copy 11031.mot to the CF card (Set properties "write protected" and "hidden").
2. Insert the CF card, but do not push it right in. Switch on the autoclave at the mains.
3. When the LED lights up, insert the CF card right in to the end.
4. Wait until the update is over.
5. LED blinks.
6. Switch off the autoclave at the mains.
7. Remove the CF card.

Autoclave update on the colour touch display

8. Look for the update files for the country concerned on the download area of the Internet page (Retrofit kit-download).
9. Save the 3 files concerned on the CF card.

<i>Firmware</i>	<i>CFW_3036.Mot</i>
<i>Device update file</i>	<i>UP40B_3032_XSa.MTK</i>
<i>Dongle file (Type specific)</i>	<i>xTyp_40B.Mtd</i>
10. Copy to the main directory of the CF card.
11. Insert the CF card into the card slot on the side of the colour touch display.
12. Switch on the autoclave at the mains, wait for the update.
13. Remove the CF card.

MELAtherm®10 – Changing the control electronics-----

Update using the boot loader:

1. Ask the MELAG staff for the respective files.
2. Copy all necessary files to the main directory of the CF card.

<i>Firmware</i>	<i>TF_1303.bin</i>
<i>Language 1</i>	<i>Z1L_1304D.mot</i>
<i>Language 2</i>	<i>Z2L_1304GB.mot</i>
<i>Parameter</i>	<i>TP10DTA_1304.mot/TP10DTB_1304</i>
<i>User interface</i>	<i>ZUS_1304.mot</i>
<i>for the update using the boot loader</i>	<i>11400.mot</i>
3. Insert the CF card into the card slot, but do not push it right in.
4. Switch on the device at the mains.
5. When the LED lights up, insert the CF card properly.
6. Wait for the update.
7. Remove the CF card after "Welcome" appears on the screen.

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E61	FW11 > no feed water the tank. Please fill up feed water tank!	To feed water- measuring turbine at the feed pump.	1. No water in the external supply tank 2. MELAdem 37/40/47 produces no water because the water tap was not opened	1. Feed pump defective 2. Solenoid valve feed water polluted/defective 3. Filter in conductivity cell clogged 4. Leakage in the intake path up to feed pump 5. Measuring turbine defective/blocked 6. Electrical connections to the measuring turbine or to the printed-circuit loose/cable breakage 7. Double jacket is under pressure and the feed pump has air sucked in, so that the pump cannot vent independently
E62	FW10: No feed water in the tank. Please fill up feed water tank!	Float-actuated level switch feed water DIN4	1. No water in the external supply tank	1. Float-actuated level switch in the supply tank defective 2. Cable connections at the float-actuated level switch or to printed-circuit board To 1. Test float-actuated level switch on digital inputs, DIN 6
E63	WW12: Conductivity insufficient! No start possible! Fill new feed water or exchange cartridge from MELAdem. Thereafter control conductivity measurement or restart.		1. Mixed-bed-resin pellets in the MELAdem 37/40 or 47 exhausted. 2. Water with inadequate water quality was poured into the supply tank	1. Short-circuit to the electrical connections of the conductance sensor 2. Tap water or other conducting liquid has dripped from above on the conductivity cell. 3. Conductance sensor in the conductivity cell is contaminated 4. Short- circuit in the wiring 5. Conductance measuring defective on the printed-circuit 6. Leakage current between transformer and electronic console because of too thin or damaged rubber plates
E64	FW12 > Conductivity insufficient	Uppermost limit value (par) at the conductance sensor	1. Mixed-bed-resin pellets in the MELAdem 37/40 or 47 exhausted. 2. Water with inadequate water quality was poured into the supply tank	1. Short-circuit to the electrical connections of the conductance sensor 2. Tap water or other conducting liquid has dripped from above on the conductivity cell. 3. Conductance sensor in the conductivity cell is contaminated 4. Short- circuit in the wiring 5. Conductance measuring defective on the printed-circuit 6. Leakage current between transformer and electronic console because of too thin or damaged rubber plates 1. 2. 4./5. withdraw cable connections from the conductance sensor, and thereafter from the AIN 11 to the printed- circuit and measure conductance with each. A measured conductance of about "0 µS" should be displayed. A raised value at removed socket plug AIN 11 indicates a faulty analogue measurement at the printed-circuit board. The printed- circuit is to be exchanged. Between cable conductivity cell and device housing no voltage may be measurable of > 1.5 V DC/AC, at greater voltages the digital input and outputs and analogue inputs are consecutively to be withdrawn with simultaneous voltage metering.
E65	WW13: Conductivity insufficient User measurement: Start no longer possible! Fill new feed water or exchange cartridge from MELAdem. Control conductivity measurement or restart.		1. Mixed-bed-resin pellets in the MELAdem 37/40 or 47 exhausted. 2. Water with inadequate water quality was poured into the supply tank	1. Short-circuit on the electrical connections of the conductance sensor 2. Tap water or other conducting liquid has dripped from above on the conductivity cell. 3. Conductance sensor in the conductivity cell is contaminated 4. Short- circuit in the wiring 5. Conductance measuring defective on the printed-circuit 6. Leakage current between transformer and electronic console because of too thin or damaged rubber plates
E66	FW15 > Leakage of distilled water	Feed water jacket during a pumping process leakage or level sensor error	None possible	1. Double jacket is flooded over, water is pressed out at the overpressure valve or by the open steam valve (below 95°C), defect at the level sensor (receptor, wiring, printed-circuit board), jacket permanently loses water by: 1. Leaky rinsing valve 2. Leaky overpressure valve 3. Leaky steam valve 4. Press-back from water by the feed pump

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E67	WW22: Rinse value insufficient. Were the prerequisites for the rinsing process truly fulfilled? If the rinsing process continues to be ignored, stains can develop and the device can deteriorate.		1. Bent waste water hose (watering units) 2.	1. Filter at the SV rinsing clogged 2. SV rinsing clogged or defective 3. Rinsing line from the jacket to the filter rinsing clogged or bent 4. Temperature sensors entry 3 defective or imprecise 5. SV rinsing is not activated by the electronics
E68	FW22 > Rinsing value insufficient TB-internal: This event should no longer come from the GS			This error no longer exists
E70	WW10: No feed water		1. No water in the external supply tank 2. Float-actuated level switch is jammed due to foreign body	1. Float-actuated level switch blocked 2. Float-actuated level switch defective 3. Wiring harness interruption 4. Intermittent contact at the digital input to the printed-circuit board
E71	WW10: No feed water in the tank .Next program start only possible after filling the storage tank!		1. No water in the external supply tank 2. Float-actuated level switch is jammed due to foreign body	1. Float-actuated level switch blocked 2. Float-actuated level switch defective 3. Wiring harness interruption 4. Intermittent contact at the digital input to the printed-circuit board
E72	WW14: Feed water quality declining! User measurement: poor feed water quality. Please guarantee materials quality as soon as possible! Starting the autoclave is still possible, but when insufficient quality this may no longer be possible.		1. Mixed-bed-resin pellets in the MELAdem 37/40 or 47 exhausted. 2. Water with inadequate water quality was poured into the supply tank	1. Short-circuit on the electrical connections of the conductance sensor 2. Tap water or other conducting liquid has dripped from above on the conductivity cell. 3. Conductance sensor in the conductivity cell is contaminated 4. Leakage current in the wiring 5. Conductance measuring defective on the printed-circuit board 6. Leakage current between transformer and electronic console because of too thin or damaged rubber plates
E73	WW15: Conductivity poor. Should the program nevertheless be continued? Improve quality as soon as possible! (Warning from START)		1. Mixed-bed-resin pellets in the MELAdem 37/40 or 47 exhausted. 2. Water with inadequate water quality was poured into the supply tank	1. Short-circuit on the electrical connections of the conductance sensor 2. Tap water or other conducting liquid has dripped from above on the conductivity cell. 3. Conductance sensor in the conductivity cell is contaminated 4. Leakage current in the wiring 5. Conductance measuring defective on the printed-circuit board 6. Leakage current between transformer and electronic console because of too thin or damaged rubber plates
E74	WW16: Conductivity poor CONFIRM: Feed water is nevertheless still used. Improve quality as soon as possible! (Warning from Standby)		1. Mixed-bed-resin pellets in the MELAdem 37/40 or 47 exhausted. 2. Water with inadequate water quality was poured into the supply tank	1. Short-circuit on the electrical connections of the conductance sensor 2. Tap water or other conducting liquid has dripped from above on the conductivity cell. 3. Conductance sensor in the conductivity cell is contaminated 4. Leakage current in the wiring 5. Conductance measuring defective on the printed-circuit board 6. Leakage current between transformer and electronic console because of too thin or damaged rubber plates

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E75	WW17: Feed water quality declining! Feed water quality is poor. Please replenish with new feed water or exchange the cartridge of the MELAdem as soon as possible! It is still possible to start the autoclave, but at insufficient quality it will no longer be possible.		1. Mixed-bed-resin pellets in the MELAdem 37/40 or 47 exhausted. 2. Water with inadequate water quality was poured into the supply tank	1. Short-circuit on the electrical connections of the conductance sensor 2. Tap water or other conducting liquid has dripped from above on the conductivity cell. 3. Conductance sensor in the conductivity cell is contaminated 4. Leakage current in the wiring 5. Conductance measuring defective on the printed-circuit board 6. Leakage current between transformer and electronic console because of too thin or damaged rubber plates
E76	WW21: Rinse value poor. Automatic cleaning has not taken place or was not sufficient. If this warning is repeated, please contact Service!			1. Clogged filter in rinsing line (all equipment) 2. Blockage in the pipe fitting rinsing below at the jacket (all equipment) 3. Defective or plugged rinsing valve (all equipment) 4. Blockage in the waste water lines for the siphon (all watering devices) 5. Blockage or calcification at the Y-piece wastewater (all watering devices) 6. Defective or plugged condensate emptying valve (all membrane equipment) 7. Clogged cooler or plugged pipe fitting at the cooler (all membrane equipment) 8. Plugged lines from the condensate drain valve to the cooling spiral (all membrane equipment) 9. Defective temperature sensor in the jacket receptor/entry 3
E78	FW31 > Waste water tank full or waste water drainage not possible! (Error at Start)	Float-actuated level switch waste water DIN5	Waste water tank was not or only incompletely emptied	1. Float-actuated level switch in the supply tank defective 2. Cable connections at the float-actuated level switch or to printed-circuit board 3. float-actuated level switch in the supply tank jammed
E79	WW33: Waste water tank full		Waste water tank was not or only incompletely emptied	1. Float-actuated level switch in the supply tank defective 2. Cable connections at the float-actuated level switch or to printed-circuit board 3. float-actuated level switch in the supply tank jammed
E80	WW32: Waste water tank full or waste water drainage not possible! (Warning from Standby)		Waste water tank was not or only incompletely emptied	1. Float-actuated level switch in the supply tank defective 2. Cable connections at the float-actuated level switch or to printed-circuit board 3. float-actuated level switch in the supply tank jammed
E81	FD1 > Blockage door motor	During locking process in gas-tight state (Z3)	Door was not pressed energetically enough, so that the thread failed to catch at entry	1. Damaged threaded inlet, usually the red bronze nut in the door beam but also the threaded spindle 2. Other blockages in the area of the thrust bearing of the tensioning spring as well as at the door motor 3. Defect of the current measurement at the printed-circuit board 4. Adjusted parameters in the parameter table for the current monitoring
E82	FD2 > Door Z2 Z1	During door opening (state Z1)	Door opening was obstructed by objects placed in front of the device. During the opening process the door was manually kept shut	1. Door switch that does not move easily 2. Door that does not move easily because: - Jamming between hinge at the beam - jamming between beam and door trim panel - jamming between door trim panel and other cladding elements 3. Chamber is still under pressure 4. Door seal stuck together 5. Threaded spindle in the nut jammed 6. Threaded spindle does not turn because wedge has fallen out on the shaft 7. defect at the motor
E83	FD3 > Door Z3 Z4	During locking process in the pressurized state Z4 (wedge, shaft)		1. Threaded spindle does not turn because wedge has fallen out at the shaft 2. Defect at motor 3. Threaded nut in the door beam turns because spiral pin has fallen out
E84	FD4 > Door Z4 Z2	During opening process for pressure compensation (state Z3)		1. Threaded spindle does not turn because wedge has fallen out at the shaft 2. Defect at motor 3. Threaded nut in the door beam turns because spiral pin has fallen out
E85	FD5 > Door shaft broken			1. Threaded spindle does not turn because wedge has fallen out at the shaft 2. Defect at motor 3. Threaded nut in the door beam turns because spiral pin has fallen out

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E86	MD11 > Start not possible: door open		For Vacuklav 13 and 14-B the door not properly closed ??? Clarify whether error is used at all	1. Defective door switch 2. Door switch wrongly adjusted 3. Cable connections at the float-actuated level switch or to the printed-circuit board
E87	MD2 > Please load device or close door immediately			
E88	MD33 > Please open door shortly and then press closed tightly			
E89	WA1 > Attention! No sterilization! Should program nevertheless be continued?		This warning comes when e.g. a disinfection program was selected	
E100	FRKP	Clarify with progpar and visio #135		
E102	FA8 > Cancel/pressure exhaust Monitoring time for pressure exhaust elapsed during the cancellation	At pressure discharge during the termination stage	Bent waste water hose, poor installation (numerous equipment attached to one siphon, waste water hose with large droop)	1. Clogged filter in pressure discharge line (all equipment) 2. Blockage in the pipe fitting of the pressure discharge at the chamber (all equipment) 3. Defective or plugged pressure relief valve (all equipment) 4. Blockage in the waste water lines for the siphon (all watering devices) 5. Blockage or calcification at the Y-piece wastewater (all watering devices) 6. Defective or plugged condensate emptying valve (all membrane equipment) 7. Clogged cooler or plugged pipe fitting at the cooler (all membrane equipment) 8. Plugged lines from the condensate drain valve to the cooling spiral (all membrane equipment) 9. Defective pressure sensor entry 7
E103	FB1 > Ventilation monitoring time	For ventilation causes: filter clogged		1. Clogged ventilation filter 2. Defective or plugged air-bleeding valve 3. Defective or plugged emergency discharge valve 4. Plugged or bent line to the chamber 5. Defective pressure sensor entry 7
E104	FB2 > Cancel/ventilation monitoring time	For ventilation during termination stage		1. Clogged ventilation filter 2. Defective or plugged air-bleeding valve 3. Defective or plugged emergency discharge valve 4. Plugged or bent line to the chamber 5. Defective pressure sensor entry 7
E106	FH1P > Under pressure Halt level	Pressure shortfall during the stage "maintain level"		Pressure drop evoked by: 1. Defective/plugged steam valve 2. Plugged/bent steam line 3. Bad heating power due to defect in the double jacket or burned cable 4. Defective steam producer relay 5. Pressure loss due to leaking/burst lines at the jacket or chamber 6. Defective pressure sensor receptor/entry 7. Defective analogue input pressure sensor on printed-circuit board
E107	FH2P > Overpressure Halt level	during the stage "maintain level"		Over-pressure evoked by: 1. Defective/pierced steam producer relay 2. Defective pressure sensor 3. Defective analogue input pressure sensor on printed-circuit board
E108	FH3T > Halt temperature difference chamber/mantle	Between temp. sensor1 (chamber) and temp. sensor3 (jacket) during the stage "maintain level"		1. Defective/plugged steam valve 2. Plugged/bent steam line 3. Defective temperature sensor in the jacket receptor/entry 3 4. Defective temperature sensor in the chamber receptor/entry 1 5. Defective analogue inputs 1 or/and 3 on printed-circuit board 6. Pressure loss due to leaking/burst lines at the jacket or chamber
E109	FRM1 > Overheating steam generator	Heating element DE (still not employed)		Is not implemented at this time
E110	FRM2: Monitoring time Temperature regulation	For attainment of the jacket temperature	Bad power supply (too small house connection to the power plug, defective wall-socket, numerous equipment connected to one wall-	1. Bad heating power due to defective heating in the double jacket or burned cable 2. Defective steam producer relay 3. Bad power supply (under voltage, bad installation) 4. Pressure loss due to leaking/burst lines at the jacket or chamber 5. Defective temperature sensor in the jacket receptor/entry 3 6. Defective analogue input temperature receptor/entry 3

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E111	FRM3: Monitoring time Pressure regulation	For attainment of the jacket pressure (not employed)	socket or fuse protection)	Is not implemented at this time
E112	FRM4: Over temperature Double jacket	Double-jacket		1. Defective/pierced steam producer relay 2. Defective temperature sensor in the jacket receptor/entry 3 3. Defective analogue input temperature receptor/entry 3
E113	FUN1 > Mains loss		1. Device was switched off during current program 2. Power connection not completely created 3. Power plug was pulled out 4. Power supply breakdown in the building supply line	1. Bad power supply (under voltage, bad installation) 2. Burned connections at the mains cable 3. Defective mains switch 4. Defective toroidal core transformer 5. Defective supply voltage monitoring on the electronic printed circuit board
E114	FR1 > Rinse monitoring time	For attainment of the rinsing temperature in the rinsing phase (criteria for discharge of the DE)		1. Clogged filter in rinsing line (all equipment) 2. Blockage in the pipe fitting rinsing below at the jacket (all equipment) 3. Defective or plugged rinsing valve (all equipment) 4. Blockage in the waste water lines for the siphon (all watering devices) 5. Blockage or calcification at the Y-piece wastewater (all watering devices) 6. Defective or plugged condensate emptying valve (all membrane equipment) 7. Clogged cooler or plugged pipe fitting at the cooler (all membrane equipment) 8. Plugged lines from the condensate drain valve to the cooling spiral (all membrane equipment) 9. Defective temperature sensor in the jacket receptor/entry 3
E115	WW3U: Cancel rinsing process Waste water tank was not completely emptied.			This message no longer exists
E116	FS1T > Sterilization Temperature difference	Between temp.sensor1 chamber and theoretical temp. (calculated from pressure of the pressure sensor) during the stage "sterilization"		1. Deviations of the theoretical and actual temperature are evoked by vacuum leakage or non condensable gases in the steam 2. Defective temperature sensor receptor/entry 1 3. Defective analogue input temperature sensor at printed-circuit board 4. Defective pressure sensor/entry 7 5. Defective analogue input pressure sensor at printed-circuit board
E117	FS2T > Temperature sensor 1 or 2 defective	between temp.sensor1 (chamber) and temp.sensor3 (jacket) during the stage "Sterilization"		1. Defective temperature sensor receptor/entry 1 2. Defective analogue input temperature sensor on printed-circuit board 3. Defective temperature sensor/entry 2 4. Defective analogue input pressure sensor on printed-circuit board
E118	FS3T > Sterilization Over temperature	During stage "sterilization" (temp.sensor1)		Over-pressure evoked by: 1. Defective/pierced steam producer relay 2. Defective temperature sensor receptor/entry 1 3. Defective analogue input temperature sensor on printed-circuit board 4. Defective pressure sensor receptor/entry 7 5. Defective analogue input pressure sensor on printed-circuit board
E119	FS4T > Sterilization Under temperature	During stage "sterilization" (temp.sensor1)		Under-temperature evoked by: 1. Shortfalls of the temperature are evoked by vacuum leakage or non condensable gases in the steam 2. Defective/pierced steam valve 3. Plugged/bent steam line 4. Bad heating power due to defective heating in the double jacket or burned cable 5. Defective steam producer relay 6. Pressure loss due to leaking/burst lines at the jacket or chamber 7. Defective pressure sensor receptor/entry 7 8. Defective analogue input pressure sensor on printed-circuit board
E120	FS5P > Sterilization Over pressure	The theoretical temp. (calculated from value of the pressure sensor)		Over-pressure evoked by: 1. Defective/pierced steam producer relay 2. Defective sensor/entry 7 3. Defective analogue input pressure sensor on printed-circuit board

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
		During stage "sterilization"		
E121	FS6P > Sterilization Under pressure	The theoretical temp. (calculated from value of the pressure sensor) during stage "sterilization"		Pressure shortfall evoked by: 1. Defective/plugged steam valve 2. Plugged/bent steam line 3. Bad heating power due to defect in the double jacket or burned cable 4. Defective steam producer relay 5. Pressure loss due to leaking/burst lines at the jacket or chamber 6. Defective pressure sensor receptor/entry 7. Defective analogue input pressure sensor on printed-circuit board
E122	FS7T > Sterilization Temperature difference Chamber/Jacket	between temp. sensor1 (chamber) and temp. sensor3 (jacket) during the stage "sterilization"		1. Defective temperature sensor chamber receptor/entry 1 2. Defective analogue input temperature sensor on printed-circuit board 3. Defective temperature sensor jacket receptor/entry 3 4. Defective analogue input temperature sensor on printed-circuit board 5. Defective/plugged steam valve 6. Plugged/bent steam line
E123	FT1 > Vacuum drying monitoring time	For attainment minimum pressure for vacuum drying.	1. Device is in too high ambient temperature 1. Device is installed and gets no or too little cooling air 3. Device is overloaded 4. Device was operated without loading rack so that e.g. textile load made direct contact with the chamber and large quantities of condensate were absorbed. This condensate vaporizes upon evacuation and forms large quantities of steam 5. Cooling air intake holes in the base plate are blocked by e.g. sucked-in paper	1. Clogged filter in vacuum line boiler SV vacuum 1 (all equipment) 2. Blockage in the pipe fitting vacuum to the chamber (all equipment) 3. Low performance/ blocked vacuum pump defective or plugged SV vacuum 1 (all equipment) 4. Blockage in the waste water lines for the siphon (all watering devices) 5. Blockage or calcification at the Y-piece wastewater (all watering devices) 6. Defective or plugged SV vacuum 1, secondary air is drawn on 3rd path (all membrane equipment) 7. Defective or plugged condensate emptying valve (all membrane equipment) 8. Clogged cooler or plugged pipe fitting at the cooler (all membrane equipment) 9. Leaky cooler (leakage) 10. Defective/blocked fan (AC output 4) at the cooler 11. Plugged lines from the condensate drain valve to the cooling spiral (all membrane equipment) 12. Defective pressure sensor entry 7
E124	FV5 > Vacuum rating (end pressure) insufficient	For attainment end pressure for evacuation/fractionating	See also FT 1 1. Device is in too high ambient temperature 1. Device is installed and gets no or too little cooling air 3. Device is overloaded 4. Device was operated without loading rack so that e.g. textile load made direct contact with the chamber and large quantities of condensate were absorbed. This condensate vaporizes upon evacuation and forms large quantities of steam 5. Cooling air intake holes in the base plate are blocked by e.g. sucked-in paper	1. Clogged filter in vacuum line boiler SV vacuum 1 (all equipment) 2. Blockage in the pipe fitting vacuum to the chamber (all equipment) 3. Low performance/ blocked vacuum pump defective or plugged SV vacuum 1 (all equipment) 4. Blockage in the waste water lines for the siphon (all watering devices) 5. Blockage or calcification at the Y-piece wastewater (all watering devices) 6. Defective or plugged SV vacuum 1, secondary air is drawn on 3rd path (all membrane equipment) 7. Defective or plugged condensate emptying valve (all membrane equipment) 8. Clogged cooler or plugged pipe fitting at the cooler (all membrane equipment) 9. Leaky cooler (leakage) 10. Defective/blocked fan (AC output 4) at the cooler 11. Plugged lines from the condensate drain valve to the cooling spiral (all membrane equipment) 12. Defective pressure sensor entry 7

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E125	FV0 > Vacuum rating (min. pressure) insufficient	For attainment minimum pressure for evacuation/fractionating	See also FT 1 1. Device is in too high ambient temperature 2.Device is installed and gets no or too little cooling air 3.Device is overloaded 4.Device was operated without tray rack so that e.g. textile load made direct contact with the chamber and large quantities of condensate were absorbed. This condensate vaporizes upon evacuation and forms large quantities of steam 5.This condensate vaporizes upon evacuation and forms large quantities of steam 5.Cooling air intake holes in the base plate are blocked by e.g. sucked-in paper	1.Clogged filter in vacuum line boiler SV vacuum I (all equipment) 2.Blockage in the pipe fitting vacuum to the chamber (all equipment) 3.Low performance/ blocked vacuum pump defective or plugged SV vacuum 1 (all equipment) 4.Blockage in the waste water lines for the siphon (all watering devices) 5.Blockage or calcification at the Y-piece wastewater (all watering devices) 6. Defective or plugged SV vacuum 1, secondary air is drawn on 3rd path (all membrane equipment) 7.Defective or plugged condensate emptying valve (all membrane equipment) 8.Clogged cooler or plugged pipe fitting at the cooler (all membrane equipment) 9.Leaky cooler (leakage) 10.Defective/blocked fan (AC output 4) at the cooler 11.Plugged lines from the condensate drain valve to the cooling spiral (all membrane equipment) 12.Defective pressure sensor entry 7
E126	FV8 > Cancel/Evacuation Monitoring time	For attainment end pressure for evacuation on termination	See also FT 1 and FV 0	See also FT 1 and FV 0
E127	FV9 > Monitoring time for vacuum drying elapsed	Only in the Quick- Program B or in the Quick-Program S, other programs run without faults, respective device types: Vacuklav 41-B, 43-B, Vacuquick 13-B		1st heavy plunger in the solenoid valve "Vacuum 2"; 2nd defective solenoid valve "condensate" 3rd leaky/ torn cooler 4th leaky connection wires between solenoid valve "Vacuum 1" and cooler and between cooler and solenoid valve "Vacuum 2" 5th general leakage of vacuum at the chamber 6th great quantities of non condensing gases inside the chamber For software version 2.3: If a mechanical fault is excluded, you can set the parameter values in the Quick-Program B or Quick-Program S (parameter type 02 and 03) in the parameter category 01 (delta, offsets, differences) in each case pd2tr from 700 to 250; only for Vacuklav® 41-B/43-B, Vacuquick® 13-B Notice: By changing the parameters the pump is wear out more, because the pump is activated earlier during the evacuation. so the pump is feeded with more heat. Usually the cooler takes on the the evacuation in that phase (until 600mbar), in which much steam and thus heat is evacuated out of the chamber. Because of that the parameter values should change only temporary, as long as the final fault cause is found .
E128	WA3 > Pressure exhaust tempo diminished! Please inform service soon and have repaired. Sterilization still possible. (area 3)		1. Bent waste water hose, 2. bad installation (numerous equipment attached to one siphon, 3.waste water hose with large droop)	1.Clogged filter in pressure discharge line (all equipment) 2.Blockage in the pipe fitting of the pressure discharge at the chamber (all equipment) 3. Defective or plugged pressure relief valve (all equipment) 4.Blockage in the waste water lines for the siphon (all watering devices) 5.Blockage or calcification at the Y-piece wastewater (all watering devices) 6. Defective or plugged condensate emptying valve (all membrane equipment) 7. Clogged cooler or plugged pipe fitting at the cooler (all membrane equipment) 8.Plugged lines from the condensate drain valve to the cooling spiral (all membrane equipment) 9.Defective pressure sensor entry 7
E129	WE3 > Steam intake tempo decreasing! Possible fault Causes: > Too low mains voltage under load > Capacity decline of the heating >		Bad power supply (too small house connection to the power plug, defective wall-socket, numerous equipment connected to one wall- socket or fuse protection)	1.Bad heating power due to defective heating in the double jacket or burned cable 2.Defective steam producer relay 3.Bad power supply (under voltage, bad installation) 4.Plugged SV steam 5.Defective solenoid coil SV steam 6.Blockage in the pipe fitting steam at the jacket 7.Blockage at the steam intake fitting at the clapper floor chamber

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
	Throughput of the steam valve too small (area 3)			
E131	LA2 > Steam exhaust performance declining		1. Bent waste water hose, 2. bad installation (numerous equipment attached to one siphon, 3. waste water hose with large droop)	1. Clogged filter in pressure discharge line (all equipment) 2. Blockage in the pipe fitting of the pressure discharge at the chamber (all equipment) 3. Defective or plugged pressure relief valve (all equipment) 4. Blockage in the waste water lines for the siphon (all watering devices) 5. Blockage or calcification at the Y-piece wastewater (all watering devices) 6. Defective or plugged condensate emptying valve (all membrane equipment) 7. Clogged cooler or plugged pipe fitting at the cooler (all membrane equipment) 8. Plugged lines from the condensate drain valve to the cooling spiral (all membrane equipment) 9. Defective pressure sensor entry 7
E132	LE2 > Steam intake performance declining		Bad power supply (too small house connection to the power plug, defective wall-socket, numerous equipment connected to one wall- socket or fuse protection)	1. Bad heating power due to defective heating in the double jacket or burned cable 2. Defective steam producer relay 3. Bad power supply (under voltage, bad installation) 4. Plugged SV steam 5. Defective solenoid coil SV steam 6. Blockage in the pipe fitting steam at the jacket 7. Blockage at the steam intake fitting at the clapper floor chamber
E133	LV2 > Vacuum performance declining		See also FT 1. Device is in too high ambient temperature 1. Device is installed and gets no or too little cooling air 3. Device is overloaded 4. Device was operated without loading rack so that i.e. textile load made direct contact with the chamber and large quantities of condensate were absorbed. This condensate vaporizes upon evacuation and forms large quantities of steam 5. Cooling air intake holes in the base plate are blocked by e.g. sucked-in paper	1. Clogged filter in vacuum line boiler SV vacuum 1 (all equipment) 2. Blockage in the pipe fitting vacuum to the chamber (all equipment) 3. Low performance/ blocked vacuum pump defective or plugged SV vacuum 1 (all equipment) 4. Blockage in the waste water lines for the siphon (all watering devices) 5. Blockage or calcification at the Y-piece wastewater (all watering devices) 6. Defective or plugged SV vacuum 1, secondary air is drawn on 3rd path (all membrane equipment) 7. Defective or plugged condensate emptying valve (all membrane equipment) 8. Clogged cooler or plugged pipe fitting at the cooler (all membrane equipment) 9. Leaky cooler (leakage) 10. Defective/blocked fan (AC output 4) at the cooler 11. Plugged lines from the condensate drain valve to the cooling spiral (all membrane equipment) 12. Defective pressure sensor entry 7

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E134	WV3 > Vacuum performance declining (area 3)		See also FT 1 1. Device is in too high ambient temperature 1. Device is installed and gets no or too little cooling air 3. Device is overloaded 4. Device was operated without loading rack so that e.g. textile load made direct contact with the chamber and large quantities of condensate were absorbed. This condensate vaporizes upon evacuation and forms large quantities of steam 5. Cooling air intake holes in the base plate are blocked by e.g. sucked-in paper	1. Clogged filter in vacuum line boiler SV vacuum 1 (all equipment) 2. Blockage in the pipe fitting vacuum to the chamber (all equipment) 3. Low performance/ blocked vacuum pump defective or plugged SV vacuum 1 (all equipment) 4. Blockage in the waste water lines for the siphon (all watering devices) 5. Blockage or calcification at the Y-piece wastewater (all watering devices) 6. Defective or plugged SV vacuum 1, secondary air is drawn on 3rd path (all membrane equipment) 7. Defective or plugged condensate emptying valve (all membrane equipment) 8. Clogged cooler or plugged pipe fitting at the cooler (all membrane equipment) 9. Leaky cooler (leakage) 10. Defective/blocked fan (AC output 4) at the cooler 11. Plugged lines from the condensate drain valve to the cooling spiral (all membrane equipment) 12. Defective pressure sensor entry 7
E135	FW40 > No cooling water! The vacuum pump has no cooling water. Please ensure sufficient water input! Open water tap or Aqua-Stop, check connections and pipelines.	To cooling water measuring turbine of the feed pump. The vacuum pump receives no cooling water. Please assure sufficient water flow! Open water tap or water-stop, verify connections and hoses. Revise error message depending on SW changes, new error FW41 if turbine too fast (air)	1. Water tap for cooling water supply not opened 2. Coolant hose bent 3. Shut-off cock in the basement turned off or filter in the basement clogged	1. Throttling nozzle in the SV cooling water is clogged 2. Solenoid valve cooling water polluted/defective 3. Filter in SV cooling water is clogged 4. Leakage in the intake path to feed pump 5. Measuring turbine defective/blocked 6. Electrical connections to the measuring turbine or to the printed-circuit board loose/cable break
E136	FA1 > Display over heated	Display	1. Ambient temperature too high 2. Ventilation slits covered 3. Door was permanently open	1. Fan behind the display defect/blocked 2. Fan to the back wall (watering units) defect/blocked 3. Defective temperature measurement on the display controller
E137	FIX1 > Short-circuit DC output 1	SV Vacuum 1		1. Short-circuit at the user (solenoid valve SV vacuum 1) 2. Short-circuit in the wiring harness to the solenoid valve 3. Short-circuit on the printed-circuit board at output 1
E138	FIX10 > Short-circuit DC output 10	SV Air injection		1. Short-circuit at the user (solenoid valve SV air injection 10) 2. Short-circuit in the wiring harness to the solenoid valve 3. Short-circuit on the printed-circuit board at output 10
E139	FIX11 > Short-circuit DC output 11	Door locking		1. Short-circuit at the user (magnet door locking) 2. Short-circuit in the wiring harness for door locking 3. Short-circuit on the printed-circuit board at output 1
E140	FIX12 > Short-circuit DC output 12	Ventilator housing ventilation at back wall		1. Short-circuit at the user (ventilator housing ventilation) 2. Short-circuit in the wiring harness to the ventilator 3. Short-circuit on the printed-circuit board at output 1
E141	FIX2: Short-circuit DC output 2	SV pressure discharge		1. Short-circuit at the user (solenoid valve SV pressure release) 2. Short-circuit in the wiring harness to the solenoid valve 3. Short-circuit on the printed-circuit board at output 1
E142	FIX3: Short-circuit DC output 3	SV steam		1. Short-circuit at the user (solenoid valve SV steam) 2. Short-circuit in the wiring harness to the solenoid valve 3. Short-circuit on the printed-circuit board at output 1
E143	FIX4: Short-circuit DC output 4	SV ventilation		1. Short-circuit at the user (solenoid valve SV ventilation) 2. Short-circuit in the wiring harness to the solenoid valve 3. Short-circuit on the printed-circuit board at output 1

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E144	FIX5: Short-circuit DC output 5	SV emergency discharge		1.Short-circuit at the user (solenoid valve SV emergency exhaust) 2.Short-circuit in the wiring harness to the solenoid valve 3.Short-circuit on the printed-circuit board at output 1
E145	FIX6: Short-circuit DC output 6	SV rinsing		1.Short-circuit at the user (solenoid valve SV rinsing) 2.Short-circuit in the wiring harness to the solenoid valve 3.Short-circuit on the printed-circuit board at output 1
E146	FIX7: Short-circuit DC output 7	SV Vacuum 2		1.Short-circuit at the user (solenoid valve SV vacuum 2) 2.Short-circuit in the wiring harness to the solenoid valve 3.Short-circuit on the printed-circuit board at output 7
E147	FIX8: Short-circuit DC output 8	SV Cooler discharge		1.Short-circuit at the user (solenoid valve SV cooler emptying) 2.Short-circuit in the wiring harness to the solenoid valve 3.Short-circuit on the printed-circuit board at output 1
E148	FIX9: Short-circuit DC output 9	SV MELAsteam		1.Short-circuit at the user (solenoid valve SV MELAsteam) 2.Short-circuit in the wiring harness to the solenoid valve 3.Short-circuit on the printed-circuit board at output 1
E149	FIY1: short-circuit AC output 1	Main heating		1.Short-circuit at the user (heating 2300 or 2000 W) 2.Short-circuit in the wiring harness to the heating 3.Short-circuit on the printed-circuit board at output 1
E150	FIY2: short-circuit AC output 2	Regular heating		1.Short-circuit at the user (heating 700 or 500 W) 2.Short-circuit in the wiring harness to the heating 3.Short-circuit on the printed-circuit board at output 2
E151	FIY3: short-circuit AC output 3	Vacuum pump		1.Short-circuit at the user (vacuum pump) 2.Short-circuit in the wiring harness to the vacuum pump 3.Short-circuit on the printed-circuit board at output 3
E152	FIY4: short-circuit AC output 4	(unfortunately not available)		This error is not implemented at this time, since the current consumption of the user is too small 1.Short-circuit at the user (fan cooler or SV cooling water) 2.Short-circuit in the wiring harness 3.Short-circuit on the printed-circuit board at output 3
E153	FIY5: short-circuit AC output 5	(unfortunately not available)		This error is not implemented at this time, since the current consumption of the user is too small 1.Short-circuit at the user (feed pump and SV feed water) 2.Short-circuit in the wiring harness 3.Short-circuit on the printed-circuit board at output 3
E154	FIZ1: Sensor error No.1	Temperature sensor 1 chamber (control)		1.Defective temperature sensor receptor/entry 1 2.Defective analogue input 1 on printed-circuit board
E155	FIZ2: Sensor error No.2	Temperature sensor 2 chamber (log)		1.Defective temperature sensor receptor/entry 2 2.Defective analogue input 1 on printed-circuit board
E156	FIZ3: Sensor error No.3	temperature sensor 3 jacket		1.Defective temperature sensor receptor/entry 3 2.Defective analogue input 2 on printed-circuit board
E157	FIZ4: Sensor error No.4	temperature sensor 4 vacuum pump		1.Defective temperature sensor receptor/entry 4 2.Defective analogue input 2 on printed-circuit board
E158	FIZ5: Sensor error No.5	Not occupied (reserved)		1.Defective temperature sensor receptor/entry 5 2.Defective analogue input 5 on printed-circuit board
E159	FIZ6: Sensor error No.6	Not occupied (reserved)		1.Defective temperature sensor receptor/entry 6 2.Defective analogue input 6 on printed-circuit board
E160	FIZ7: Sensor error No.7	Pressure sensor chamber		1.Defective temperature sensor receptor/entry 7 2.Defective analogue input 7 on printed-circuit board

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E161	FIZ8: Sensor error No.8	Not occupied (reserved)		1. Defective voltage input standby input 8 2. Defective analogue input 8 on printed-circuit board
E162	FIZ9: Sensor error No.9	Not occupied (reserved)		1. Defective voltage input standby input 9 2. Defective analogue input 9 on printed-circuit board
E163	FOX1: Opener DC output 1	SV Vacuum 1 (DCOUT1)		1. electr. interruption at the user (solenoid valve SV vacuum 1) 2. electr. interruption in the wiring harness to the solenoid valve 3. electr. interruption at the printed-circuit board at output 1
E164	FOX10: Opener DC output 10	SV Airinjection (DCOUT10)		1. electr. interruption at the user (solenoid valve SV Airinjection 10) 2. electr. interruption in the wiring harness to the solenoid valve 3. electr. interruption at the printed-circuit board at output 10
E165	FOX11: Opener DC output 11	Door locking (DCOUT11)		1. electr. interruption at the user (solenoid door locking) 2. electr. interruption in the wiring harness to the solenoid valve 3. electr. interruption at the printed-circuit board at output 11
E166	FOX12: Opener DC output 12	Ventilator housing ventilation at back wall (DCOUT 12)		1. electr. interruption at the user (solenoid door locking) 2. electr. interruption in the wiring harness to the solenoid valve 3. electr. interruption at the printed-circuit board at output 12
E167	FOX2: Opener DC output 2	SV pressure discharge (DCOUT2)		1. electr. interruption at the user (solenoid valve SV pressure discharge 10) 2. electr. interruption in the wiring harness to the solenoid valve 3. electr. interruption at the printed-circuit board at output 2
E168	FOX3: Opener DC output 3	SV steam (DCOUT3)		1. Electr. interruption at the user (solenoid valve SV steam) 2. Electr. interruption in the wiring harness to the solenoid valve 3. Electr. interruption on the printed-circuit board at output 1
E169	FOX4: Opener DC output 4	SV ventilation (DCOUT4)		1. electr. interruption at the user (solenoid valve SV ventilation 10) 2. electr. interruption in the wiring harness to the solenoid valve 3. electr. interruption at the printed-circuit board at output 4
E170	FOX5: Opener DC output 5	SV emergency discharge (DCOUT5)		1. electr. interruption at the user (solenoid valve SV emergency discharge 10) 2. Electr. interruption in the wiring harness to the solenoid valve 3. electr. interruption at the printed-circuit board at output 5
E171	FOX6: Opener DC output 6	SV rinsing (DCOUT6)		1. electr. interruption at the user (solenoid valve SV rinsing 10) 2. electr. interruption in the wiring harness to the solenoid valve 3. electr. interruption at the printed-circuit board at output 6
E172	FOX7: Opener DC output 7	SV Vacuum 2 (DCOUT7)		1. electr. interruption at the user (solenoid valve SV vacuum 2) 2. electr. interruption in the wiring harness to the solenoid valve 3. electr. interruption at the printed-circuit board at output 7

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E173	FOX8: Opener DC output 8	SV cooler discharge (DCOUT8)		1.electr. interruption at the user (solenoid valve SV cooler emptying 10) 2.Electr. interruption in the wiring harness to the solenoid valve 3.electr. interruption at the printed-circuit board at output 8
E174	FOX9: Opener DC output 9	SV MELAsteam (DCOUT9)		1.electr. interruption at the user (solenoid valve SV MELAsteam 10) 2.electr. interruption in the wiring harness to the solenoid valve 3.electr. interruption at the printed-circuit board at output 9
E175	FOY1: Opener AC output 1	Main heating (ACOUT1)		1.electr. interruption at the user (heating 2300 or 2000 W) 2.electr. interruption in the wiring harness to the heating 3.Electr. interruption at the printed-circuit board at output 1
E176	FOY2: Opener AC output 2	Regular heating (ACOUT2)		1.electr. interruption at the user (heating 700 or 500 W) 2.electr. interruption in the wiring harness to the heating 3.Electr. interruption at the printed-circuit board at output 2
E177	FOY3: Opener AC output 3	feed pump (ACOUT3)		1.electr. interruption at the user (feed pump/ SV feed water) 2.electr. interruption in the wiring harness to the feed pump/ SV feed water 3.Electr. interruption at the printed-circuit board at output 3
E178	FOY4: Opener AC output 4	(unfortunately not available)		1.electr. interruption at the user (Ventilator/SV cooling water) 2.electr. interruption in the wiring harness to the ventilator /SV cooling water 3.electr. interruption at the printed-circuit board at output 4
E179	FOY5: Opener AC output 5	Vacuum pump (ACOUT5)		1.electr. interruption at the user (vacuum pump) 2.electr. interruption in the wiring harness to the door locking mechanism 3.electr. interruption at the printed-circuit board at output 5
E180	FD14 > Door locking (14- B, 13-B) Move the hand wheel so that the locking pin is released. Please avoid touching the hand wheel after locking in the future!	For locking/unlocking the door	The locking pin was jammed due to twisting of the hand wheel	1. Locking pin does not move easily 2. Pull solenoid does not move easily 3.Defect to the coil of the pull solenoid 4.Switch for door locking is defective 5.Switch for door locking wrongly adjusted 6.interruption in the wiring harness between switch and printed-circuit board (digital input.... 7.Defect at the digital input..... on the printed-circuit board
E181	W-AT > Display too warm Please keep the door locked!		1.Ambient temperature too high 2. Ventilation slits covered 3. Door was permanently open	1.Fan behind the display defect/blocked 2.Fan to the back wall (watering units) defect/blocked 3.Defective temperature measurement on the display controller
E182	WMNU > Shortfall mains voltage		1. Device was not properly installed (too weak installation configuration) 2.Numerous equipment connected to one wall-socket	1.Bad power supply (under voltage, bad Installation) 2.Burned connections at the mains cable 3.Defective mains switch 4.Defective toroidal core transformer 5.Defective supply voltage monitoring on the electronic printed circuit board

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E183	WV4 > Vacuum performance low (area 4)		See also FT 1 1. Device is in too high ambient temperature 1. Device is installed and gets no or too little cooling air 3. Device is overloaded 4. Device was operated without loading rack so that i.e. textile load made direct contact with the chamber and large quantities of condensate were absorbed. This condensate vaporizes upon evacuation and forms large quantities of steam 5. Cooling air intake holes in the base plate are blocked by e.g. sucked-in paper	1. Clogged filter in vacuum line boiler SV vacuum 1 (all equipment) 2. Blockage in the pipe fitting vacuum to the chamber (all equipment) 3. Low performance/ blocked vacuum pump defective or plugged SV vacuum 1 (all equipment) 4. Blockage in the waste water lines to the siphon (all watering devices) 5. Blockage or calcification at the Y-piece wastewater (all watering devices) 6. Defective or plugged SV vacuum 1, secondary air is drawn on 3rd path (all membrane equipment) 7. Defective or plugged condensate emptying valve (all membrane equipment) 8. Clogged cooler or plugged pipe fitting at the cooler (all membrane equipment) 9. Leaky cooler (leakage) 10. Defective/blocked fan (AC output 4) at the cooler 11. Plugged lines from the condensate drain valve to the cooling spiral (all membrane equipment) 12. Defective pressure sensor entry 7
E185	FH1T > Monitoring time exceeded during heating	For heating up to target temperature/level	Bad power supply (too small house connection to the power plug, defective wall-socket, numerous equipment connected to one wall-socket or fuse protection)	1. Bad heating power due to defective heating in the double jacket or burned cable 2. Defective steam producer relay 3. Bad power supply (under voltage, bad installation) 4. Pressure loss due to leaking/burst lines at the jacket or chamber 5. Defective temperature sensor in the jacket receptor/entry 3 6. Defective analogue input temperature receptor/entry 3
E186	FE5 > Steam input speed insufficient	For attainment of goal pressure in the pressure build-up/steam intake in the chamber (conditioning, fractionation, heating up)	Bad power supply (too small house connection to the power plug, defective wall-socket, numerous equipment connected to one wall-socket or fuse protection)	1. Bad heating power due to defective heating in the double jacket or burned cable 2. Defective steam producer relay 3. Bad power supply (under voltage, bad installation) 4. Plugged SV steam 5. Defective solenoid coil SV steam 6. Blockage in the pipe fitting steam at the jacket 7. Blockage at the steam intake fitting at the clapper floor chamber
E187	FA5 > pressure exhaust tempo insufficient	For attainment goal pressure During pressure discharge chamber (conditioning, fractionation, pressure discharge)	Bent waste water hose, bad installation (numerous equipment attached to one siphon, waste water hose with large droop)	1. Clogged filter in pressure discharge line (all equipment) 2. Blockage in the pipe fitting of the pressure discharge at the chamber (all equipment) 3. Defective or plugged pressure relief valve (all equipment) 4. Blockage in the waste water lines for the siphon (all watering devices) 5. Blockage or calcification at the Y-piece wastewater (all watering devices) 6. Defective or plugged condensate emptying valve (all membrane equipment) 7. Clogged cooler or plugged pipe fitting at the cooler (all membrane equipment) 8. Plugged lines from the condensate drain valve to the cooling spiral (all membrane equipment) 9. Defective pressure sensor entry 7
E192	WA2 > Rinsing necessary For service-dependent rinsing process please completely fill feed water tank with distilled water!			
E193	WA3 > Rinsing necessary For service-dependent rinsing process, please completely drain waste water tank!			

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E195	FDXU > Error with the door system has been triggered by the automatic or manual door process, when no plausible switch settings (Door contact 1 and 2) are registered	Door contact 1 or door contact 2 (not plausible switch setting)		1. Defective door switch 1 2. Defective door switch 2 3. Cable break or loose contact in the connection lines between switch and electronics 4. Contact problems to the terminal connector between printed-circuit board 5. Contact switch 2 loose, causing regular unwanted setting of the switch point
E203	f.e203 Program log not activated.			
E204	f.e204 no memory for program log recording			
E205	f.e205 Program log recording not inactive			
E206	f.e206 Output device for program log unknown			
E207	f.e207 Data log recording not activated.			
E208	f.e208 No memory space for data log recording			
E209	f.e209 Data log recording not inactive			
E210	f.e210 Output device for data log unknown			
E211	f.e211 No memory space for system log recording			
E212	f.e212 System log recording not activated.			
E213	f.e213 Output device for system log unknown			
E214	w.e214 Other error (e.g. internal CFC under power, ...)			
E215	w.e215 File not available			
E216	w.e216 Subdirectory does not exist.			
E217	w.e217 File can not be opened (maximum number has been reached)			

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E218	FD 218 : Log already exists An attempt was made to overwrite an existing log with this name (write-protected) file. Empty data medium and new OUTPUT.			
E219	w.e219 Invalid file descriptor			
E220	w.e220 Filename too long or invalid chars...			
E221	w.e221 No more memory space on the CF-Card			
E222	w.e222 File / Directory with the same name exists.			
E223	w.e223 Unknown format or not formatted			
E224	w.e224 Invalid open mode			
E225	w.e225 Function not implemented			
E226	w.e226 CF-Card write buffer overflow			
E227	w.e227 Invalid drive letter			
E228	w.e228 For CFC identify after reset failed, e.g. if internal CFC is under power, ...			
E229	w.e229			
E230	w.e230			
E231	e 231: Warning No CF-Card in the slot. For immediate print=YES is required			
E232	w.e 232 Card still in Initializing or is still processing write operations ...			
E233	w.e 233			

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E234	w.e 234			
E235	w.e 235			
E236	w.e 236			
E237	w.e 237			
E238	w.e 238			
E239	w.e 239			
E240	w.e 240			
E248	MX1R > Repeat vacuum test Because of large temperature deviations, the vacuum test is repeated. Otherwise an exact test result would not be possible. These deviations come e.g. from the residual moisture in the vessel. Vessel must be empty and dry for the vacuum test!			Wavering temperature sensors chamber entry 1 or jacket entry 3
E249	MX1P : System leaking Maximum permitted subatmospheric pressure was exceeded. Vacuum test was not even attempted. Extensive leakage!			See also FT 1 1.Large leakage in the door area 2.Large leakage in the boiler area 3.Possibly steam streams on the leaking SV steam in the chamber 4.Leaky SV ventilation 5.Leaky SV emergency discharge 6.Leaky SV pressure discharge 7.Defective pressure sensor entry 7
E250	MX7N > MELAsteam is READY Required minimum temperature is achieved. MELAsteam can be used.			
E251	MX7U > MELAsteam: fallen below minimum temperature MELAsteam-SV is locked. MELAsteam is only then ready to run when the temperature is attained once again. A message will be displayed.		MELAsteam was employed for too long a period of time or run with screwed-off jet to the pistol	1.Bad heating power due to defective heating in the double jacket or burned cable 2.Defective steam producer relay 3.Bad power supply (under voltage, bad installation) 4.Plugged SV steam 5.Defective solenoid coil SV steam 6.Blockage in the pipe fitting steam at the jacket 7.Blockage at the steam intake fitting at the clapper floor chamber

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E253	FC1			
E254	M254: Calibration aborted			
E255	e255 EVT_ERR_SENSOR10 / Thermo element defective			1.Receptor error at the entry 10 Thermo element 2.error on the printed-circuit board at the entry 10 3. Thermo element at present not installed, error message should not come
E256	e256 EVT_ERR_SENSOR11 / Thermo element comparator defective			1.Receptor error at the entry 11 thermo element comparison measuring point 2.Error on the printed-circuit board at the entry 11
E257	e257 EVT_ERR_SENSOR12 / Conductivity sensor defective	at measurement	Water with very inferior water quality (tap water) was employed, possibly another liquid (for example disinfectant or developer) was employed	1.Short-circuit to the electrical connections of the conductance sensor 2.Tap water or other conducting liquid has dripped from above on the conductivity cell. 3. Conductance sensor in the conductivity cell is contaminated 4.Short-circuit in the wiring 5.Conductance measuring defective on the printed-circuit board 6.Leakage current between transformer and electronic console because of too thin or damaged rubber plates
E258	e258 EVT_ERR_SENSOR13 / Voltage monitoring defective			1.Receptor error at the entry 13 voltage metering 2.Error on the printed-circuit board at the entry 13 3.Receptor for the voltage monitoring is component of the printed-circuit board
E259	e259 EVT_ERR_SENSOR14 / DC Motor current measurement defective			1.Receptor error at the entry 14 current measurement door motor 24 VDC 2.Error on the printed-circuit board at the entry 14 3.Receptor for the current measurement is component of the printed-circuit board
E260	e260 EVT_ERR_SENSOR15 / AC Current measurement defective			1.Receptor error at the entry 15 current measurement AC outputs 1 to 5 heating to feed pump 230 VAC 2.Error on the printed-circuit board at the entry 15 3.Receptor for the current measurement is component of the printed-circuit board
E261	e261 EVT_ERR_SENSOR16 / Temp. heat sink monitoring defective			1.Receptor error at the entry 16 temperature measurement of the heat sink of the electronic relay AC outputs 1 to 5 heating to feed pump 230 VAC 2.Error on the printed-circuit board at the entry 16 3.Receptor for the temperature measurement is component of the printed-circuit board
E262	e262 EVT_ERR_SENSOR1DI FF / ADC defective	Temperature sensor 1 chamber (control)		Defect at the 14 or 10 bit AD converter in connection with sensor inputs 1 PT1000 chamber temperature
E263	e263 EVT_ERR_SENSOR2DI FF / ADC defective	Temperature sensor 2 chamber (log)		Defect at the 14 or 10 bit AD converter in connection with sensor inputs 2 PT1000 chamber temperature
E264	e264 EVT_ERR_SENSOR3DI FF / ADC defective	temperature sensor 3 jacket		Defect at the 14 or 10 bit AD converter in connection with sensor inputs 3 PT1000 jacket temperature
E265	e265 EVT_ERR_SENSOR5DI	temperature sensor 4 vacuum pump		Defect at the 14 or 10 bit AD converter in connection with sensor inputs 4 PT1000 vacuum pump

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
	FF / ADC defective			
E266	e266 EVT_ERR_SENSOR5DI FF / ADC defective	Not occupied (reserved)		Reserve not occupied
E267	e267 EVT_ERR_SENSOR6DI FF / ADC defective	Not occupied (reserved)		Reserve not occupied
E268	e268 EVT_ERR_SENSOR7DI FF / ADC defective	Pressure sensor chamber		Defect at the 14 or 10 bit AD converter in connection with sensor inputs 7 pressure sensor chamber
E269	e269 EVT_ERR_SENSOR8DI FF / ADC defective	Not occupied (reserved)		Reserve not occupied
E270	e270 EVT_ERR_SENSOR9DI FF / ADC defective	Not occupied (reserved)		Reserve not occupied
E271	e271 EVT_ERR_SENSOR10D IFF / ADC defective	??		1. Defect at the 14 or 10 bit AD converter in connection with sensor inputs 10 thermo element steam producer 2. Receptor for the temperature measurement is component of the printed-circuit board
E272	e272 EVT_ERR_SENSOR11D IFF / ADC defective	??		1. Defect at the 14 or 10 bit AD converter in connection with sensor inputs 11 Comparison measuring point on the printed-circuit board at entry 10 for thermo element steam producer 3. Sensor for the temperature measurement is a component of the printed-circuit board
E273	e273 EVT_ERR_SENSOR12D IFF / ADC defective	??		Defect at the 14 or 10 bit AD converter in connection with sensor input 12 conductance measuring
E274	e274 EVT_ERR_SENSOR13D IFF / ADC defective	??		1. Defect at the 14 or 10 bit AD converter in connection with sensor inputs 13 Voltage metering on the printed-circuit board 2. Receptor for the voltage metering is a+ component of the printed-circuit board
E275	e275 EVT_ERR_SENSOR14D IFF / ADC defective	??		1. Defect at the 14 or 10 bit AD converter in the connection with sensor inputs 14 current measurement of the 24 V DC output for the door motor 2. Receptor for the current measurement is a component of the printed-circuit board
E276	e276 EVT_ERR_SENSOR15D IFF / ADC defective	??		1. Defect at the 15 or 10 bit AD converter in connection with sensor inputs 15 current measurement to the electronic relay AC outputs 1 to 5 heating to feed pump 230 VAC 2. Receptor for the current measurement is a component of the printed-circuit board
E277	e277 EVT_ERR_SENSOR16D IFF / ADC defective	??		1. Defect at the 16 or 10 bit AD converter in connection with sensor inputs 16 temperature measurement at the heat sink of the electronic relay AC outputs 1 to 5 heating to feed pump 230 VAC 2. Receptor for the current measurement is a component of the printed-circuit board
E302	F302: EEPROM-Error (e302)			
E303	f.e303 EVT_ERR_RTC Clock			

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E304	f.e304 EVT_ERR_CSROM Checksum error	Checksum error		
E305	f.e305 EVT_ERR_DISPLAY Display data transfer	Display data transfer		
E306	f.e306 EVT_ERR_DATABASE1 / Error 1 in operator database.	Operating database error 1		
E307	f.e307 EVT_ERR_DATABASE2 / Error 2 in operator database.	Operating database error 2		
E308	f.e308 EVT_ERR_DATABASE3 / Error 3 in operator database.	Operating database error 3		
E309	f.e309 EVT_ERR_DATABASE4 / Error 4 in operator database.	Operating database error 4		
E310	f.e310 EVT_ERR_DATABASE5 / Error 5 in operator database.	Operating Database error 5		
E311	f.e311 EVT_ERR_DATABASE6 / Error 6 in operator database.	Operating database error 6		
E312	f.e312 EVT_ERR_DATABASE7 / Error 7 in operator database.	Operating database error 7		
E315	M315: Parameter write protected			
E316	FUDP > Cancel pressure compensation	For attainment pressure compensation between chamber and jacket during termination		1.Defective or plugged SV steam (all equipment) 2.Defective pressure sensor entry 7 3.defective temperature sensors in the jacket entry 3.
E350	f.e350 EVT_ERR_TIMEBASE			

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E351	W351: Please maintain! For conservation of value and failure-free operation this autoclave should be regularly maintained by a specialist dealer or from MELAG customer service. You are still able to start the autoclave.			
E353	W353.e: A blackout during change of parameters. Please do not turn off during SETUP. Please notify technical service if incidence recurs.			
E354	e354 EVT_ERR_SRAM			
E355	e355 EVT_ERR_HARDWARE	Hardware error		
E357	e357 EVT_ERR_CSROMBOO TLOADER / Boot loader storage defective. Please inform Service!	Boot loader memory		
E358	e358 EVT_ERR_CSROMAPP LIKATION / Program storage defective. Please inform Service!	Program memory		
E359	e359 EVT_ERR_CSROMCALI BDATA / Calibration data defective. Please inform Service!	Calibration data		
E360	e360 EVT_ERR_CSROMPRO GPAR / program/Parameter storage defective. Please inform Service!	Program/parameter memory		

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E361	e361 EVT_ERR_CSROMUSE R / User interface storage defective. Please inform Service!	User interface memory		
E362	e362 EVT_ERR_CSROMLAN GUAG1 / Language 1 memory defective. Please inform Service!	Language 1 memory		
E363	e363 EVT_ERR_CSROMLAN GUAG2 / Language 2 memory defective. Please inform Service!	Language 2 memory		
E364	e364 EVT_ERR_CSROMPRO GPROT / Program log storage defective. Please inform Service or delete program log!	Program protocol memory		
E365	e365 5EVT_ERR_CSROMST OERPROT / Program log storage defective. Please inform Service or delete error log!	Error logging memory		
E366	e366 EVT_ERR_STOERPRO TNOTACTIVE / Error log not active. Data cannot be stored.	Error logging		
E367	e367 EVT_ERR_STOERPRO TOVERFLOW / Malfunction log storage full. data cannot be stored.	Error logging memory		
E368	e368 EVT_ERR_STOERPRO TNOTINACTIVE / New error log can not be created. Logging still running.			

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E369	e369 EVT_ERR_STOERPRO TUNKNOWNDVICE / error log output not possible. Output device unknown.			
E377	W377: No printer is connected! Please connect a printer and configure it via the SETUP-menu respectively			
E378	W378.e : Battery voltage too low! Please notify technical support!			1. Battery worn out (normal fatigue after about 5 years) 2. Battery must continually be replaced (leakage current on the printed-circuit board lead to rapid discharging) 3. Contact resistance to the contacts to the battery
E378	W378: Battery voltage too low Please inform technical services.		1. Battery exhausted (normally after c. 5 years) 2. Battery must be replaced constantly (creepage current on the main board results in rapid discharge) 3. Transfer resistance on the battery contacts	1. The lithium battery (CR 2032 3V) buffers the voltage provision upon disconnection of the device from the network. If the battery is exhausted, the message E378 will appear. Information such as the date and time are set to the year 2000 and the maintenance message E351 will be displayed upon the next program start. The device can be operated further. Check or reset the date, time and printer settings after replacing the battery. 2. Creepage current on the electronics can result in frequent battery discharge. Upon repeated battery failure, replace the electronics. 3. Transfer resistance or a bad connection on the battery contacts can result in battery malfunction. Check the contacts on the battery bracket and clean with alcohol or spirit where necessary. Disconnect the device from the network before doing so.
E386	WMP8: Internal log memory almost full. This means that the memory in the autoclave (not the CF-Card) will soon be unable to accept any more programs. Ensure the availability of the output media for the charge documentation.			
E387	WMS8: Internal malfunction log memory almost full. This means that the memory in the autoclave (not the CF-Card) will soon be unable to save any more malfunction logs. Ensure the availability of the output media for the charge			

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
	documentation so that important logs can be archived for the technicians.			
E394	WCF8: Should all logs not issued on the CF-Card be saved on the CF-Card?			
E395	WPR8: Should all the logs not printed be printed?			
E396	WPC8: Should all the logs not saved on the computer now be saved?			No device fault, A sterilization process was completed without a connection to the pre-set medium having been established.
E397	WPC0: Connection to the computer not possible or aborted			
E402	FD08: Door was blocked by a foreign object and thus was re-opened. Remove the foreign object and clear.	1. foreign object (e.g. hand) between the door and front plate prevents the door from being closed correctly	1.tight door 2.defective door contact 3 3.incorrectly set door contact 3 4.Cable break to the door contact 3	1. The door spindle and the door nut not optimally aligned or are stiff, possibly due to a lack of oil / grease in the ball bearing. The plastic disk on the spindle end in the cage also requires lubrication (noise)
E403	FD30: Door malfunction No plausible switch setting from door contact 3		Only with T13/T14	
E404	F404: Sleep mode activated due to lack of cable bridge on DIN between "-" und "0"		1. Cable bridge on main board interrupted or entirely missing. 2.Loose contact in the connector	1. DIN "-" and "0" are connected with each other independently of software (white cable). This cable bridge is essential. 2. Check the connector DIN.
E407	FD23: The door does not reach gas-tight state Z3 (from Z2 - Z3)	1. The door was not pressed against the spindle beginning tightly enough. The spindle was unable to catch.	1. The thread spindle does not twist through; the shaft pin is missing. 2.Malfunction on the motor / pressure switch 3.Threaded nut in the door bar rotates because the spiral pin has fallen out / is defective	Switch settings: Z1 – door open; Z2 – Door beaked , door switch 1 activated / closed, door switch 2 not activated / open; Z3 - Door closed steam tight, door switch 1 and 2 activated / closed; Z4 - Door closed pressure tight, door contact switches 1 and 2 activated / closed and power off (int2 5A). 1. The thread spindle and the motor shaft are connected via a pin. If this pin is defective or is missing, the pin cannot turn. The pin should be replaced. 2. Check the power supply DCOUTB2, The LED illuminates when the door closes and a voltage of c. 24-28 VDC can be measured at X30. For safety reasons, the door is secured with a 0.2 bar pressure switch. If the pressure switch is open or defective, the voltage supply to the motor is interrupted. The door motor can only be operated with a chamber pressure under 0.2 bars. The pressure switch may be bypassed for measurement purposes only. The voltage can be measured directly at the motor on the cable ties. A defective motor can be removed from its guide after removing the retaining ring and disconnecting the cable ties. 3. The threaded nut / red brass nut is secured in the door bar with a spiral pin (first remove the door cladding). If the pin is defective or has fallen out, the threaded nut rotates in unison. A new spiral pin needs to be fitted.

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E408	FW41: No cooling water. Pump draws in air. The vacuum pump does not receive any cooling water. Please ensure sufficient water supply. To do so, open the tap or check the aquastop, connections and hose.	1. The tap is not opened: the vacuum pump draws in air from the duct aerator.	1. The measurement turbines signals too high a flow rate by bleed air in front of the measurement turbine. 2. The wrong calibration factor has been set in the parameter table. 3. Faulty turbine 4. Faulty counting input on the main board. 5. Cooling water valve flow restrictor defective.	1. The measurement turbine also registers turbine "wheel spin." With measurement values over 70ml/min it is assumed that there is air in the turbine. Defective cooling water feed, bleed air or no cooling water. 2. 3./4. The turbine registers too many impulses, although water is flowing through the device. Let the water outflow flow into a bucket and compare the manually measured value with the counter display. If the value displayed is considerably higher, e.g. > 70 and the measured value lies at e.g. 2.4 L/min (counter value 40), the turbine or the analogue input on the main board is defective. 5. If the flow restrictor in the cooling water valve (2 L/min green) is defective, or not present, too much water flows through the measurement turbine and counts too high a through flow, the value lies at over 70ml/min, air will be assumed. This fault is rare, as the pump also restricts the through flow. Readings produced a flow of 60-65 ml/min. This varies between vacuum pumps.
E414	FV1: Evacuation aperture blocked (gradient). With type 40-B, 44-B, 14-B, check effluent hose for kinks or pinching. Check for packaging residue in the chamber. Vacuum system outage.	1. Kinked effluent hose with water ring device 2. The surrounding temperature of the device is too high (membrane device) 3. Device is fitted and receives too little or no cooling air (membrane device) 4. Device is overloaded 5. Device was operated without a tray mount so that e.g. the textile load has come into direct contact with the chamber and has absorbed large quantities of condensate. This condensate evaporates upon evacuation and forms large quantities of steam 6. The cooling air ingress holes in the base plate are covered e.g. with paper sucked on to the holes	1. Blocked filter in chamber vacuum conduit SV vacuum I (all devices) 2. Blockage in the vacuum nozzle on the chamber (all devices) 3. Inefficient / blocked vacuum pump, motor protection, mech. defective or blocked SV vacuum 1 (all devices) 4. Blockage in the effluent conduit to siphon (all water-ring devices) 5. Blockage / calcination on the effluent Y-piece (all water-ring devices) 6. Defective or blocked SV vacuum 2, bleed air is drawn via the 3. route (all membrane devices) 7. Electric. Defective / blocked condensate drain valve (all membrane devices) 8. Blocked cooler or blocked nozzles on the cooler (all membrane devices) 9. Leaky cooler (leakage) 10. Defective / blocked ventilator (ACOUT4 output 4) on the cooler 11. Blocked conduits from the condensate drain valve up to cooling coil (all membrane devices) 12. Defective pressure sensor input 7, see E160 13. Defective vacuum pump kickback valve (Water ring device)	NOTE: The abbreviation FV0 has been introduced from version 1.023. (previously only: Insufficient vacuum rating (minimum pressure)). 1. Open the filter from the SV vacuum and clean. With a free vacuum conduit, the program step - evacuation time SF12 for achieving a pressure of 500mbar - amounts to 20 sec. 2. The vacuum nozzle on the rear of the chamber must be kept free of residue. In particular, packaging residue can hinder pressure release. 3. With water ring devices, a blocked vacuum pump can trigger the motor protector and the message E179 FOY5 open output. Resetting the motor protector and activating the vacuum pump returns the device to operation (control the amount of cooling water). A blocked membrane pump is not subject to electrical evaluation and can only be detected by visual or acoustic inspection or indicated by unchanged vacuum pressure. Blockage / a blocked plunger means that sufficient evacuation can no longer be guaranteed. The SV vacuum 1 can be opened in the device using a 24 mm spanner. Check the plunger and the position of the SV and clean. It could be necessary to change the SV. 4. The hose between the effluent exit on the Vacuklav and the siphon must be kink-free and taut. Obstruction resulting from blockage must be checked and removed. The hose connection on the siphon must be kept free from residue. A backflow preventer can be fitted. However, with siphons with two parallel hose connections and two fitted backflow preventers, it is necessary to check for reciprocal influence during pressure-release. Example: A thermo-disinfector and Vacuklav with simultaneous pressure release / water discharge. 5. Visual control of the connections at the Y-piece for possible visible blockage. Remove the drain hose on the rear of the device. Examine the screw connection with a torch. Guide a flexible cable section of approx. 2.5mm² through the screw connection from outside until it is visible inside the PTFE hose. In this way it is possible to recognize a blockage and possibly rectify it. Alternatively, you can replace the Y-piece. 6. The solenoid valve vacuum 2 is a 3/2 way valve. Upon being switched on, the vacuum pump will be supplied with bleed air via the open current-free path for a short time and can then start freely. A blockage or a blocked plunger means that sufficient evacuation can no longer be guaranteed. The SV can be opened in the device with a 24 mm spanner. Check the plunger and the position of the SV and clean. It could be necessary to change the SV. 7. Blockage / a blocked plunger means that sufficient condensate drain can no longer be guaranteed. The SV condensate drain can be opened in the device with a 24 mm spanner. Check the plunger and the position of the SV and clean. It could be necessary to change the SV. 8. A permeable cooler is necessary for error-free vacuum generation and pressure release (rinsing). Blocked coolers initially result in extended operating times, especially whilst generating the vacuum. Remove the cooler and rinse with water. A clogged cooler should be changed. 9. Perform a normal vacuum test and an extended vacuum test (from version 3.0 and above). If the leakage rates vary strongly, (e.g. 0.3 and 1.5) the cooler has a leak. During an extended vacuum test, the vacuum valve 1 DCOUT1 is also

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
				activated. The cooler is thus included in the test. A leaky cooler can also be controlled in excess pressure. The steam or condensate leaves (e.g. through the screw connection) upon pressure-release or rinsing. A leaky screw connection can be sealed using Loctite. A defective / torn cooler must be replaced. 10. A defective ventilator means that the vacuum pump is not cooled sufficiently. The high temperature has an extreme effect on the performance of the pump especially after the second evacuation. Its functionality should be checked via activating ACOUT 4 in the SERVICE menu. Proceeding from the electronics, LED output ACOUT4, the 230V voltage and the cable connection should be checked up to the ventilator. 11. To check the conduits, open the hose connection to the hose filter at the SV condensate drain. Activate the DC OUT 08 SV cooler drainage in the SERVICE menu. You can check for permeance by blowing lightly into the cooling coil in the effluent tank. The air movement can be felt at the SV exit. When blocked, the corresponding screws or hoses may need to be changed or unblocked. 12. Imprecise measurement by the pressure sensor (sensor defective outside the area) can result in the minimum evacuation pressure of 600 mbar not being reached. The P sensor is checked for plausibility in the SERVICE menu under analogue inputs AIN 07. Normally, the surrounding pressure is around 1000 mbar (+/- 50 mbar). Upon considerable deviation, the sensor should be replaced. 13. The vacuum kick-back valve is located on the suction conduit of the vacuum pump. Upon an error, the plunger jams and the under-pressure remains unchanged during evacuation, remaining e.g. at 800 mbar. A light blow to the kick-back valve returns the device to evacuation immediately. The kick-back valve was jammed. The kick-back valve is dismantled, the plunger taken apart and re-fitted without a plunger (cement with Loctite). Alternatively, the kickback valve can be replaced by an 8 mm straight socket-end fitting.
E416	WCFD: No CF card for log			Insert CF card. This could also be unreadable. In such a case, save the log, re-format the CF card and perform a test print.
E417	WFTD: No computer connection for the log			
E422	WDT: Invalid date or time			Reset date and time

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E425	FTUN: Unknown type			
E426	FONA: Object not active			
E428	FG1: Pressure-release inlet blocked (gradient). With type 40-B, 44-B, 14-B, check effluent hose for kinks or pinching. Check for packaging residue in the chamber,	1. Kinked effluent hose 2. Bad installation (multiple devices on a single siphon, sagging effluent hose) 3. Pressure release clips blocked by packaging residue	1. Filter in pressure release pipe is blocked (all devices) 2. Blockage in clips pressure release on the chamber (all devices) 3. Mechanically defective or blocked pressure-release valve (all devices) 4. Blockage in the effluent conduit to siphon (all water-ring devices) 5. Blockage / calcinations on the effluent Y piece (all water-ring devices) 6. Mechanically defective / blocked condensate drain valve (all membrane devices) 7. Blocked cooler / clips on the cooler blocked (all membrane devices) 8. Blocked conduits from the condensate drain valve to cooling coil (all membrane devices) 9. Defective pressure sensor input 7	1. Open the hose filter SV pressure release and clean. 2. The pressure-release clips on the floor of the chamber must be free of residue or deposits. In particular, packaging residue can hinder pressure release. 3. Blockage / a blocked plunger means that sufficient pressure-release can no longer be guaranteed. The SV pressure release can be opened in the device with a 24 mm spanner. Check the positions of the plunger and SV and clean. It could be necessary to change the SV. 4. The hose between the effluent exit on the Vacuklav and the siphon must be kink-free. Obstruction resulting from blockage must be checked and removed. The hose connection to the siphon must be kept free from residue. A backflow preventer can be fitted. However, with siphons with two parallel hose connections and two fitted backflow preventers, it is necessary to check for reciprocal influence during pressure-release. Example: A thermo disinfectant and Vacuklav with simultaneous pressure / water discharge. 5. Visual control of the connections at the Y piece for possible visible blockage. Remove the drain hose on the rear of the device. Examine the screw connection with a torch. Guide a flexible cable section of c. 2.5mm ² through the screw connection from outside until it is visible inside the PTFE hose. In this way it is possible to recognize a blockage and possibly rectify it. Alternatively, you can replace the Y piece. 6. Blockage / a blocked plunger means that sufficient condensate drain can no longer be guaranteed. The SV condensate drain can be opened in the device with a 24 mm spanner. Check the positions of the plunger and SV and clean. It could be necessary to change the SV. 7. A permeable cooler is necessary for error-free vacuum generation and pressure release (rinsing). Blocked coolers initially result in extended operating times, especially whilst generating the vacuum. Remove the cooler and rinse with water. A clogged cooler should be changed. 8. To check the conduits, open the hose connection to the hose filter at the SV condensate drain. Activate the DC OUT 08 SV cooler drainage in the service menu. You can check for permeance by blowing lightly into the cooling coil in the effluent tank. Air movement can be felt at the SV exit. When blocked, the corresponding screws or hoses may need to be changed or unblocked. 9. Imprecise measurement by the pressure sensor (sensor defective outside the area of the sensor) can result in the target pressure pressure-release of 1300mbar not being reached. In the service menu under analogue inputs AIN 07, the P sensor is checked for plausibility. Normally, the surrounding pressure is around 1000mbar (+/- 50mbar). Upon considerable deviation, the sensor should be replaced
E429	FUNF: Unknown fault			
E430	FPAR: Fault in the parameter data			
E433	FS8T Sterilisation temp.- difference between the theoret. temp. and temp sensor2 (AIN02=Display-sensor) is too large		1. Deviation between the theoretical and actual temperatures are produced by vacuum leakage or the presence of non-condensable gasses in the steam. 2. defective temperature sensor input	The recorded temperatures input 1 "chamber control", input 2 "chamber log" and the theoretical temperature, calculated from the pressure through inputs 7 "pressure sensor" are located during sterilization within the so-called Temperature band (2K). If one of the three temperatures falls out of this temperature band, it will produce the error message FS8T. The presence of air in the chamber (leakage) or non-condensable gasses (feed water) can result in temperature differences. 1. The CO ₂ level (non-condensable gas) in the feed water varies and in the rule is not relevant. In the conditioning process, the water is subjected to additional degassing before evacuation and is improbable as a cause of malfunction.

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
			3. defective analogue input temperature sensor on the electronics board 4. defective pressure sensor input AIN 07 5. defective analogue input pressure sensor on the electronics board	Air can enter the chamber following a vacuum leakage. The leakage should be determined by a vacuum test. 2./3. Faulty temperature sensor in the chamber display or imprecise temperature measurement outside the message "faulty sensor". In the service menu analogue inputs (AIN 02) a temperature of 160°C (infinite) is displayed upon a cable break and 0.0°C for a short circuit. Select service menu analogue inputs AIN, Disconnect sensor 02 chamber display from the electronic strip terminal X45 (lower position). A temperature of c. 160°C should be displayed. If this temperature is displayed, disconnect the sensor 01 chamber log from X44 (middle position) and plug it on X45 for comparison. If a plausible temperature is registered on the sensor AIN 02 according to the temperature previously displayed on sensor AIN 01, the temperature sensor input 02 chamber display is defective and needs to be replaced. If a temperature other than 160°C e.g. 70°C is displayed when removing X45, the analogue measurement on the electronics is faulty. The electronics must be replaced. WARNING After troubleshooting, the temperature sensor 03, jacket heating may not be swapped with another sensor. 4./5. pressure sensor is not connected or "open", cable interruption, sensor short circuit. In the service menu analogue inputs (AIN 07) a pressure of 0-10mbar is displayed upon a cable break and c. 3000mbar for a short circuit. The plausible absolute surrounding temperature with an open door is c. 1000mbar (+/- 50mbar), depending on the location and weather conditions. Select the service menu input AIN. Disconnect sensor 07 chamber pressure sensor from the electronic strip terminal X50. A pressure of c. 0-10mbar should be displayed. If this value is displayed, the analogue processing in the electronics is fully functional. The pressure sensor must be replaced. If the value displayed is higher, e.g. 200mbar or higher, when the X50 has been disconnected, the analogue measurement on the electronics is defective. The electronics must be replaced.
E434	FS5T: Sterilisation overheat temp.sensor2 (AIN02=display-sensor)		Overheat produced by: 1. defective/ penetrated steam generator relays 2. defective temperature sensor input 2 "chamber control" 3. defective analogue input temperature sensor on electronics board 4. defective pressure sensor / input 7 5. defective analogue input pressure sensor on electronics board	1. Heating is controlled via semi-conductor relay outputs. The main heater produces 2300W at 230V, which corresponds to 10A, the regular heating 700W at 230V, which corresponds to 3A. Faulty steam generator relays (semi-conductor) can result in the output steam generator AC-OUT 01 or AC-OUT 02 being supplied with current. Low voltage results in overheating. In the service menu, select outputs AC-OUT. The steam generator main heater 01 is switched at "0". Measure the voltage on the electronics AC OUT 01 (X39/1-bl/red). You should not be able to measure any voltage. Switch on AC OUT 01, switch state "1". A voltage of c. 230V AC is measured. After renewed switching off, no more voltage will be measured. If voltage can be measured e.g. > 3 V when switched at "0", this would indicate defective semi-conductor relays. The electronics must be replaced. Similar to the main heating, the steam generator regular heating 02 should also be checked AC OUT 02 (X41/2-bl/red). 2./3. The temperature sensor is not connected or "open", current break, short circuit at sensor. In the service menu analogue inputs (AIN 01) a temperature of 160°C (infinite) is displayed with a cable break and 0.0°C for a short circuit. Select service menu analogue inputs AIN, Disconnect sensor 02 chamber display from the electronic strip terminal X45 (lower position). A temperature of c. 160°C should be displayed. If this temperature is displayed, remove the sensor 01 chamber control from X44 (middle position) and plug it on X45 for comparison. If a plausible temperature is registered on the sensor AIN 02 according to the temperature previously displayed on sensor AIN 01, the temperature sensor input 02 chamber display is defective and needs to be replaced. If a temperature other than 160°C e.g. 70°C is displayed when removing X45, the analogue measurement on the electronics is faulty. The electronics must be replaced. WARNING After troubleshooting, the temperature sensor 03, jacket heating may not be swapped with another sensor. 4./5. pressure sensor is not connected or "open", cable interruption, sensor short circuit. In the service menu analogue inputs (AIN 07) a pressure of 0-10mbar is displayed upon a cable break and c. 3000mbar for a short circuit. The plausible absolute surrounding temperature with an open door is c. 1000mbar (+/- 50mbar), depending on the location and weather conditions. Select the service menu input AIN. Disconnect sensor 07 chamber pressure sensor from the electronic strip terminal X50. A pressure of c. 0-10mbar should be displayed. If this value is displayed, the analogue processing in the electronics is fully functional. The pressure sensor must be replaced. If the value displayed is higher, e.g. 200mbar or higher, when the X50 has been disconnected, the analogue measurement on the electronics is defective. The electronics must be replaced.

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E435	FS6T: Sterilisation insufficient temp.sensor2 (AIN02=display-sensor)		Insufficient temperature produced by: 1. Insufficient temperature is produced by vacuum leakage or the presence of non-condensable gasses in the steam. 2. mech. defective / blocked steam valve 3. blocked / kinked steam conduit 4. Bad heating performance resulting from defective heating in the double jacket (tubular heating element) or charred cable 5. defective steam generator relays 6. pressure loss via leaky / burst conduits on the jacket or chamber 7. defective pressure sensor / input7 8. defective analogue input pressure sensor on electronics board.	The recorded temperatures input 1 "chamber control", input 2 "chamber log" and the theoretical temperature, calculated from the pressure through inputs 7 "pressure sensor" are located during sterilization within the so-called Temperature band (2K). If one of the three temperatures falls out of this temperature band, it will produce the error message FS1T. The presence of air in the chamber (leakage) or non-condensable gasses (feed water) can result in temperature differences. 1. The CO₂ level (non-condensable gas) in the feed water varies and in the rule is not relevant. In the conditioning process, the water is subjected to additional degassing before evacuation and is improbable as a cause of malfunction. Air can enter the chamber following a vacuum leakage. The leakage should be determined by a vacuum test. ATTENTION: WHEN PERFORMING ANY WORK ON THE DOUBLE JACKET (FEED PUMP; SOLENOID VALVE; MEASUREMENT TURBINE; CAPILLARY. SENSOR, TUBULAR HEATING ELEMENT) IT IS NECESSARY TO ENSURE THAT IT IS DEPRESSURIZED -DANGER OF SUFFERING BURNS - 2. Steam flows over the SV steam from the double jacket to the chamber. Blockage or a blocked plunger means that steam conduction can no longer be guaranteed. The SV can be opened in the device with a 24 mm spanner. WHEN UNDERTAKING ANY WORK ON THE STEAM VALVE; ENSURE THAT THE DOUBLE JACKET IS FIRST DEPRESSURIZED: SEE E 192 / WA2. Check the positions of the plunger and SV and clean. It could be necessary to change the SV. 5. Heating is controlled via semi-conductor relay outputs. The main heater produces 2300W at 230V, which corresponds to 10A, the regular heating 700W at 230V, which corresponds to 3A. Faulty steam generator relays (semi-conductor) can result in the output steam generator AC-OUT 01 or AC-OUT 02 being supplied with current. Low voltage results in overheating. In the service menu, select outputs AC-OUT. The steam generator main heater 01 is switched at "0". Measure the voltage on the electronics AC OUT 01 (X39/1-bl/red). You should not be able to measure any voltage. Switch on AC OUT 01, switch state "1". A voltage of c. 230V AC is measured. After renewed switching off, no more voltage will be measured. If voltage can be measured e.g. > 3 V when switched at "0", this would indicate defective semi-conductor relays. The electronics must be replaced. Similar to the main heating, the steam generator regular heating 02 should also be checked AC OUT 02 (X41/2-bl/red). 7.8.. The pressure sensor is not connected or "open", current break, short circuit at sensor. In the service menu analogue inputs (AIN 07) a pressure of 0-10mbar is displayed upon a cable break and c. 3000mbar for a short circuit. The plausible absolute surrounding temperature with an open door is c. 1000mbar (+/- 50mbar), depending on the location and weather conditions. Select the service menu input AIN. Disconnect sensor 07 chamber pressure sensor from the electronic strip terminal X50. A pressure of c. 0-10mbar should be displayed. If this value is displayed, the analogue processing in the electronics is fully functional. The pressure sensor must be replaced. If the value displayed is higher, e.g. 200mbar or higher, when the X50 has been disconnected, the analogue measurement on the electronics is defective. The electronics must be replaced.
E436	FS1A: The Air-Detector 2/3 detected air in the chamber during sterilization. Therefore, the started programme was interrupted.	1. A foreign object between the door seal and the flange on the chamber. 2. Door seal is damaged	1. Air detector temperature sensor AIN 5 is defective 2. blocked pressure-release filter in chamber	1. The temperature sensor is not connected or "open", current break, short circuit at sensor. In the service menu analogue inputs (AIN 01) a temperature of 160°C (infinite) is displayed with a cable break and 0.0°C for a short circuit. Select service menu analogue inputs AIN, sensor 01. Disconnect the chamber control from electronic strip terminal X44 (lower position). A temperature of c. 160°C should be displayed. If this temperature is displayed, disconnect the sensor 02 chamber log from X45 (middle position) and plug it on X44 for comparison. If a plausible temperature is registered on the sensor AIN 01 according to the temperature previously displayed, the temperature sensor input 01 chamber control is defective and needs to be replaced. If a temperature other than 160°C e.g. 70°C is displayed when removing X44, the analogue measurement on the electronics is faulty. The electronics must be replaced. WARNING After troubleshooting, the temperature sensor 03, jacket heating may not be swapped with another sensor. 2. The pressure-release filter connected or screwed on the rear base of the chamber is impermeable, resulting e.g. from packaging residue, or is choked. To determine the fault, perform cleaning or sterilization without the filter.

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
E437	M437: Please wait Jacket must first be allowed to cool for the vacuum test with the solenoid valve "steam."	1. If the special vacuum test with an open steam valve should be carried out, the jacket temperature may not be too high so that re- evaporation does not occur during the test.		
E438	WVAL: Validation necessary	Warning message in a situation similar to the maintenance message, after a specified number of charges has been reached, you will be requested to perform the next validation. This message can be issued in relation to the number of charges and the number of operating hours (with or without standby) depending on the mode set.		Perform a validation procedure. The autoclave can still be started. The validation message will continue to appear until the procedure has been carried out successfully.
E439	FK1: Blockage pressure- release during conditioning. With type 40-B, 44-B, 14-B, check effluent hose for kinks or pinching. Check for packaging residue in the chamber,	1. Kinked effluent hose 2. Bad installation (multiple devices on a single siphon, sagging effluent hose) 3. Pressure release clips blocked by packaging residue	1. Filter in pressure release pipe is blocked (all devices) 2. Blockage in clips pressure release on the chamber (all devices) 3. Mechanically defective or blocked pressure-release valve (all devices) 4. Blockage in the effluent conduit to siphon (all water-ring devices) 5. Blockage / calculations on the effluent Y piece (all water-ring devices) 6. Mechanically defective / blocked condensate drain valve (all membrane devices) 7. Blocked cooler / clips on the cooler blocked (all membrane devices) 8. Blocked conduits from the condensate drain valve to cooling coil (all membrane devices) 9. Defective pressure sensor input 7	1. Open the hose filter SV pressure release and clean. 2. The pressure-release clips on the floor of the chamber must be free of residue or deposits. In particular, packaging residue can hinder pressure release. 3. Blockage / a blocked plunger means that sufficient pressure-release can no longer be guaranteed. The SV pressure release can be opened in the device with a 24 mm spanner. Check the positions of the plunger and SV and clean. It could be necessary to change the SV. 4. The hose between the effluent exit on the Vacuklav and the siphon must be kink-free. Obstruction resulting from blockage must be checked and removed. The hose connection to the siphon must be kept free from residue. A backflow preventer can be fitted. However, with siphons with two parallel hose connections and two fitted backflow preventers, it is necessary to check for reciprocal influence during pressure-release. Example: A thermo disinfectant and Vacuklav with simultaneous pressure / water discharge. 5. Visual control of the connections at the Y piece for possible visible blockage. Remove the drain hose on the rear of the device. Examine the screw connection with a torch. Guide a flexible cable section of c. 2.5mm ² through the screw connection from outside until it is visible inside the PTFE hose. In this way it is possible to recognize a blockage and possibly rectify it. Alternatively, you can replace the Y piece. 6. Blockage / a blocked plunger means that sufficient condensate drain can no longer be guaranteed. The SV condensate drain can be opened in the device with a 24 mm spanner. Check the positions of the plunger and SV and clean. It could be necessary to change the SV. 7. A permeable cooler is necessary for error-free vacuum generation and pressure release (rinsing). Blocked coolers initially result in extended operating times, especially whilst generating the vacuum. Remove the cooler and rinse with

E***	Incident text Display (BO v2.4x)	Possible operating cause	Possible device cause	Suggestion / Corrective
				water. A clogged cooler should be changed. 8. To check the conduits, open the hose connection to the hose filter at the SV condensate drain. Activate the DC OUT 08 SV cooler drainage in the service menu. You can check for permeance by blowing lightly into the cooling coil in the effluent tank. Air movement can be felt at the SV exit. When blocked, the corresponding screws or hoses may need to be changed or unblocked. 9. Imprecise measurement by the pressure sensor (sensor defective outside the area of the sensor) can result in the target pressure pressure-release of 1300mbar not being reached. In the service menu under analogue inputs AIN 07, the P sensor is checked for plausibility. Normally, the surrounding pressure is around 1000mbar (+/- 50mbar). Upon considerable deviation, the sensor should be replaced
E442	WEU0: Increased deviation of the eutectic temperature sensor		1. The temperature measured by the temperature sensor chamber control AIN 01 exhibits greater deviation from the melting temperature than the tolerance permits 2. Temperature sensor AIN 01 is imprecise 3. Analogue input AIN 01 on the electronics board is imprecise 4. Transfer resistance by the signal transmission and from the temperature sensor at the electronics input.	
E443	MEU1: Calibration procedure started			

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MELAG Vacuklav - Comfort Line (Input + Output Anschlüsse / Connections)


Type	PIN	Function	Bezeichnung/Description	Farbe/Color:
Alle / All	X 2/1-2	+24V-	Displaylüfter / Fan display	rt , bl
Alle / All	X 3/1-3	S	Eingang / Input 24V AC	sw
		+	Eingang / Input + 24V DC	rt
		-	Eingang / Input - 24V DC	rt-sw
Alle / All	X 15/1-6	Display 1	Bedienpanel / Controlpanel	gr
Alle / All	X 16/1-2	+24V-	LED CF-Card	rt , ws
Alle / All	X 25/1-6	EN 1	Ethernet 1 (Front)	gr
Alle / All	X 26/1-6	EN 2	Ethernet 2 (Rückseite / Rear)	gr
Alle / All	X 27/1-6	EN 3	Ethernet 3 (Rückseite / Rear)	gr
Alle / All	X 29/+	+	DCOUT + 24 V Versorgung / Supply	bn
Alle / All	X 29/1	DCOUT 1	MV Vakuum 1 / SV Vacuum 1	gr
Alle / All	X 29/2	DCOUT 2	MV Druckablass / SV Pressure release	rt
Alle / All	X 29/3	DCOUT 3	MV Dampf / SV Steam	vio-sw
Alle / All	X 29/4	DCOUT 4	MV Belüften / SV Ventilation	ge
Alle / All	X 29/5	DCOUT 5	MV Notablass / SV Emergency release	or
Alle / All	X 29/6	DCOUT 6	MV Spülen / SV Rinse	rs
13-B/41-B/43-B	X 29/7	DCOUT 7	MV Vakuum 2 / SV Vacuum 2	sw-gn
13-B/41-B/43-B	X 29/8	DCOUT 8	MV Kühlerentleerung / SV Cooler drain.	sw
Alle / All	X 29/9	DCOUT 9	MV Melasteam / SV Melasteam	rt-gr
Alle / All	X 29/10	DCOUT 10	MV Air Injektion / SV Air Injection	ws-rs
13-B/14-B	X 29/11	DCOUT 11	Türverriegelung manuell / doorlock manually	rt-gn
14-B/40-B/44-B	X 29/12	DCOUT 12	Gehäuselüfter / Fan housing	sw
Alle / All	X 30/1	1	Türmotor über Druckschalter Tür /	sw
	X30/2	2	Doormotor through pressure switch door	rt
Alle / All	X 37/1-2	230 V Trafo	Versorgung Ringtrafo / Supply Transformator	bl , bn
Alle / All	X 39/1	ACOUT 1	Phase Hauptheizung / Heating-main (2300 W)	rt
Alle / All	X 40/N	N	Neutral Haupt-Heizung / Heating-main (2300 W)	bl
Alle / All	X 41/2	ACOUT 2	Regelheizung / Heating regulation (700 W)	rt
Alle / All	X 41/3	ACOUT 3	Speisepumpe + MV / Feedpump +SV	gn
14-B/40-B/44-B	X 41/4	ACOUT 4	MV Kühlwasser / SV Cooling water	ws
13-B/41-B/43-B			Kühlerlüfter / Cooler fan	
14-B/40-B/44-B	X 41/5	ACOUT 5	Wasserringpumpe / Waterring pump	vio
13-B/41-B/43-B			Membran-Vakuumpumpe / Membrane vacuum pump	
Alle / All	X 41/N	N	ACOUT 2-5 Neutral Versorgung / Supply	bl
Alle / All	X 42/1	DINZ 1	Turbine Speisewasser / feedwater	bn , or , gn
14-B/40-B/44-B	X 42/2	DINZ 2	Turbine Kühlwasser / Cooling water	bn , ge , gn
Alle / All	X 42/3	DINZ 3		
Alle / All	X 43/1	DIN 1	Türschalter 1 / Doorswitch 1	rt
	X 43/2	DIN 2	Türschalter 2 / Doorswitch 2	ws
	X 43/3	DIN 3	Niveausensor Mantel / Jacket	rs
	X 43/4	DIN 4	Schwimmer / Float switch Aqua-Dem	bl
	X 43/5	DIN 5	Schwimmer Abwasser / Float switch drain	gr
	X 43/6	DIN 6	Schalter / Switch Melasteam (optional)	vio
	X 43/-	-	DIN - Versorgung / Supply	
	X 43/+	+	DIN + Versorgung / Supply	bn
Alle / All	X 44/1-2	AIN 1	Temp.-sensor Kammer Steuerung / Chamber control	bl , sw
Alle / All	X 45/1-2	AIN 2	Temp.-sensor Kammer - Protokoll / Chamber protocol	bl , sw
Alle / All	X 46/1-2	AIN 3	Temp.-sensor Mantel / Jacket	bl , sw
14-B/40-B/44-B	X 47/1-2	AIN 4	Temp.-sensor Vakuumpumpe / Vacuum pump	bl , br
Alle / All	X 50/1-4	AIN 7	Drucksensor Kammer / Pressuresensor chamber	gn , ws , bn
Alle / All	X 54/1-2	AIN 11	Leitwertsensor / Conductivity sensor	sw , sw

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