

GALIT

EXPERTS IN DENTISTRY



GALIT CONSOLE

EXPERTS IN DENTISTRY

GALIT is the premier manufacturer of dental equipment in Ukraine. It is distinguished by its integration of high-quality European components, commitment to universality, and emphasis on ergonomic design.

EXPERTISE IS OUR PRIORITY

For a quarter of a century, our mission has been to craft high-quality dental and medical equipment, delivering comprehensive solutions with an expert focus on addressing client needs. We prioritize the durability of our solutions and value enduring relationships with our customers and partners.

BRAND VALUES:

QUALITY – We adhere to international quality standards for medical devices.

TRUST – Respect, integrity, and sincerity guide all our actions.

CARE – Our commitment to our people and customers drives us to act with dedication and passion.

INNOVATION – We embrace the latest technologies to achieve the best outcomes.

PARTNERSHIP – We collaborate with healthcare professionals and partners to achieve mutual success.

DEDICATION – We work meticulously to deliver exceptional results while providing the highest customer service.



Why do dentists rely on the GALIT CONSOLE in their daily practice?

ERGONOMICS AND COMFORT

An integrated chair within the dental unit improves ergonomics and maximizes working space for both the dentist and the assistant. The open space beneath the chair allows easy interaction with the patient during treatment. The 180° swiveling cuspidor provides effortless access for the assistant while keeping the patient comfortable.

FUNCTIONALITY AND FLEXIBILITY

The robust and stable design of the unit allows safe mounting of additional equipment such as a microscope, X-ray, or monitor for efficient and precise work. The long arms of the dental light provide optimal illumination and convenience during treatment.

QUALITY AND RELIABILITY

The GALIT CONSOLE is made in Europe from premium European components and ensures lasting durability and reliable performance.

GALIT CONSOLE — an ergonomic and reliable dental unit that combines patient comfort with convenient working conditions for both the dentist and the assistant.



PATIENT CHAIR

Maximum comfort, minimum effort

The patient's chair is an essential and integral part of any dental unit. At GALIT, we prioritize patient comfort at every stage of treatment while enabling the dentist and assistant to focus entirely on their work.

All GALIT dental units feature premium ergonomic chairs manufactured by the Italian company Tecnodent. They stand out for their intelligent design, proven reliability, and comfort certified to European standards.



Customizable upholstery colours from a palette with more than 25 options. In addition to standard upholstery, Soft Foam or Memory Foam options are available for increased comfort.



The chair can feature a left-side and/or right-side armrest for maximum patient comfort.



Additional Trendelenburg function (simultaneous 15° leg lifting and 90° backrest reclining) for proper circulation.



Precision joystick for easy and effective hands-free chair control.



The chair can lift patients up to 210kg




9002 White	300 40 20 Lilac	3012 Peach	5003* Dark blue
5011 Steel blue	6001 Taupe	7502 Anthracite	4009* Mauve
3493 Turquoise	3016 Coral red	6033 Mint green	8017 UM33** Suede brown
1300T Zinc yellow	1011 Light gray	3491 Periwinkle blue	7021** Black gray satin
1414 Red	5515 Apple green	1604T Burgundy	15141*** Sparkle blue
3494 Silk gray	5023 Pastel blue	1007 Black	14966*** Sparkle green
1009 D Electric blue	5510 Purple		21366*** Platinum

*The product colour may differ from the actual item



DENTIST'S CONSOLE

Ergonomics based on experience



DENTIST'S CONSOLE WITH UPPER HOSE DELIVERY

Soft counterbalance springs ensure easy and comfortable instrument operation, preventing hand fatigue even during lengthy procedures.



DENTIST'S CONSOLE WITH LOWER HOSE DELIVERY

Easy instrument position and angle control for maximum comfort and ergonomics during the treatment.

CONTROL PANEL

The dentist's console is available with either a standard panel or an advanced touchscreen panel.

Standard control panel

Easy power adjustment for all instruments and patient chair. Everything is within reach for full control.



Touchscreen control panel

Intuitive interface, advanced functionality and real time display of all the parameters enable instant response to any changes, as well as more precise and efficient work.



INSTRUMENTS

The dentist's console of the GALIT dental unit can be equipped with up to five custom-tailored instruments, designed to meet modern dentistry standards. Arranged in a convenient and familiar order, the instruments ensure comfort and efficiency during procedures.

With a wide selection of instruments, dentists can tailor their workspace to any treatment, from prophylaxis to surgery. GALIT provides flexible options that enhance both functionality and ergonomics.



CUSPIDOR UNIT



The ceramic cuspidor bowl rotates 180°, improving assistant access and ensuring maximum patient comfort.

ASSISTANT'S CONSOLE

DÜRR Comfort assistant's console

Enables efficient interaction between the dentist and the assistant, which is the key to high-quality and comfortable treatment.

Main advantages:

- Instrument holder with adjustable angle.
- Easy horizontal movement of the console.
- Available with 2–5 instruments.
- Smooth and precise instrument positioning.



PRO assistant's console

The smart ergonomic design of the console allows for convenient instrument arrangement and easy patient chair control.

Main advantages:

- Available with up to 4 instruments.
- Independent operation of the aspiration system — the saliva and dust aspirators operate separately to enhance performance.
- Easy horizontal movement of the console.





DENTAL LIGHT


Craftsmanship in details

All GALIT dental units are equipped with lights from the market leader, the Italian company FARO. Each model delivers high-precision, even illumination, simple and intuitive control, and removable handles for autoclave sterilization. Designed for maximum reliability and patient comfort, FARO lights ensure optimal visibility for every procedure.

Operating temperature:




For accurate dental colour matching, prosthetics, and conservative dentistry




For orthodontic therapy, conservative dentistry, and children's dentistry



For surgery



Composave: prevents premature curing



MAIA LED

Light field: 170 x 95 mm.

Light flux intensity: 3,000–45,000 LUX at a distance of 700 mm.

5 000 K

EVA LED TUNABLE WHITE

Light field: 190 x 100 mm.

Light flux intensity: 3,000–50,000 LUX at a distance of 700 mm.

5 700 K

5 000 K

4 000 K

2 700 K



B75. Dentistry mode:

Light field: 200 x 100 mm.

Light flux intensity: 3,000–50,000 LUX at a distance of 700 mm.

5 700 K

5 000 K

4 000 K

2 700 K

B75. Surgery mode:

Light field: 250 mm.

Light flux intensity: Up to 160,000 LUX at a distance of 800 mm.

5 000 K

2 700 K





FOOT CONTROL

GALIT dental units can be equipped with foot controls featuring various functions, giving dentists the flexibility to select the setup that meets their specific needs.



WIRED MULTIFUNCTIONAL FOOT CONTROL:

- Switching dental instruments on/off.
- Switching spray and air/water functions of the instruments on/off .
- Adjusting instrument speed and power.



WIRELESS MULTIFUNCTIONAL FOOT CONTROL:

- Advanced foot control with wireless Bluetooth connection and high-capacity battery, providing extended service life of up to 6 months.



PNEUMATIC FOOT CONTROL:

- Switching dental instruments on/off.

Impressive quality. Inspiring innovations



At GALIT, every stage of a project is carefully managed to ensure outstanding results. For non-standard options, the process begins with a tailored technical design solution, ensuring that each client receives a product of fully guaranteed quality.

GALIT operates on 5,000 m² of modern production space. Our advanced facilities ensure maximum precision and product quality, while enabling fast and efficient fulfillment of orders.



With a professional design division, state-of-the-art equipment, top manufacturing experts, and partnerships with over 30 renowned dental equipment brands, GALIT provides expert solutions that meet the needs and expectations of even the most sophisticated clients.

We place particular emphasis on the durability of our dental units, as we highly value the longevity of our relationships with our clients.

GALIT – EXPERTS IN DENTISTRY



GALIT

Manufacturer

PC Galit
6E, 15 Kvitnya Str., Baikivtsi, Ternopil region,
47711, UKRAINE

+380-800-502-998

office@galit.te.ua
service@galit.te.ua

European representative

Galdent CZ s.r.o.
K.Vapence 448, 41503 Teplice, Czech Republic

+420-775-471-378

export@galit.te.ua
galdentcz@gmail.com



Version: 20251223

The technical specification or appearance of the equipment is subject to change without prior notice